



# Protekt® Aire Geri-Chair Overlay System

## Alternating Pressure/Low Air Loss Recliner Mattress & Pump Overlay System

This 100% USA made alternating pressure recliner geri-chair air mattress system allows for residents who are at high risk for pressure ulcers to be positioned in a variety of mobile devices. Residents are no longer confined to bed due to risks of sitting without appropriate pressure redistribution. Protekt Aire Geri-Chair Overlay Systems can be integrated into a geriatric chair, recliner, or standard high-back wheelchair. Thirty one bladders gently alternate, providing low interface pressures in addition to comfortably supporting the resident.

- Designed to deliver proven alternating pressure relief therapy to patients while in the reclining geri-chair.
- Low air loss top sheet has a special multi-chambered air distribution layer which administers low air loss relief directly underneath the patient.
- User friendly quiet pump provides 8 Liters of air per minute.
- On demand low air loss.
- Chair pad is covered with a four way stretch quilted nylon comforter specifically designed to prevent friction and shear.
- Fits most standard geriatric recliners and is easily transportable from room to room, 69"x20"x3" (XL system is 69"x28"x3").
- Three adjustable placement straps with buckles to secure to any stretcher or chair.
- 1 Year warranty.

Item #	Description	Size of Chair Pad	Weight Capacity
80018	Protekt Aire Geri-Chair Overlay System – 20"	69"x20"x3"	250 lb
80019	Protekt Aire Geri-Chair Overlay System – 28"	69"x28"x3"	450 lb

[www.proactivemedical.com](http://www.proactivemedical.com)



270 Washington Street, Mount Vernon, NY 10553

Ph: 855-237-7622 • Fax: 845-205-6005