


# Protekt® Supreme Support

Self-Adjusting Air/Foam Mattress with  
Optional Alternating Pressure Pump

HCPCS  
E0373/  
E0277



## FEATURES AND BENEFITS:

- The Protekt® Supreme Support utilizes a series of 10 interconnected air cells and pressure valves to continually equalize pressure in response to patient movement and position changes. The Equalization Technology automatically redistributes weight to relieve pressure by increasing body surface contact.
- The top foam layer of Visco Gel Polymer Technology (VGPT™) is an open cell visco memory foam infused with gel polymers. This VGPT™ provides rapid heat dissipation, excellent pressure redistribution and vastly increases breathability compared to traditional foam.
- The sloped heel section reduces pressure and increases blood circulation by redistributing the load from the heel to the less vulnerable lower leg.
- Stretch nylon top cover is fluid resistant, breathable and fire retardant.
- The optional pump easily converts the Protekt® Supreme Support to a powered air mattress.
- Warranty: Mattress; 5 Year Non-Prorated.
- Optional Protekt® 3000 Pump; 18 Month.

94001	Non-Powered Self Adjusting	35"x80"x6"	500 lb.
94001P	Powered Self Adjusting	35"x80"x6"	500 lb.
94003	Non-Powered Self Adjusting	42"x80"x7"	600 lb.
94003P	Powered Self Adjusting	42"x80"x7"	600 lb.
94004	Non-Powered Self Adjusting	48"x80"x7"	650 lb.
94004P	Powered Self Adjusting	48"x80"x7"	650 lb.
94005	Non-Powered Self Adjusting	54"x80"x7"	700 lb.
94005P	Powered Self Adjusting	54"x80"x7"	700 lb.



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com