


4-WHEEL SCOOTER



OUTSTANDING
MANEUVERABILITY,
OUTDOOR PERFORMANCE
AND LUXURY FEATURES

The Pursuit® Sport 36 Volt Mobility Scooter offers the maneuverability of a 4-wheel mobility scooter along with high-performance features like 36 volts of power, full suspension and a hydraulic sealed brake system. The Pursuit Sport 36 Volt features sleek, aggressive style and a wealth of standard luxury touches including a full steering wheel with integrated digital display, making it the ultimate outdoor scooter.



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

www.pridemobility.com

SPECIFICATIONS

PURSUIT® SPORT 36 VOLT 4-wheel Model # MV714

Weight capacity:	400 lbs.
Maximum speed ¹ :	Up to 9 mph
Ground clearance ² :	5.25" at mid frame, 2.5" at rear frame
Turning radius:	83"
Overall length ² :	56"
Overall width ² :	27.25"
Front tires:	13" pneumatic
Rear tires:	13" pneumatic
Suspension:	Front and Rear
Range per charge ^{1,4} :	Up to 22.5 miles
Total weight without batteries ⁵ :	271 lbs.
Heaviest piece when disassembled:	216 lbs.
Standard seating:	Type: Deluxe High-back pillow top seat with sliders and headrest Weight: 55 lbs. Material: Black Vinyl Dimensions: Width: 19.5", Depth: 17.5"
Drive system:	Rear-wheel drive, sealed transaxle, 36-volt DC motor
Brakes:	Electro-magnetic with drum brake system
Battery requirements ³ :	(3) 12-volt, deep cycle Size: (3) 55 AH Weight: 38 lbs. each
Battery charger:	36V 8-amp, off-board

Color:



Silver



Red

- ¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- ² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- ³ AGM or gel cell type required.
- ⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- ⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

FEATURES

- Full suspension for optimal handling
- Deluxe reclining high-back seat with headrest and sliders
- Full steering wheel with integrated digital display
- 36 volt battery system with (3) 12-volt 55 AH batteries
- Full lighting package that includes directional signals and two lower headlights
- Throttle control selection switch allows control by the left or right hand
- Parking brake
- Increased horn volume
- Water bottle holder standard
- Two rearview mirrors standard



Full steering wheel with integrated digital display



Water bottle holder standard

Pride
Mobility Products Corp.
— LIVE YOUR BEST —

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113
www.pridemobility.com