


4-WHEEL SCOOTER



ENGINEERED FOR
POWERFUL PERFORMANCE
AND **EXCEPTIONAL** STYLE

The **Pursuit® XL** is an exciting blend of power and precision that features a 24-volt DC motor, full suspension, 13" low-profile tires and a hydraulic-sealed brake system. Standard luxury touches like an infinite adjustable wraparound delta tiller, high-back pillow top seat and LED lighting make it the perfect outdoor scooter.



Full Suspension



facebook.com/PrideMobilityProductsCorp



youtube.com/PrideMobility

PRIDE® **MOBILITY** SCOOTERS
LIVE YOUR BEST®

pridemobility.com

SPECIFICATIONS

PURSUIT® XL 4-wheel Model # SC714

Weight capacity:	400 lbs.
Maximum speed ¹ :	Up to 8.5 mph
Ground clearance ² :	5.25"
Turning radius:	82.5"
Overall length ² :	56"
Overall width ² :	27.25"
Front tires:	13" pneumatic
Rear tires:	13" pneumatic
Suspension:	Front and Rear
Range per charge ¹ :	Group 24 (70 AH): Up to 22.13 miles Group 27 (100 AH): Up to 29.51 miles
Total weight without batteries:	259.9 lbs. (without batteries)
Heaviest piece when disassembled:	N/A
Standard seating:	Type: High-back pillow top with sliders and headrest Weight: 55 lbs. Material: Black Vinyl Dimensions: 19.5" x 17.5"
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes:	Regenerative and electromechanical
Battery requirements ³ :	Size: Group 24 (75 AH), Group 27 (100 AH) Weight ⁴ : Group 24: 53 lbs. each, Group 27: 60 lbs. each
Battery charger:	8-amp, off-board

Color:



Silver



Candy Apple Red

FEATURES

- Infinitely-adjustable tiller angle with conveniently located handle
- Powerful drivetrain; features 200 amp microprocessor-based controller
- Deluxe reclining high-back seat with headrest and sliders
- Full lighting package that includes directional signals and lower headlight
- Delta tiller with wraparound handles
- Large 13" low-profile, pneumatic tires for excellent outdoor performance
- Front and rear suspension
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front basket standard
- Rearview mirror standard
- Hydraulic handbrake

OPTIONS

- Rear basket
- Safety flag
- Durable weather cover
- Cup holder
- Saddlebag
- Lap belt
- Weather cover
- Weather breaker canopy
- Cell phone holder

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com