




# MANUAL WHEELCHAIR

## REFERENCE GUIDE

A COMPREHENSIVE  
GUIDE TO SELECTING  
THE RIGHT QUICKIE  
MANUAL WHEELCHAIR



LIVE  
WITH  
OUT  
LIMITS



		RIGID					TILT-IN-SPACE			FOLDING							
		K0005					E1161			K0005					K0004	K0009	
		7 <sup>RS</sup>	7 <sup>R</sup>	5 <sup>R</sup>	QRi	GP / GPV	GPSA	IRIS / IRIS SE	SR45	TS / TS SE	Quickie 2 Lite	Quickie 2	Quickie 2 HP	QXi	LXI	LX	M6
																	
Design	Frame Style	Rigid - Open	Rigid - Open	Rigid - Open	Rigid - Open	Rigid - Traditional	Rigid - Swing-Away	Rotation in Space	Rotation in Space	Tilt	Folding	Folding	Folding	Folding	Folding	Folding	Folding
	Standard Weight†	17.7 lbs (Std. Config.) 9.9 lbs (Transport)	18.9 lbs (Std. Config.) 11.1 lbs (Transport)	21.5 lbs (Std. Config.) 12.9 lbs (Transport)	22 lbs (Std. Config.) 13.4 lbs (Transport)	GP: 23.75 lbs GPV: 23.6 lbs	26 lbs	39 lbs (w/footrests, 18x18)	39 lbs	49 lbs	25.5 lbs (fixed) 27.75 (swing-away)	28.9 lbs	27.9 lbs	28.8 lbs	29.6 lbs	29.6 lbs	53 lbs (w/o footrests)
	User Weight Capacity	265 lbs	265 lbs	265 lbs	265 lbs	250 lbs*	250 lbs*	250 lbs (Std) 350 lbs (HD)	265 lbs	250 lbs (Std) 350 lbs (HD)	265 lbs	265 lbs - STD 350 lbs - HD	250 lbs	265 lbs.	265 lbs	250 lbs	650 lbs
Adjustability	Seat Width	12" to 20"	12" to 20"	12" to 20" *	12" to 20"	12" to 22" *	12" to 22" *	14" to 22"	14" to 20"	14" to 20"	14" to 20"	12" to 22" *	11" to 22"	14" to 20"	12" to 20"	12" to 20"	22" to 30"
	Seat Depth	12" to 20"	12" to 20"	12" to 20" *	12" to 20"	12" to 18" * (GP) 12" to 22" * (GPV)	12" to 22"	14" to 22"	14" to 22"	14" to 22"	14" to 20"	12" to 20" *	10" to 20"	14" to 20"	12" to 20" *	12" to 18"*	18" to 22"
	Back Angle Adjustment	Fixed	Fore (-) 27° to Aft (+) 8°	Fore (-) 17° to Aft (+) 8°	Fore (-) 17° to Aft (+) 8°	Fore (-) 8° to Aft (+) 4°	Fore (-) 8° to Aft (+) 4°	Fore (-) 5° to Aft (+) 30°	Fore (-) 5° to Aft (+) 30°	Fore (-) 12° to Aft (+) 30°	Fore (-) 4° to Aft (+) 14°	Fore (-) 3° to Aft (+) 15°	Fore (-) 3° to Aft (+) 15°	Fore (-) 5° to Aft (+) 14°	Fore (-) 5° to Aft (+) 14°	Fore (-) 5° to Aft (+) 14°	N/A
	Camber	0°, 3°, 6°	0°, 3°, 6°	0°, 3°, 6°	0° to 3°	0° to 11°	0° to 11°	0°	0°	0°, 4°	0°, 2°, 4°, 6°	0°, 2°, 4°, 6°	0° to 7°	0°, 3°	1° to 5°	1° to 5°	0°
	Tension Adj. Upholstery	Yes	Yes	Yes	Yes	Yes	Yes	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	No
Clinical Applications	Front Seat-to-Floor Ht.	16" to 21"	16" to 21"	16" to 21"	15" to 20"	15.5" to 23"	15.5" to 23"	12.5" to 19"	14" to 20.5"	—	15.5" to 21"	13.5" to 22"	17" to 22"	12.5" to 19.5"	16" to 20.5"	16" to 20.5"	17" to 20"
	Rear Seat-to-Floor Ht.	13" to 20"	13" to 20"	13" to 20"	14" to 20"	14.5" to 20.5" (GP) 14.5" to 20.6" (GPV)	14.5" to 20.5"	12.5" to 19"	14" to 20.5"	—	13" to 19.5"	13" to 22"	15.5" to 21"	12" to 19"	15" to 19.5"	15" to 19.5"	17" to 20"
	Footrest/Frame Angle	75°, 80°, 85°	70°, 80°, 85°	75°, 80°, 85°	75°, 80°, 85°	70°, 80°, 85°	Swing-Away	60°, 70°, 80°, 90°, Elevating	70°, 80°, Elevating	60°, 70°, 80°, 90°, Elevating	70°, 80°	Swing-Away	60°, 70°	60°, 70°, 80°, Elevating	Swing-Away	Swing-Away	Swing-Away
	Fold-Down Back	No	Yes	Yes	Yes	Yes	Yes	Yes	Optional	No	No	No	No	No	No	No	No
	Adjustable Back Height	Fixed	Adj. 8" to 20"	Adj. 10" to 20"	Adj. 14" to 20"	Adj. 8" to 20"	Adj. 8" to 20"	Adj. 15" to 21" or 18" to 24"	Adj. 15" to 21" or 18" to 24"	No	9" to 21"	9" to 21"	8.5" to 20"	13" to 17" or 16" to 20"	14" to 19"	14" to 19"	16.75" to 20.13"
	Center of Gravity (COG)	0.5" to 4"	0" to 4"	-1.5" to 4"	1" to 3"	0.5" to 3.5"	0.5" to 3.5"	N/A	N/A	N/A	0.5" to 3"	-1" to 4.5"	0.5" to 3.5"	1" to 3"	0.5" to 3.5"	0"	N/A
	Tilt Range	—	—	—	—	—	—	40° or 55°	45°	53°	—	—	—	—	—	—	—
Lifestyle Options	Xtender™ Power Assist	No	Yes	Yes	No	Yes	Yes	No	No	No	No*	Yes	Yes	No	No	No	No
	Front Suspension	Yes	Yes	Yes	Yes	Built-4-Me	Built-4-Me	Built-4-Me	Built-4-Me	Built-4-Me	Yes	Yes	Yes	Built-4-Me	Built-4-Me	Built-4-Me	No
	Rear Suspension	No	Yes	Yes	No	No	No	No	No	No	No	No	No	No	No	No	No
	Transit	Yes°	Yes°	Yes°	Yes°	No	No	Yes°	Yes	No	Yes°	Yes°	Yes	Yes	Yes	No	No
	Warranty	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime	Limited Lifetime

\* Additional options available through Quickie Built-4-Me.

† Transport Weight = all items that can be removed without tools to depict transporting into or out of a vehicle.  
Standard Weight = 16" x 16" chair with wheel locks.

° Conforms with ANSI/RESNA WC-19 standards. Please refer to Quickie® & Zippie® WC-19 Transit Securement Supplement (MK-100072) for additional information. Supplement is available for download at [www.SunriseMedical.com](http://www.SunriseMedical.com)

## Quickie Colors

IMPORTANT NOTE:  
Print processes cannot accurately represent paint finishes. Actual colors may vary.




Aztec Gold




Black



Black Cherry



Black Opal



Blue Opal



Candy Blue



Candy Purple




Candy Red



Electric Blue



Evergreen




Glow




Green Apple



Hot Sparkle Pink



Mauve Pink



Rootbeer



Sparkle Silver



Sunrise Orange



Yellow



Titanium Color Paint



Matte Black



Matte Black Cherry



Matte Purple



Matte Electric Blue



Matte Evergreen



Desert Camo¹



Mossy Oak Camo¹



Stars and Stripes¹



Zebra¹²

IMPORTANT NOTE: The color options shown are for reference only. They are not intended as accurate representations. Refer to the Quickie Color Brochure or a Quickie Color Ring for representations of Quickie colors.

1. Kolorfusion™ available on select Quickie chairs.  
2. Black Zebra Kolorfusion™ pattern is available on any base color. Shown here on Candy Blue.



# JAY<sup>®</sup> Seating Selection

## JAY BACKS



**Basic**  
E2611



**GO™**  
E2611, E2612



**Care™**  
E2615



**J2™**  
E2613



**J2™ Deep Contour**  
E2620



**Focus Point**  
E2620



**J3™ PA**  
E2613, E2614




**J3™ PL**  
E2615




**J3™ PD, PDL  
& PDC**  
E2620

## JAY CUSHIONS



**Basic**  
E2601, E2602




**Soft Combi™ P**  
E2605, E2606



**Care™**  
E2607, E2608



**Zip™**  
E2607



**GO™**  
E2603, E2604



**Ion™**  
E2607, E2608



**Easy™**  
E2607, E2608



**Lite™**  
E2603, E2604,  
E2607, E2608




**Union™**  
E2607, E2608



**X2™**  
E2622



**J2™**  
E2622, E2623,  
E2624, E2625



**J2™ Deep Contour**  
E2622, E2623,  
E2624, E2625



**Fusion™**  
E2622, E2623



**J3™**  
E2622, E2623  
E2624, E2625

built4me

Our Built-4-Me Team is dedicated to creating product customizations beyond what is listed on our product order forms. Experts in the engineering and manufacturing of custom rehab products will listen to your specific requests and build a chair as individual as your fingerprint!