



Quickie GT

January 2018

EIR3

Mark For:

Submitting for:

☐ Quote

☐ Order

Date: _____

PO#: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: ☐ Fax ☐ Email

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HPCPS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions.

The use of **HPCPS CODES** does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

✓	Option #	HPCPS	Definition	Price	✓	Option #	HPCPS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. FRAME STYLE & SETUP

Base Model

For sizes 14" or less please use HPCPS E1235

☐	145WM13	K0005	GT (weight capacity 265 lbs)	\$2,800.00
--------------------------	---------	-------	------------------------------	------------

Frame Angle (see Figure 1)

☐	145FA1	75°	NC
☐	145FA8	80°	NC
☐	145FA2	85°	NC

Frame Inset Per Side (See Figure 2)

☐	145FF1	0"	NC
☐	145FF2	1" (N/A w/ 12" wide)	NC
☐	145FF3	2" (N/A w/ 12", 13", or 14" wide)	NC

Frame Inset	Footrest Width
0	Subtract 0 from Frame Width
1"	Subtract 2" from Frame Width
2"	Subtract 4" from Frame Width

Attitude Bike (Optional)

☐	145IN04	Attitude Bike bundle	See Order Form
--------------------------	---------	----------------------	----------------

Frame Style (See Figure 3)


XTR Suspension N/A with Xtender

☐	145FS34	Classic Frame	NC
☐	145FS35	Performance Frame	NC
☐	145FS36	XTR Suspension Frame Only available in Performance Frame Set Up	\$895.00
☐	145SK1	Light Shock (up to 180 lbs.)	NC
☐	145SK3	Heavy Shock (up to 265 lbs.)	NC


Depth Adjustable Backrest (Optional)

See **Figure 5** Depth Adjust Backrest extends the frame rear of the seat sling, select a factory setting. N/A on 19" and 20" Seat Depths.


☐	145DP02	0 Setting Backrest at standard position, field adjustable to -1" or -2"	\$165.00
☐	145DP03	-1 Setting Backrest -1" offset, field adjustable to flush or -2"	\$165.00
☐	145DP04	-2 Setting Backrest -2" offset, field adjustable to flush or -1"	\$165.00



0 POSITION



-1" POSITION



-2" POSITION

Figure 1




Figure 2





Figure 3



Classic Frame - Traditional GT seat sling to front frame bend angle dimensions.

Performance Frame - Seat Sling to front frame bend angle is 1.5" shorter than the classic position. Front of back cane post to front edge of seat sling will remain the same for both Classic and Performance.

1. FRAME STYLE & SETUP (CONT)

Seat Width

See Figure 2. Measure from outside of seat tube to outside of opposite seat tube. XTR and J3 Back N/A w/ 12" and 13" widths

☐	145S3		12"		NC
☐	145S4		13"		NC
☐	145S5		14"		NC
☐	145S6		15"		NC
☐	145S7		16"		NC
☐	145S8		17"		NC
☐	145S9		18"		NC
☐	145S10		19"		NC
☐	145S11	E2201	20"		\$500.00

Seat Depth

See Figure .1 Measures from front of back post to front edge of seat sling. Seat sling will be 1" shorter than depth ordered, but actual seat depth will remain the same as ordered.

☐	145SD5		12"		NC
☐	145SD6		13"		NC
☐	145SD7		14"		NC
☐	145SD8		15"		NC
☐	145SD9		16"		NC
☐	145SD10		17"		NC
☐	145SD11		18"		NC
☐	145SD12		19"		NC
☐	145SD13	E2203	20"		\$550.00

2. FRONT SEAT HEIGHT SET UP

Front Seat Height

See Figures 4 and 5

☐	145SH14	16"	NC	☐	145SH20	18.5"	NC
☐	145SH15	16.5"	NC	☐	145SH21	19"	NC
☐	145SH16	17"	NC	☐	145SH23	19.5"	NC
☐	145SH17	17.5"	NC	☐	145SH24	20"	NC
☐	145SH18	18"	NC	☐	145SH26	20.5"	NC
☐				☐	145SH27	21"	NC

Front to Rear Seat Height Differences

Frame Depth	Max Seat Squeeze
12"/13"	2.5"
14"/15"	3.0"
16"/17"	3.5"
18"/19"	4.0"
20"	4.5"

Caster Wheels

☐	145CW1		3" Micro		NC
☐	145CW3		4" Low Profile Polyurethane		NC
☐	145CW31		4" Micro Lighted		\$70.00
☐	145CW72		4" x 1.25" Aluminum Polyurethane		\$110.00
☒	145CW44	E2219	4" x 1.25" Semi Pneumatic		\$55.00
☐	145CW64		4" x 1.5" Aluminum Soft Roll		\$185.00
☐	145CW5		5" Low Profile Polyurethane		NC
☐	145CW51		5" Micro Lighted		\$85.00
☐	145CW6	E2219	5" Soft Roll (Grey) N/A w/ Frog Legs		\$90.00
☐	145CW65		5"x 1" Aluminum Polyurethane		\$110.00
☐	145CW7	E2219	5"x 1.5" Semi Pneumatic		\$85.00
☐	145CW66		5"x 1.5" Aluminum Soft Roll		\$195.00
☐	145CW34		6" Polyurethane		NC
☐	145CW35	E2219	6"x 1.5" Semi Pneumatic		\$85.00
☐	145CW68		6"x 1.5" Aluminum Soft Roll		\$205.00

Caster Fork

See Figures 5 and 6

☐	145CF27		Short Fork		NC
☐	145CF45		Short Wide Fork (w/ 1.25" or 1.5" caster)		NC
☐	145CF32		Tall Wide Fork		NC

Caster Fork Suspension

☐	145CF53	E1015	Frog Legs 0 weight to 100 lbs		\$425.00
☐	145CF54	E1015	Frog Legs 101 lbs weight to 190 lbs		\$425.00
☐	145CF55	E1015	Frog Legs 191 lbs weight to 265 lbs		\$425.00

Figure 4




Figure 5

Front STFH	Fork			
	Short/Single	Short Wide	Tall Wide	Frog Legs
3"	16"-20"	16"-20.5"	16.5"-21"	17.5"-21"
4"	16"-20.5"	16.5"-21"	16.5"-21"	17.5"-21"
5"	17", 18", 19", 20", 21"	17"-21"	17"-21"	18.5"-21"
6"	NA	NA	18"-21"	19", 20", 21"

Figure 6

Caster Size	Fork	Tubular/Flat			Flared			Platform		
		Inset			Inset			Inset		
		0"	1"	2"	0"	1"	2"	0"	1"	2"
75° Frame Angle	3"	Short, Short W, Tall W			NA					
		Single Sided			NA					
		Froglegs			NA					
	4"	Short, Short W, Tall W			NA					
		Single Sided			NA					
		Froglegs			NA	2"				
	5"	Short, Short W, Tall W	NA		NA	NA				
		Single Sided	NA		NA	NA				
		Froglegs	NA		NA	NA		5"		
	6"	Tall Wide	NA		NA	NA		2"		
		Froglegs	NA	NA	NA	NA	1.5"	NA	.5"	
80° Frame Angle	3"	Short, Short W, Tall W			NA					
		Single Sided			NA					
		Froglegs			NA					
	4"	All NC Forks			NA					
		Single Sided			NA					
		Froglegs	NA		NA	NA				
	5"	Short, Short W, Tall W	NA		NA	NA		1"		
		Single Sided	NA		NA	NA		1"		
		Froglegs	NA		NA	NA		2"		
	6"	Wide	NA	1"	NA	NA		NA	1"	
		Froglegs	NA	NA	NA	NA	NA	NA	2"	
85° Frame Angle	3"	All NC Forks			NA					
		Single Sided			NA					
		Froglegs			NA					
	4"	Short, Short W, Tall W			NA			.5"		
		Single Sided			NA			.5"		
		Froglegs	NA		NA	NA		.5"		
	5"	Short	NA	.5"	NA	NA		2"	.5"	
		Short Wide/Tall Wide	NA	1"	NA	NA		2"	.5"	
		Single Sided	NA	1"	NA	NA		2"	.5"	
		Froglegs	NA	1.5"	NA	NA		NA	1.5"	
	6"	Tall Wide	NA	2"	NA	NA		NA	2"	
		Froglegs	NA	NA	1.5"	NA	NA	NA	NA	1.5"

☐ = Go
 ☐ = No Go
 # = "Soft" Go with limited seat angle
 Min seat angle listed in cell

3. FOOTPLATE OPTIONS

Footrest

Heel loop or leg strap not included, See Footplate Options Section

☐	145F228		Aluminum Tubular	
☐	145F229		Aluminum Tubular Long	NC
☐	145F230		Aluminum Tubular w/ Flat Cover	NC
☐	145F231		Aluminum Tubular Long w/ Flat Cover	nc
☐	145F234	K0040	Aluminum Platform (Ang. Adj.)	\$95.00
☐	145F235	K0040	Aluminum Platform Long (Ang. Adj.)	\$95.00
☐	145F278	K0040	High Mount Platform N/A w/ 12", 13" W with 1" or 2" Inset or 14" or 15" W with 2" insets.	\$95.00

Adaptive Platform Spacer

N/A with Tubular Footrest

☐	145AD1		Adaptive Platform 2" Spacer	\$95.00
☐	145AD2		Adaptive Platform 4" Spacer	\$95.00
☐	145AD3		Adaptive Platform 6" Spacer	\$95.00

Footrest Options

☐	145HL5	K0038	Leg Strap	\$35.00
☐	145DI01		Removable Snap On Cover N/A w/ Platform Footrests	\$60.00
☐	145DI02		Wear Pads N/A w/ Platform Footrests	\$50.00
☐	145A179	K0108	Positioning Plates for Platform Footrest	\$50.00
☐	145A54	K0038	Extra Leg Strap	\$30.00

4. REAR SEAT HEIGHT SET UP

Rear Seat Height

14" & 14.5" only available with 20" rear wheel. N/A with Amputee Axle Plate.

13" & 13.5" only available with 20" or 22" rear wheel. N/A with Amputee Axle Plate.

☐	145RH11		13"	NC
☐	145RH12		13.5"	NC
☐	145RH13		14"	NC
☐	145RH14		14.5"	NC
☐	145RH15		15"	NC
☐	145RH16		15.5"	NC
☐	145RH17		16"	NC
☐	145RH18		16.5"	NC
☐	145RH19		17"	NC
☐	145RH26		17.5"	NC
☐	145RH20		18"	NC
☐	145RH27		18.5"	NC
☐	145RH21		19"	NC
☐	145RH37		19.5"	NC
☐	145RH22		20"	NC

Rear Wheel Size

Ultra Lightweight Wheels only available with aluminum anodized handrim

Contact Built 4 Me for Spinerger color options

☐	145RW3		20" Mag Wheel	NC
☐	145RW8		22" Mag Wheel	NC
☐	145RW1		24" Mag Wheel	NC
☐	145RW90		20" Lite Spoke Wheel	\$200.00
☐	145RW91		22" Lite Spoke Wheel	\$200.00
☐	145RW92		24" Lite Spoke Wheel	NC
☐	145RW93		25" Lite Spoke Wheel	NC
☐	145RW94		26" Lite Spoke Wheel	\$80.00
☐	145RW77		22" Ultra Lightweight Spoke Wheel	\$250.00
☐	145RW78		24" Ultra Lightweight Spoke Wheel	\$250.00
☐	145RW28		22" Spinerger SPOX Black	\$820.00
☐	145RW42		24" Spinerger SPOX Black	\$820.00
☐	145RW96		25" Spinerger SPOX Black	\$820.00
☐	145RW44		26" Spinerger SPOX Black	\$820.00
☐	145RW32		22" Spinerger LX Black	\$1,025.00
☐	145RW88		24" Spinerger LX Black	\$1,025.00
☐	145RW99		25" Spinerger LX Black	\$1,025.00
☐	145RW40		26" Spinerger LX Black	\$1,025.00
☐	145RW50		24" Mountain	\$250.00

Xtender Power Assist

RSH 17.5"-19.5" with 22" wheels, RSH 18.5"-20" with 24" wheels 0 camber, N/A w/ Amputee Axle. 22" & 24" wheels only. Not available with XTR or standard rear suspension. Anti-Tips no charge with Xtender. Dual Drive mode available w/B4ME

☐	145XT6	E0986	22" Xtender Power Assist JWX-2 NICKEL	\$7,895.00
☐	145XT7	E0986	22" Xtender Power Assist JWX-2 LITHIUM	\$8,695.00
☐	145XT4	E0986	24" Xtender Power Assist JWX-2 NICKEL	\$7,895.00
☐	145XT5	E0986	24" Xtender Power Assist JWX-2 LITHIUM	\$8,695.00

Please select one of the following handrim options

☐	145XH5		Plastic Coated	NC
☐	145XH6		Stainless Steel	NC

Xtender - Other Options

☐	145XA10		Spare NiMH Power Pack - NICKEL	\$1,395.00
☐	145XA11		Spare Li-ion Power Pack - LITHIUM	\$1,795.00

Xtender Wheel Options

☐	145XA8	E2211 E2212	SPARE 22" Lite Spoke Pneu Alum	NC
☐	145XA9	E2211 E2213	SPARE 22" Lite Spoke Pneu Plastic	\$90.00
☐	145XA2	E2211 E2212	SPARE 24" Lite Spoke Pneu Alum	NC
☐	145XA4	E2211 E2213	SPARE 24" Lite Spoke Pneu Plastic	\$90.00

Rear Wheel Tire

See Availability Matrix

☐	145RT1		Pneumatic	NC
☐	145RT2	E2211 E2213	Pneumatic w/ Airless Insert N/A with Xtender	\$80.00
☐	145RT3		Full Poly	NC
☐	145RT6	E2211 E2212	High-Pressure Clincher	\$95.00
☐	145RT27	E2211 E2212	Knobby V- Trak	\$95.00
☐	145RT21	E2211 E2212	Schwalbe Marathon Plus Evolution	\$135.00
☐	145RT38	E2218	SHOX Performance Poly	\$95.00
☐	145RT12	E2211 E2212	Mountain Tire	\$95.00

4. REAR SEAT HEIGHT SET UP (CONT)

Rear Tire & Wheel Matrix

Tire	Wheel	Mag	Lite Spoke	Ultra Lightweight	Spinergy
Pneumatic		20,22,24	20,22,24,26		22,24,26
Pneumatic w/ Airless Insert		20,22,24	20,22,24		22,24
Full Poly		20,22,24	20,22,24		
High-Pressure Clincher			22,24,25,26	22,24	22,24,25,26
Knobby V Trak by Kenda			24		24
Schwalbe Marathon			24,25,26		24,25,26
SHOX Performance Poly		24	24,25		24,25
22" = 489 mm 24" = 540 mm 25" = 559 26" = 590					

Handrim

Aluminum and Plastic are not available on all wheels in all sizes. Natural Fits are N/A w/ 26".
Handrim Selection not applicable with Xtender Power Assist

☐	145HR1	Aluminum Anodized	NC
☐	145HR4	Plastic Coated N/A w/ 25" & 26" Lite Spoke Wheels	\$95.00
☐	145HR22	Natural Fit - Standard Grip	\$395.00
☐	145HR23	Natural Fit - Super Grip	\$425.00
☐	145HR30	Natural Fit LT - Standard Grip	\$415.00
☐	145HR31	Natural Fit LT - Super Grip	\$475.00
☐	145HR24	Spinergy Flex Rim Only available on Spinergy LX wheels 24" & 25"	\$650.00
☐	145HR8	8 Oblique Projections N/A w/ Kik tires, only available with Mag or Spoke	\$325.00
☐	145HR7	Omit Handrim N/A w/ Ultra Lightweight Wheels	NC

Axle

☐	145AX1	Quick Release Steel	NC
☐	145AX13	Quad Quick Release Only available with Mag wheels	\$125.00

Axle Plates

Selection not required for XTR Suspension

☐	145AL40	Standard Short (COG 1" to 4.5") Some limitations with STFH	NC
☐	145AL43	Standard Long (COG 1" to 4.5")	NC
☐	145AL37	E0959 Amputee (COG -1" to 2.5")	\$175.00
☐	145AL38	E1015 Suspension (COG-.5 to +3) N/A with XTR Suspension Frame	\$395.00
☐	145AL39	Xtender (COG 1")	NC

Rear Seat to Floor Height	Short Axle Plate- Wheel Size:				
	20	22	24	24 (559)	26
	13	14	15	15.5	16
	13.5	14.5	15.5	16	16.5
	15	16	17	17.5	18
	15.5	16.5	17.5	18	18.5
	16	17	18	18.5	19

Rear Suspension Rating

Selection not required for XTR Suspension

☐	145SV04	From 0 weight to 100 lbs (Yellow)	NC
☐	145SV05	From 101 lbs weight to 190 lbs (Red)	NC

Camber

0, 3, 6 available with std frame. 0, 4, 8 available with XTR frame

☐	145CI1	0°	NC
☐	145CI2	3°	NC
☐	145CI3	4°	NC
☐	145CI4	6°	NC
☐	145CI5	8°	NC

Center of Gravity

Measure from front of back post to center of rear axle. XTR Suspension Frame will be factory set at 1.5", Xtender 1" - 4"

☐	145CC11	-1"	NC
☐	145CC10	-.5"	NC
☐	145CC12	0"	NC
☐	145CC1	.5"	NC
☐	145CC2	1"	NC
☐	145CC3	1.5"	NC
☐	145CC4	2"	NC
☐	145CC5	2.5"	NC
☐	145CC6	3"	NC
☐	145CC7	3.5"	NC
☐	145CC8	4"	NC
☐	145CC9	4.5"	NC




Figure 7

5. BACKREST

Back Type

☐	145BT3	Non-Folding Backrest	NC
☐	145BT82	Folding Lock Down Angle Adj	NC

Height Adjustable Backrest

☐	145B110	10" - 14" Back Height	NC
☐	145B127	12" - 16" Back Height	NC
☐	145B111	14" - 18" Back Height	NC
☐	145B137	16" - 20" Back Height	NC

Back Upholstery

☐	145BU1	Standard Upholstery	NC
☐	145BU8	Adjustable	\$175.00
☐	145BU57	Vented	\$200.00
☐	145BU14	Omit Upholstery (Use w/ JAY Backs)	-\$15.00

JAY J3 Back

☐	145JB01	JAY J3 Back	See OF
--------------------------	---------	-------------	--------

Military Service Patches

Only available on tension adjustable upholstery and Jay 3 Back

☐	145VU1	US Marines Patch	\$90.00
☐	145VU2	US Army Patch	\$90.00
☐	145VU3	US Navy Patch	\$90.00
☐	145VU4	US Air Force Patch	\$90.00
☐	145VU12	US Coast Guard	\$90.00
☐	145VU5	Purple Heart	\$90.00
☐	145VU8	Purple Heart with Marines Patch	\$145.00
☐	145VU9	Purple Heart with Army Patch	\$145.00
☐	145VU10	Purple Heart with Navy Patch	\$145.00
☐	145VU11	Purple Heart with Air Force Patch	\$145.00
☐	145VU13	Purple Heart with Coast Guard	\$145.00

Push Handles

☐	145PH6	Integral Push Handles	NC
☐	145PH5	Fold-Down Push Handles N/A w/ 10-14" Back Height or Standard Upholstery.	\$100.00
☐	145PH2	Omit Push Handles	NC

6. ARMREST

Armrest - Padded Swing Away

Includes angle adjustable receiver hardware. Some limitations may exist w/ integral push handles, some RSHs and wheel sizes

☐	145AA18	E0973	Padded Swing-Away - Standard	\$210.00
--------------------------	---------	-------	------------------------------	----------

Armrest - Height Adjustable

☐	145AA11	E0973	Single Post Ht. Adjust. Std (8.5"-13.5")	\$285.00
☐	145AA29	E0973	Single Post Ht. Adjust. Low (7"-11")	\$285.00

Armpad - For Height Adjustable

Suggest using classic model for tray tables

☐	145AM23		Desk Length (10") Classic	NC
☐	145AM21		Desk Length (10") Contour	NC
☐	145AM24		Full Length (14") Classic	NC
☐	145AM22		Full Length (14") Contour	NC

7. Key Options & Accessories

Wheel Lock

Not all wheel locks available in all seat depths, wheel sizes and COG settings. XTR suspension only available w/ high mount push

☐	145WL1		High Mount Push	NC
☐	145WL2		High Mount Pull	NC
☐	145WL27		Ergo Scissor	NC
☐	145WL30		Compact	\$125.00
☐	145WL10		Omit Wheel Lock	NC
☐	145WL12		Do Not Mount	NC

Extension Handle

Only available w/ High Mount Push or Pull Wheel Locks

☐	145EH7	E0961	Extension handles (Pair)	\$70.00
--------------------------	--------	-------	--------------------------	---------

Grade Aid

Only available w/ high mount push or pull. N/A w/ clincher, poly, Schwalbe, or Xtender

☐	145GA1	E0974	Grade Aid	\$150.00
--------------------------	--------	-------	-----------	----------

Anti Tips

☐	145AT1	E0971	Anti Tips (Pair)	\$140.00
☐	145AT3		Omit Rear Anti Tip	NC

Seat Sling

☐	145SL1		Standard Seat Sling	NC
☐	145SL6		Omit Seat Sling	-\$15.00

Cushion

☐	145CU1	E0669	2" Cushion	\$120.00
--------------------------	--------	-------	------------	----------

JAY Cushions (Optional)

To order any JAY cushions on-chair, please see the JAY price list.

Cushion name or part number must be filled out with desired width and depth

☐	145JC20	JAY Bundled Cushion		
		Cushion Name or Part #		Price List
		Width:	Depth:	

Side Guards

All side guards N/A w/ Height Adj. Arms. Some interference may exist w/ padded SA arms at lowest settings

☐	145SG33	K0108	Plastic - Standard	\$140.00
☐	145SG29	K0108	Plastic - Low	\$140.00
☐	145SG12	K0108	Carbon Fiber N/A w/ push or pull wheel locks	\$450.00

Pack Accessories

☐	145A1		Quickie Backpack Black	\$65.00
☐	145A2		Quickie Kids Backpack Black	\$60.00
☐	145A8		Seat Pouch Black	\$50.00
☐	145A500		Gel Cover Hanger/Frame	\$105.00
☐	145A994		Gel Cover Swing Away Arm	\$105.00

Positioning Belt

☐	145A142	E0978	Belt 1.5" Auto Buckle Universal w/clamps	\$70.00
☐	145A12	E0978	Belt 2" Aircraft Buckle	\$110.00
☐	145A13	E0978	Belt 2" Aircraft Buckle Padded	\$120.00

Spoke Guards

N/A w/ Xtender Power Assist

☐	145A23	K0065	Quickie Black Spoke Guard	\$115.00
☐	145A302	K0065	Quickie Clear Spoke Guard	\$115.00

Miscellaneous

☐	145A32	E0950	Wheelchair Tray	\$250.00
☐	145IG6		Neoprene Impact Guards	\$85.00
☐	145A51	E2207	Crutch Holder Only available with Integral or Fold Down push handles	\$195.00
☐	145A85		Bolt-On Push Handles	\$115.00
☐	145A83		Tool Kit	\$45.00

7. Key Options & Accessories (Cont)

Standard Frame Colors

☐	145C140	Aztec Gold	NC
☐	145C2	Black <i>(Not avail as base w/Zebra)</i>	NC
☐	145C141	Black Cherry	NC
☐	145C10	Black Opal	NC
☐	145C38	Blue Opal	NC
☐	145C36	Candy Blue	NC
☐	145C14	Candy Purple	NC
☐	145C35	Candy Red	NC
☐	145C142	Electric Blue	NC
☐	145C143	Evergreen	NC
☐	145C126	Glow	NC
☐	145C144	Green Apple	NC
☐	145C147	Hot Sparkle Pink	NC
☐	145C148	Mauve Pink	NC
☐	145C103	Rootbeer	NC
☐	145C114	Sparkle Silver	NC
☐	145C145	Sunrise Orange	NC
☐	145C113	Titanium Color Paint	NC
☐	145C8	Yellow	NC

Matte Frame Colors

☐	145C116	Matte Black	NC
☐	145C151	Matte Black Cherry	NC
☐	145C150	Matte Electric Blue	NC
☐	145C149	Matte Evergreen	NC
☐	145C152	Matte Purple	NC

Custom Frame Color

☐	145C112	Desert Camouflage	\$295.00
☐	145C43	Mossy Oak Camouflage	\$295.00
☐	145C42	Stars & Stripes	\$295.00
☐	145C108	Zebra <i>(Choose above base color) Std with White unless above base paint is chosen N/A w/ Black</i>	\$295.00

Anodized Color Package

Camber tubes are titanium grey except for the black option. Frog Legs black only. Black only w/ amputee, Xtender, suspension axle plates.

☐	145EB02	Black	NC
☐	145EB06	Blue	\$160.00
☐	145EB05	Green	\$160.00
☐	145EB04	Purple	\$160.00
☐	145EB03	Red	\$160.00
☐	145EB01	Ti Grey	\$160.00



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA