

The Invacare® line of Reliant™ lifts was conceived to help reduce the possibility of caregiver back injury and to ensure dignity in patient handling. These high-quality lifts are comfortable, reliable and integral tools for staff and caregivers. Using the lift for patient handling can help reduce staff injuries, save you money and reinforce safety as a priority in your health care facility. Invacare Reliant lifts support up to 450 pounds, making potentially challenging transfers smooth for the resident and caregiver. The high lift range allows easy access to tubs, chairs, commodes, floors and beds (with at least 4.5" under-bed clearance). Along with functional and innovative lifts, Invacare also offers facility assessment plans to further aid in your safety management program. Our all-encompassing lift system also includes in-service training, videos and an ergonomic safety plan. Our goal is to partner with you in lowering costs and protecting your staff from injury. The sturdy and practical design of Reliant lifts provides the solution you need for success in your quest to provide safe and effective health care.



Ask us about our injury prevention program.



Invacare® Reliant™ 450 Hydraulic Lift

Model no. RHA450-1 Hydraulic Adjustable Base (shown), RHL450-1 Hydraulic Low Base Adjustable

- Pump mechanism allows transfers of up to 450 lb. with minimal effort
- Comfortable handle placement for ease of use
- Wide base for maximum stability
- Low maintenance
- Limited warranty includes three years on lift, one year on pump



Invacare® Reliant™ 450 Battery-Powered Lift with Low Base

Model no. RPL450-1 Power with Low Base

- Removable battery easily charges overnight on wall-mounted charging station
- Equipped with a quiet, long-life motor
- Wide lift range for easier transfer
- Wide base for maximum stability
- 4.5" Under bed clearance
- Limited warranty includes three years on lift, one year on battery

Invacare® Reliant™ 450

Model no. RPA450-1 Power with Adjustable Base



Reliant Lift Features

- Height range allows lifting from floor and high surfaces
- Low friction casters significantly improve rollability and contribute to caregiver safety
- Removable battery pack eliminates downtime and simplifies the charging process
- Manual lowering for instant response to power loss
- Electronic system includes two 24-volt batteries: one for lift, one for battery-charging station; a manual emergency lowering device; and an anti-entrapment feature if boom meets any resistance while lowering
- Available with power opening base option (model no. RPL450-2)
- Multi-grip push handle is angled to minimize wrist strain
- Base legs adjust easily and lock securely into open position with the padded shifter handle or optional power option
- Total pinch-point protection covers all moving parts for added safety
- Padded swivel bar with 360° rotation and six-point hook-up allows use of Invacare chain-free slings as well as older style slings with chains or straps
- Attractive neutral design blends well into any environment
- Covered base provides protection against wear, dirt and moisture
- Manufactured in an ISO13485 facility



Invacare® Reliant™ 450 Power Lift with Power Opening Low Base

Model no. RPL450-2

Available Parts Kits for Reliant™ Lifts

Description	Part no.	Description	Part no.
Reliant touch-up paint	1075222	Reliant low base replacement kit (manual)	1079810
Reliant replacement batteries	1078276	Reliant sling holder bag (four-pack)	1088911
Reliant swivel bar kits	1143629	Reliant training video: <i>Lifting with Ease</i>	1090823
Extra-wide Reliant swivel bars	1143630	Reliant training video: <i>Standing with Ease</i>	1090824