

Pride®

RAPTOR



**Now
Available**

Soar into Adventure!

Soar into Adventure with the all new **Pride Raptor!** With speeds up to 14 mph and a 400 lb. weight capacity, it's the perfect scooter for anyone and everyone. Travel greater distances with a per charge range up to 31 miles and experience a smoother ride over varied terrain with front and rear suspension.



Digital Display



Front LED Light Package



Rear LED Light Package



Under Seat Storage



**Fast speed, with
advanced safety features.**

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

pridemobility.com



SPECS

PRIDE RAPTOR

Weight Capacity	400 lbs.
Front Tire	3" x16"
Rear Tires	3" x16"
Maximum Speed ¹	14 mph
Ground Clearance ³	3.5" at center of frame
Turning Radius ³	60"
Length/Width ³	59"/24.75"
Seating ³	High back 20w x 18d (100° Limited Recline)
Seat-to-Ground	28"
Weight (w/o Batteries)	208 lbs.
Seat-to-Deck	19"
Battery Weight	53 lbs. each
Battery Requirements ^{4,5}	75 (Ah)
Per Charge Range ^{1,2} (Max.)	19 miles (up to 400 lbs. and full speed) 31 miles (up to 220 lbs. and at low speed)
Battery Charger	8 amp
Suspension	Front and rear
Weight of Heaviest Piece ^{3,6}	(unit does not disassemble)
Warranty	2 year limited

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or 3%).

³ AGM or gel cell type required.

⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

FEATURES AND BENEFITS

- Front and rear suspension
- Two rear-view mirrors
- Full LED lighting package with directional signals and brake lights that offer full visibility for increased safety
- Hand-brake with regenerative braking system
- Secure, lockable under seat storage compartment with seat support rod
- Easy to read LCD console which shows estimated speed, mileage and battery display
- Limited recline fold-flat seat with manual adjustment and angle/width adjustable armrest
- Sleek and solid design for enhanced stability
- Front storage compartment located on the tiller
- Tubeless pneumatic tires for a smooth ride

Accessories

- Cup holder
- Saddle bag
- Safety flag

*Military Patches Available



Color



Silver

(US) 800-800-8586

(Canada) 888-570-1113

pridemobility.com

Pride
Mobility Products Corp.
— LIVE YOUR BEST —