

2023 order form

Use dimension chart to select appropriate size.



User dimensions (inches)		E420 small	E430 large
Height		30–50	46–72
 Key user dimension: height Select the appropriate stander by the user's overall height. Choose the model that allows for growth. Important: User's weight must not exceed the maximum working load.			
Item dimensions (inches)		E420 small	E430 large
Board length and width		53 x 16½	73 x 19
Distance between laterals		7–13	9–16
Height when horizontal		28½	21–30
Height when vertical		54½	74
Base length and width		37½ x 27	46½ x 29½
Max. working load (lb)		100	250

		E420 small		E430 large	
		E40	\$3,130 ☐	E41	\$4,900 ☐
* Supine Stander includes padding, footboard, adjustable armrests, hip strap, chest strap, leg straps, one pair of laterals, one pair of head blocks and double-locking casters.					
Additional laterals (pair)	CHOOSE Small	E504	\$160 ☐	N/A	
	Large	N/A		E439	\$235 ☐
Tray		E404	\$315 ☐	E414	\$390 ☐
CHOOSE	Adjustable abduction wedge	E945	\$260 ☐	E945	\$260 ☐
	Round abductor	E933	\$110 ☐	E933	\$110 ☐
Collar (requires round abductor)		E934	\$90 ☐	E934	\$90 ☐
Sandals (pair)	CHOOSE Small	E422	\$100 ☐	N/A	
	Medium	E433	\$110 ☐	E433	\$110 ☐
	Large	N/A		E432	\$110 ☐
Wedges (pair, require sandals)	CHOOSE Small	E424	\$70 ☐	N/A	
	Medium	E435	\$75 ☐	E435	\$75 ☐
	Large	N/A		E434	\$80 ☐
Hand Anchor		K821	\$165 ☐	K821	\$165 ☐
Total		E420	\$ _____	E430	\$ _____

* Required