


Support Station dimensions

User dimensions (inches)		K710
Height		46 – 80
 Important: User height and weight must not exceed the maximum number.		
Item dimensions (inches)		K710
Trunk board size		24 x 24
Maximum trunk board height		45
Minimum trunk board height		26
Height of frame		65½
Width of frame		27½
Depth when trunk board horizontal		31
Depth when collapsed against wall		8
Max. working load (lb)		250

2023 order form



* Frame

		K710	
		K715	\$3,610 <input type="checkbox"/>
Support strap		K716	\$115 <input type="checkbox"/>
Kneeboard		K717	\$720 <input type="checkbox"/>
Handholds	Pair	K718	\$90 <input type="checkbox"/>
	Single	K719	\$45 <input type="checkbox"/>
Total		K710	\$ _____

* Required