

2023 order form


Use dimension chart to select appropriate size.



		K310 Standard (better for gait training)		K320 Low-base (fits under low beds)	
* Frame includes one battery, one battery charger, accessories tote and charger power cord. (International power cords available on request.)	CHOOSE	Frame with scale	K301 \$6,400 <input type="checkbox"/>	K321 \$6,400 <input type="checkbox"/>	
		Frame	K302 \$5,230 <input type="checkbox"/>	K322 \$5,230 <input type="checkbox"/>	
* Thigh straps (pair)	CHOOSE	Narrow	K308 \$170 <input type="checkbox"/>	K308 \$170 <input type="checkbox"/>	
		Wide	K338 \$170 <input type="checkbox"/>	K338 \$170 <input type="checkbox"/>	
Switch pole		K666 \$245 <input type="checkbox"/>	K666 \$245 <input type="checkbox"/>		
Arm supports (pair)	CHOOSE	Small arm prompts	K305 \$520 <input type="checkbox"/>	K305 \$520 <input type="checkbox"/>	
		Large arm prompts	K315 \$595 <input type="checkbox"/>	K315 \$595 <input type="checkbox"/>	
		Arm platforms	K397 \$390 <input type="checkbox"/>	K397 \$390 <input type="checkbox"/>	
		Arm platforms w/ handgrips	K396 \$440 <input type="checkbox"/>	K396 \$440 <input type="checkbox"/>	
Pelvic support	CHOOSE	Small	K306 \$265 <input type="checkbox"/>	K306 \$265 <input type="checkbox"/>	
		Medium	K316 \$298 <input type="checkbox"/>	K316 \$298 <input type="checkbox"/>	
		Large	K326 \$315 <input type="checkbox"/>	K326 \$315 <input type="checkbox"/>	
Hip positioner	CHOOSE	Small	K667 \$325 <input type="checkbox"/>	K667 \$325 <input type="checkbox"/>	
		Large	K668 \$355 <input type="checkbox"/>	K668 \$355 <input type="checkbox"/>	
Hip positioner pad	CHOOSE	Small	K541 \$110 <input type="checkbox"/>	K541 \$110 <input type="checkbox"/>	
		Large	K542 \$150 <input type="checkbox"/>	K542 \$150 <input type="checkbox"/>	
Additional battery (recommended for institutional use)		K307 \$155 <input type="checkbox"/>	K307 \$155 <input type="checkbox"/>		
Additional thigh straps	CHOOSE	Narrow single	K318 \$85 <input type="checkbox"/>	K318 \$85 <input type="checkbox"/>	
		Narrow pair	K308 \$170 <input type="checkbox"/>	K308 \$170 <input type="checkbox"/>	
		Wide single	K328 \$85 <input type="checkbox"/>	K328 \$85 <input type="checkbox"/>	
		Wide pair	K338 \$170 <input type="checkbox"/>	K338 \$170 <input type="checkbox"/>	
Swivel locks (pair)		K319 \$196 <input type="checkbox"/>	K324 \$196 <input type="checkbox"/>		
Total		K310 \$ _____	K320 \$ _____		

* Required

TRAM dimensions

User dimensions (inches)		K310 and K320
Girth		22 min – 55 max
 Key user dimension: girth Important: User's weight or height must not exceed 350 lb or 76 inches		
Item dimensions (inches)	K310 TRAM	K320 TRAM
Overall length	45½	Same
Overall width	27½–40	Same
Overall height	43½–58	Same
TRAM weight (lb)	70	Same
Turning diameter	50	Same
Floor to top of base leg for furniture clearance	6¾	4¾
Min. user armpit height (when standing or ambulating)	30½	Same
Max. user height (when standing or ambulating)	76	Same
Max. working load (lb)	350	Same

Pelvic support

User dimensions (inches)

Indicated by the red arrow below

Key user dimension:

half of thigh circumference

measured at the groin.

K306 small


K316 med.

K326 large

6 – 10

8 – 14

10 – 18



Product dimensions (inches)

K306 small

K316 med.

K326 large

a

8½

12

13½

b

16

20

25


c

9

12

14

Hip positioner

Product dimensions (inches) – without pad		
		
	K667 small	K668 large
a Crotch length	2½	3½
b Width	1¼	1¾
c Length	9	11