

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

[Submit Form](#)

Effective: 5/22/2023

Rogue2 K0005, E1235 **\$2745**
Rogue Standard Frame. 275 lb weight capacity.

Frame

Seat Width

See reference A on page 7

12"	N/C
13"	N/C
14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	N/C
20" (E2201)	\$550

Seat Depth

See reference B on page 7

14"	N/C
15"	N/C
16"	N/C
17"	N/C
18"	N/C
19"	N/C
20" (E2203)	\$550

Transit Option

Transit System	\$320
Cannot be retrofitted. Requires Expanded Width Caster Position.	
Omit	N/C

Adjustable Depth

Adds the capability to grow seat in depth

No Rear Adjustment	N/C
1/2"	N/C
1"	\$85
2"	\$165

Sling Position

In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 7.

0" (Performance)	N/C
1"	N/C
2" (Classic)	N/C

Seat Option

Standard Seat Sling	N/C
Aluminum Solid Seat Pan (E2231)	\$265

Front Frame Bend

Angle between front seat tube and ground. See reference D on page 7.

90°	\$150
85°	N/C
80°	N/C
75°	N/C
70°	N/C

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	Height Adjustable Short Back Post	11" ²	12"	13"	14"	
		Height Adjustable Medium Back Post	13"	14"	15"	16"	17"
		Height Adjustable Tall Back Post	17"	18"	19"	20"	
Fixed Height ¹	N/C	Fixed Height Straight Back Post	9" ²	10" ²	11" ²	12"	13"
			14"	15"	16"	17"	18"

¹ Not available with back upholstery

² Not available with Fold Down Push Handle

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1017.1 Pricing and specifications subject to change without notice.

Push Handles

Omit Push Handle	N/C
Standard Push Handle	\$95
Low Profile Push Handle	\$150
2" shorter than standard push handles	
Bolt-On Push Handle	\$175
Handles will be 1" - 2" above back height	
Fold Down Push Handle	\$110

Stroller Handle

Stroller Handles are removable and height adjustable

Single Stroller Handle	\$165
Double Stroller Handle	\$330
Not available with Push Handles	

Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

Standard +3° +6° +9° -3° -6° -9° -12° _____ N/C

Back Upholstery Options

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

Padded Velcro Adjustable with Standard Cover	N/C
Padded Velcro Adjustable with Air Mesh Breathable Cover	\$75
Tension Adjustable with Standard Cover (E2611)	\$350
Tension Adjustable with Air Mesh Breathable Cover (E2611)	\$390
Omit	(\$10)

Rigidizer Bar

Adjustable Height Rigidizer Bar	\$75
Non-Adjustable Height Rigidizer Bar	N/C
Not Available with Fixed Height Back Post	

US Military Patches

Army	Navy	Marines	Air Force	
Coast Guard	US Purple Heart			\$90
Combo Army/Purple Heart		Combo Navy/Purple Heart		
Combo Marines/Purple Heart		Combo Air Force/Purple Heart		
Combo Coast Guard/Purple Heart				\$145

Front Seat Height and Casters

Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend.

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 7.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 7.

Footrest Taper	Price	Caster Size	Front Seat Height													
			3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
Straight	N/C	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
		4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
		5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
		6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
1"	N/C	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"					
		4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
		5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
		6"	18"	18.5"	19"	19.5"	20"	20.5"								
2"	N/C	3"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"					
		4"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"						
		5"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"							
		6"	18"	18.5"	19"	19.5"	20"	20.5"								
3.5° V-Front	\$100	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
		4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
		5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
		6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		
7° V-Front	\$100	3"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
		4"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"
		5"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	
		6"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"		

Casters

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C				
	.75" Lighted Rollerblade	\$70				
	1" Poly	N/C				
	1" Poly Aluminum	\$110				
	1" Pneumatic	\$80				
	1.5" Poly	\$65				
	1.5" Soft Roll Aluminum	\$200				

Caster Position

See reference L on page 7

Optimized	N/C
Caster Placement minimizes the footprint of the chair	
Expanded Width	N/C
Caster arm is positioned approximately .75" wider than Optimized position	

Fork Options

Standard Fork	N/C
Single Sided Fork	\$185
Frog Legs II Suspension System	\$450
Lightweight shock absorbing fork. Not included in Colored Anodizing Package.	

Footrest

Footplate

Tubular "Open" Footplate	N/C
Not available with V-Taper	
Tubular Footplate with ABS Cover	\$30
Not available w/footrest widths greater than 15.5". See Reference G on page 7.	
Not available with V-Taper.	
Aluminum Angle Adjustable (K0040)	\$98
Aluminum Angle Adjustable with Riser (K0040)	\$165
High Mount Aluminum Angle Adjustable (K0040)	\$165
Not available with V-Taper	
Angle Adjustable Flip Under (K0040, K0037)	\$225
Not available with footrest widths less than 9.5". Not available with V-Taper.	
See Reference G on page 7.	
High Mount Angle Adjustable Flip Under (K0040, K0037)	\$270
Not available w/footrest widths greater than 15.5". See Reference G on page 7.	
Not available with V-Taper.	
Hybrid Angle Adjustable	\$90
Hybrid Angle Adjustable with Wings	\$130

Footrest Protector

Footrest Protector	N/C
Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable	
Flip Under or High Mount Aluminum Angle Adjustable Footplates	

Seat to Footrest Length

Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 7.

8"	8.5"	9"	9.5"	10"	10.5"
11"	11.5"	12"	12.5"	13"	13.5"
14"	14.5"	15"	15.5"	16"	16.5"
17"	17.5"	18"	18.5"	19"	19.5"

Rear Seat Height and Wheels

Wheel Size / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 7.

Wheel	Rear Seat Height															
20"	12.5" ¹	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"
22"	13.5" ¹	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"		
24"	14.5" ¹	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"				
25"	15" ¹	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"					
26"	15.5" ¹	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"						

¹ Available only with Low Profile tires

One Arm Drive

One Arm Drive - Quick Release (E0958)	\$1495
Available only with 20", 22" and 24" spoke wheels and quick release axles.	
Spoke wheel style provided is exclusive to OAD. Standard handrim will be Aluminum Anodized. Limitations exist based on rear wheel spacing, camber and armrest selections.	

Left Side Drive	N/C
Right Side Drive	N/C
Aluminum Drive Handrim	N/C
Poly Coated Drive Handrim	\$200

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Black Spokes	Red Spokes	Blue Spokes	White Spokes	N/C
--------------	------------	-------------	--------------	-----

Wheels

No selection needed with One Arm Drive

Rear Wheel Type		Rear Wheel Size							
		Price	20"	22"	24"	25"	26"		
			RCW	RCW	RCW	RCW	RCW		
MAXX Spoke		N/C	H11 RCW	H11 RCW	H11 RCW	H11 RCW	H11 RCW		
Mag		N/C	H12 RCW	H12 RCW	H12 RCW	H12 RCW	H12 RCW		
Mag (Color) ²		\$105	H13 RCW	H13 RCW	H13 RCW	H13 RCW	H13 RCW		
Superlight Spoke		\$315			HE2 RCW	HE2 RCW	HE2 RCW		
Spinergy Spox		\$920		HE4 RCW	HE4 RCW	HE4 RCW	HE4 RCW		
Spinergy LX		\$1150		HE5 RCW	HE5 RCW	HE5 RCW	HE5 RCW		

¹ Available only with Superlight and Plastic Coated Handrims

² Color selection on page 6

Rear Tire Type

Full Profile

Street Pneumatic 1 3/8" (E2211, E2212)	N/C
<i>Not available on 25" wheel size</i>	
Full Poly 1 3/8" (E2218)	N/C
<i>Not available on Superlight wheel or 25" or 26" wheel sizes</i>	
Pneumatic with Airless Insert 1 3/8" (E2211, E2213)	\$85
<i>Not available on Superlight wheel or 25" wheel size</i>	
Schwalbe HS-127 1 3/8" (E2211, E2212)	\$125
<i>Available only on 24" wheel size. Not available on Mag wheel</i>	
Mountain 2" Black (E2211, E2212)	\$105
<i>Available only on 22", 24" and 25" MAXX wheels and Spinergy wheels</i>	

Low Profile

Low Poly 1" (E2218)	\$75
<i>Not available on Superlight wheel, 25" or 26" wheel sizes</i>	
Sentinel Solid Tire 1"	\$95
<i>Not available on Superlight wheel</i>	
High Pressure 1" (E2211, E2212)	\$95
<i>Not available on Mag wheel</i>	
Schwalbe One 1" (E2211, E2212)	\$125
<i>Available only on 24" and 25" wheel sizes. Not available on Mag wheel</i>	
Schwalbe Marathon Plus Evolution 1" (E2211, E2212)	\$140
<i>Not available on Mag wheel, 20" or 22" wheel sizes</i>	

Handrims

No selection needed with One Arm Drive

Standard Handrims

Aluminum Anodized	N/C
Superlight Aluminum Anodized	\$50
<i>Available only on 24" and 25" wheel sizes. Available only on MAXX Spokes, Superlight or Spinergy wheels.</i>	
Plastic Coated	\$115
Projection (8 Vertical Projections)	\$350
<i>Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels</i>	
Omit Handrim	N/C

Ergonomic Handrims

Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.

Lite Regular Grip	\$465
Standard Regular Grip	\$445
Lite Extra Grip	\$515
Standard Extra Grip	\$495
Lite Omit Grip	\$395
Standard Omit Grip	\$375

Handrim Options

Handrim Mount

If no selection is made, Standard Mount will be used

Close Mount	\$25
-------------	------

Micro-Grip

3M Micro Grip	\$65
<i>Not available on Plastic Coated or Projection Handrims</i>	

Axles

Quick Release	N/C
Quad Quick Release	\$185
<i>Not available with Mag or Spinergy wheels</i>	

Center of Gravity Presets

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 7.

Amputee Positions				
-1.5"	-1"	-.5"	0"	.25"
Standard Positions				
.5"	.75"	1"	1.25"	1.5"
1.75"	2"	2.25"	2.5"	2.75"
3"	3.25"	3.5"	3.75"	4"
4.25"	4.5"	4.75"	5"	

Rear Wheel Spacing

1" minimum with 0° camber. 1.5" minimum with T-Arms. The following apply if One Arm Drive is selected: 1.75" minimum with T-Arms. 1.25" minimum with select Rear Seat Heights. See Reference K on page 7.

0.75"	1"	1.25"	1.5"	1.75"	N/C
-------	----	-------	------	-------	-----

Rear Wheel Camber

The following apply if One Arm Drive is selected: 6° and 8° not available. 4° not available with T-Arms.

0°	2°	4°	6°	8°	N/C
----	----	----	----	----	-----

Wheel Locks

Wheel Lock

Push to Lock	N/C
Push to Lock Flush Mount	N/C
Pull to Lock	N/C
Under Seat Scissor Lock	\$35
Low Profile Scissor Lock	\$65
<i>Not available with Low Poly, Full Poly or Sentinel Solid tires</i>	
Short Thro Scissor Lock	\$160
Omit	N/C

Wheel Lock Options

Extension Handles (E0961)	\$84
<i>Available only with Push to Lock or Pull to Lock Wheel Locks</i>	
Grade Aids (E0974)	\$150
<i>Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires</i>	

Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk	\$310
	Full	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk	\$340
	Full	
Height adjustable T-Arm Low (6.5"-10.5") E0973 ¹	Desk	\$340
	Full	

¹ If requested configuration is unavailable, standard height T-Arm will be used. 14" minimum seat depth with Offset Frame.
² Includes direct mount receiver

Armrest Type	Price
Swing Away (8.5"-12.5") E0973	\$225
Swing Away with Side Guard (8.5"-12.5") E0973 ²	\$285
Omit Armrests	N/C

Side Guards

Direct Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Composite Side Guards - Direct Mount Receiver (K0108)	\$145
Aluminum Straight - Direct Mount Receiver - Black Matte (K0108)	\$195
Aluminum Straight - Direct Mount Receiver - Match Frame Color (K0108)	\$245
<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver - Black Matte (K0108)	\$275
<i>Not available with Adjustable Depth greater than 1/2"</i>	
Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver - Match Frame Color (K0108)	\$325
<i>Not available with Adjustable Depth greater than 1/2". Not available with Hydrographics.</i>	
Carbon Fiber Adjustable Fender - Direct Mount Receiver (K0108)	\$485
<i>Not available with Adjustable Depth greater than 1/2"</i>	

Clamp-On Mount Side Guards

Not available with T-Arms or Swing Away with Side Guard Armrests

Composite Side Guards - Clamp-On Receiver (K0108)	\$175
Aluminum Straight - Clamp On Receiver - Black Matte (K0108)	\$225
Aluminum Straight - Clamp-On Receiver - Match Frame Color (K0108)	\$275
<i>Not available with Hydrographics</i>	
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black Matte (K0108)	\$305
Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Match Frame Color (K0108)	\$355
<i>Not available with Hydrographics</i>	
Carbon Fiber Adjustable Fender - Clamp-On Receiver (K0108)	\$515

Accessories

Anti-Tips

Rear Anti-Tips (E0971)	\$165
User Activated - Pair	\$340
User Activated - Left	\$170
User Activated - Right	\$170
Omit	N/C

Calf Strap

Velcro Adjustable (K0038)	\$42
Padded Velcro Adjustable (K0038)	\$68
Bodypoint Padded AeroMesh (K0038)	\$118

Pelvic Positioning Belt

1 1/2" Auto Buckle Pelvic Positioning Belt (E0978)	\$85
1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978)	\$170
Small Bodypoint Evoflex (16-24") (E0978)	\$188
Medium Bodypoint Evoflex (19-29") (E0978)	\$188
Large Bodypoint Evoflex (25-37") (E0978)	\$188

Belt Mounting Kit

Ki Mobility Belt Mounting Clamp	\$35
Bodypoint Belt Mounting Bracket Kit	\$108

Other Accessories

Tool Kit	\$50
Removable Underseat Pouch	\$50
Ki Mobility Backpack	\$65
Luggage Carrier	\$85
Black Neoprene Impact Guards	\$90
Under Chair Carrier and Removable Bag	\$85
<i>Not available with Seat Widths less than 14" or scissor locks</i>	
Cane and Crutch Holder	\$195

Color

Options	Price	Colors																					
		Gloss												Matte									
		Black	Blue Skies	Bubblegum Pink	Candy Blue	Candy Purple	Candy Red	Copper	Granny Smith	MRGallic Orange	Silver	Teal	Yellow	Charcoal Gray	Emerald Green	Iridescent Burgundy	Iridescent Navy	Military Green	Retro Red	Royal Blue	Shadow Black	Smooth White	
Frame	N/C																						
Backrest Release Bar	\$50																						
	N/C	Black Anodized																					
MAXX Mag Color																							

Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

Color		Black ¹	Blue	Bright Green	Gold	Orange	Pewter	Pink	Purple	Red
Retail		N/C	\$25	\$25	\$25	\$25	\$25	\$25	\$25	\$25
Option	Back Plates ²									
	Camber Tube ³									
	Towers									
	Side Guard Receiver									
	Caster Fork ³									

¹ Black is no charge

² Not available with Transit Option

³ Available only with Standard or Single Sided

References

A. Seat Width *See Diagram 1*

Measured from outside of frame tube on one side to outside of frame tube on other side.

B. Seat Depth *See Diagram 3*

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

C. Sling Position *See Diagram 2*

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1" or 2" seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

D. Front Frame Bend *See Diagram 3*

Angle between front seat tube and ground.

E. Footrest Taper *See Diagram 1*

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

F. Front Seat Height *See Diagram 1*

Measured from floor to top of seat tube at front of seat sling.

G. Footrest Width *See Diagram 1*

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1" Y Taper	2" Y Taper	3.5° V Taper	7° V Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
12"	9.5"	7.5"		7.5"	
13"	10.5"	8.5"		8.5"	
14"	11.5"	9.5"		9.5"	
15"	12.5"	10.5"		10.5"	
16"	13.5"	11.5"	9.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"	15.5"	13.5"

H. Seat to Footrest Length *See Diagram 3*

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

I. Rear Seat Height *See Diagram 3*

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

J. Center of Gravity Preset *See Diagram 3*

Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post.

K. Rear Wheel Spacing *See Diagram 1*

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

L. Caster Position *See Diagram 4*

The Optimized Position keeps the casters tighter to the frame for enhanced maneuverability. Expanded Position pushes the casters out .75" past the Optimized Position for increased stability.

Diagram 1




Diagram 2




Diagram 3





Diagram 4



Axiom G

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
	Retail	\$135	\$84	\$135	\$84	\$135	\$84	\$135	\$84	\$135	\$84	\$135	\$84	\$135	\$84	\$223	\$84
	HCPC	E2601		E2601		E2601		E2601		E2601		E2601		E2601		E2602	
Depth	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
22"																	

Axiom S

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
	Retail	\$229	\$96	\$229	\$96	\$229	\$96	\$229	\$96	\$229	\$96	\$229	\$96	\$229	\$96	\$343	\$96
	HCPC	E2603		E2603		E2603		E2603		E2603		E2603		E2603		E2604	
Depth	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
22"																	

Axiom P

		Width															
		14"		15"		16"		17"		18"		19"		20"		22"	
		Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover	Cushion	Outer Cover
	Retail	\$397	\$107	\$397	\$107	\$397	\$107	\$397	\$107	\$397	\$107	\$397	\$107	\$397	\$107	\$515	\$107
	HCPC	E2605		E2605		E2605		E2605		E2605		E2605		E2605		E2606	
Depth	14"																
	15"																
	16"																
	17"																
	18"																
	19"																
	20"																
	21"																
22"																	

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0236.7 Pricing and specifications subject to change without notice.

Axiom SP - Visco

		Width																	
		14"			15"			16"			17"			18"			19"		
		Cushion		Cover	Cushion		Cover	Cushion		Cover	Cushion		Cover	Cushion		Cover	Cushion		Cover
		Inner	Outer		Inner	Outer		Inner	Outer		Inner	Outer		Inner	Outer		Inner	Outer	
Depth	Retail	\$481		\$118	\$481		\$118	\$481		\$118	\$481		\$118	\$481		\$118	\$481		\$118
	HCPC	E2607			E2607			E2607			E2607			E2607			E2607		
	14"																		
	15"																		
	16"																		
	17"																		
	18"																		
	19"																		
	20"																		
	21"																		
	22"																		

Axiom SP - Fluid

		Width																	
		14"			15"			16"			17"			18"			19"		
		Cushion		Cover	Cushion		Cover	Cushion		Cover	Cushion		Cover	Cushion		Cover	Cushion		Cover
		Inner	Outer		Inner	Outer		Inner	Outer		Inner	Outer		Inner	Outer		Inner	Outer	
Depth	Retail	\$492		\$118	\$492		\$118	\$492		\$118	\$492		\$118	\$492		\$118	\$492		\$118
	HCPC	E2607			E2607			E2607			E2607			E2607			E2607		
	14"																		
	15"																		
	16"																		
	17"																		
	18"																		
	19"																		
	20"																		
	21"																		
	22"																		

Solid Seat Insert

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	21"
		Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion	Cushion
Depth	Retail	\$68	\$68	\$68	\$68	\$68	\$68	\$68	\$68
	14"	XI1414	XI1514	XI1614					
	15"	XI1415	XI1515	XI1615	XI1715				
	16"	XI1416	XI1516	XI1616	XI1716	XI1816	XI1916	XI2016	XI2216
	17"	XI1417	XI1517	XI1617	XI1717	XI1817	XI1917	XI2017	XI2217
	18"	XI1418	XI1518	XI1618	XI1718	XI1818	XI1918	XI2018	XI2218
	19"	XI1419	XI1519	XI1619	XI1719	XI1819	XI1919	XI2019	XI2219
	20"	XI1420	XI1520	XI1620	XI1720	XI1820	XI1920	XI2020	XI2220
	21"		XI1521	XI1621	XI1721	XI1821	XI1921	XI2021	XI2221
	22"		XI1522	XI1622	XI1722	XI1822	XI1922	XI2022	XI2222

Star Stabil-Air Cushions

Star Stabil-Air Cushion (3" Tall Cells) w/ Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
Retail		\$549	\$69	\$80	\$549	\$69	\$80	\$549	\$69	\$80	\$549	\$69	\$80	\$549	\$69	\$80	\$675	\$91	\$103
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Height	13.5"																		
	15.5"																		
	17.25"																		
	19.25"																		
	21"																		

Star Stabil-Air Cushion (4" Tall Cells) w/ Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
Retail		\$549	\$69	\$80	\$549	\$69	\$80	\$549	\$69	\$80	\$549	\$69	\$80	\$549	\$69	\$80	\$675	\$91	\$103
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Height	13.5"																		
	15.5"																		
	17.25"																		
	19.25"																		
	21"																		

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0903.1 Pricing and specifications subject to change without notice.

Star Standard Air Cushions

Star Standard Air Cushion (2" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
Retail		\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$580	\$91	\$103
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Height	13.5"																		
	15.5"																		
	17.25"																		
	19.25"																		
	21"																		

Star Standard Air Cushion (3" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"					
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover				
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent			
Retail		\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80			
HCPCS		E2622			E2622			E2622			E2622			E2622					
Height	13.5"																		
	15.5"																		
	17.25"																		
	19.25"																		
	21"																		

Star Standard Air Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
Retail		\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$580	\$91	\$103
HCPCS		E2622			E2622			E2622			E2622			E2622			E2622		
Height	13.5"																		
	15.5"																		
	17.25"																		
	19.25"																		
	21"																		

Star Standard Air Cushions Dual Valve

Star Standard Air Cushion Dual Valve (4" Tall Cells) and Extra Standard / Incontinent Cover

		Width														
		13.5"			15.5"			17.25"			19.25"			21"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
			Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent		Standard	Incontinent
Retail		\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80	\$475	\$69	\$80
HCPCS		E2622			E2622			E2622			E2622			E2622		
Height	13.5"															
	15.5"															
	17.25"															
	19.25"															
	21"															

Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" -17" (40.5cm -43.0cm)	18" -20" (45.5cm -50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)



Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	DEX
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	BTX
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	FXT

Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

Axiom Posterior Back (PB)

Includes Lumbar Pad

		Width					
		12"	14"	16"	18"	20"	22"
Retail		\$610	\$610	\$610	\$610	\$610	\$918
HCPC		E2613	E2613	E2613	E2613	E2613	E2614
Height	Low	8"	AXPB1208	AXPB1408	AXPB1608	AXPB1808	AXPB2008
		10"		AXPB1410	AXPB1610	AXPB1810	AXPB2010
		12"	AXPB1212	AXPB1412	AXPB1612	AXPB1812	AXPB2012
	Mid	14"		AXPB1414	AXPB1614	AXPB1814	AXPB2014
		16"	AXPB1216	AXPB1416	AXPB1616	AXPB1816	AXPB2016
		18"		AXPB1418	AXPB1618	AXPB1818	AXPB2018

Back Hardware

DEX Hardware BTX Hardware FXT Hardware _____ N/C

Lateral Support for AXPB Backs

Pair	Description	Item Code	HCPC	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PBLF3x4PP	E0956	\$305
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PBLF4X5PP	E0956	\$305
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLS3X4PP	E0956, E1028	\$650
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLS4X5PP	E0956, E1028	\$650
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLA53X4PP	E0956, E1028	\$790
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLA53X4PP	E0956, E1028	\$790

Axiom Posterior Lateral Back (PL)

Includes Lumbar Pad

		Width					
		12"	14"	16"	18"	20"	22"
Retail		\$704	\$704	\$704	\$704	\$704	\$978
HCPC		E2615	E2615	E2615	E2615	E2615	E2616
Height	Low	8"	AXPL1208	AXPL1408	AXPL1608	AXPL1808	AXPL2008
		10"		AXPL1410	AXPL1610	AXPL1810	AXPL2010
		12"	AXPL1212	AXPL1412	AXPL1612	AXPL1812	AXPL2012
	Mid	14"		AXPL1414	AXPL1614	AXPL1814	AXPL2014
		16"	AXPL1216	AXPL1416	AXPL1616	AXPL1816	AXPL2016
		18"		AXPL1418	AXPL1618	AXPL1818	AXPL2018
	Tall	20"	AXPL1220	AXPL1420	AXPL1620	AXPL1820	AXPL2020
		22"		AXPL1422	AXPL1622	AXPL1822	AXPL2022

Back Hardware

DEX Hardware BTX Hardware FXT Hardware _____ N/C

Axiom Lateral Back (LB)

Does not include Lumbar Pad

		Width					
		12"	14"	16"	18"	20"	22"
Retail		\$851	\$851	\$851	\$851	\$851	\$989
HCPC		E2620	E2620	E2620	E2620	E2620	E2621
Height	Low	8"	AXLB1208	AXLB1408	AXLB1608	AXLB1808	AXLB2008
		10"		AXLB1410	AXLB1610	AXLB1810	AXLB2010
		12"	AXLB1212	AXLB1412	AXLB1612	AXLB1812	AXLB2012
	Mid	14"		AXLB1414	AXLB1614	AXLB1814	AXLB2014
		16"	AXLB1216	AXLB1416	AXLB1616	AXLB1816	AXLB2016
		18"		AXLB1418	AXLB1618	AXLB1818	AXLB2018
	Tall	20"	AXLB1220	AXLB1420	AXLB1620	AXLB1820	AXLB2020
		22"		AXLB1422	AXLB1622	AXLB1822	AXLB2022

Back Hardware

DEX Hardware BTX Hardware FXT Hardware _____ N/C

Lateral Support for AXPL & AXLB Backs

Pair	Description	Item Code	HCPC	Pair
	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PLLF3x4PP	E0956	\$305
	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PLLF4x5PP	E0956	\$305
	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLS3x4PP	E0956, E1028	\$650
	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLS4x5PP	E0956, E1028	\$650
	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLAS3x4PP	E0956, E1028	\$790
	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLAS4x5PP	E0956, E1028	\$790

Axiom Deep Lateral Back (DL)

Does not include Lumbar Pad

		Width					
		12"	14"	16"	18"	20"	22"
Height	Retail	\$851	\$851	\$851	\$851	\$851	\$989
	Mid	14"	AXDL1214	AXDL1414	AXDL1614	AXDL1814	AXDL2014
		16"		AXDL1416	AXDL1616	AXDL1816	AXDL2016
		18"	AXDL1218	AXDL1418	AXDL1618	AXDL1818	AXDL2018
	Tall	20"	AXDL1220	AXDL1420	AXDL1620	AXDL1820	AXDL2020
		22"		AXDL1422	AXDL1622	AXDL1822	AXDL2022

Back Hardware

DEX Hardware BTX Hardware FXT Hardware _____ N/C

Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	CS08SS	\$339	CS08MA	\$360	CS08PL	\$506	E0955 / E1028
Comfort Soft Headrest 10"	CS10SS	\$339	CS10MA	\$360	CS10PL	\$506	E0955 / E1028
Comfort Soft Headrest 14"	CS14SS	\$339	CS14MA	\$360	CS14PL	\$506	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	CP08SS	\$339	CP08MA	\$360	CP08PL	\$506	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	CP09SS	\$339	CP09MA	\$360	CP09PL	\$506	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	CP10SS	\$339	CP10MA	\$360	CP10PL	\$506	E0955 / E1028
Tri-soft Headrest 10"	TS10SS	\$339	TS10MA	\$360	TS10PL	\$506	E0955 / E1028
Tri-soft Headrest 14"	TS14SS	\$339	TS14MA	\$360	TS14PL	\$506	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	UMP	\$50	K0108

Not required for Axiom Headrests

Anterior Thoracic Supports

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

Anterior Thoracic Supports		
Description	Item Code	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	14941	\$121
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	14942	\$121
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	14943	\$121
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	14924	\$100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	14927	\$100
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	14928	\$100
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	14851	\$107.80
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	14852	\$107.80
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	14853	\$107.80
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	14854	\$107.80
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	14856	\$107.80
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	14857	\$107.80
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	14858	\$107.80
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	14859	\$107.80
AirLogic Anterior Trunk Support: Standard, Structure, Small	14282	\$164.17
AirLogic Anterior Trunk Support: Standard, Structure, Medium	14283	\$164.17
AirLogic Anterior Trunk Support: Standard, Structure, Large	14284	\$164.17
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	14285	\$164.17
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	14257	\$164.17
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	14258	\$164.17
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	14259	\$164.17
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	14260	\$164.17

Accessories

Adjustable Shoulder Guides (K0108)	\$95
<i>Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.</i>	
Extra Lumbar Pad (K0108)	\$50
<i>Available only with Axiom Lateral Back and Axiom Deep Lateral back</i>	
Privacy Flap (K0108)	\$50

Axiom Back Modifications

Axiom Back Modifications

Extended lead times will apply for all modifications. Makes back non-refundable.

Anterior Thoracic Mounting Modification
Provides holes for Anterior Thoracic Supports and Shoulder Guides

N/C

Gloss Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Black Back Shell	N/C
Blue Skies Back Shell	N/C
Bubblegum Pink Back Shell	N/C
Candy Blue Back Shell	N/C
Candy Purple Back Shell	N/C
Candy Red Back Shell	N/C
Copper Back Shell	N/C
Granny Smith Back Shell	N/C
Metallic Orange Back Shell	N/C
Yellow Back Shell	N/C
Silver Back Shell	N/C
Teal Back Shell	N/C

Matte Finish

If no selection is made, standard Shadow Black will be used with no additional lead time

Military Green Back Shell	N/C
Charcoal Gray Back Shell	N/C
Emerald Green Back Shell	N/C
Iridescent Burgundy Back Shell	N/C
Iridescent Navy Back Shell	N/C
Retro Red Back Shell	N/C
Royal Blue Back Shell	N/C
Smooth White Back Shell	N/C