



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



Company:

Account Number:  Date:

Quote:  PO Number:

ATP:

Marked For:

Associated Quote/Order#:

Ship To:

Address:

City:  State:

Zip:  Phone:

Email/Fax:

**SubmitButton**

Effective: 1/10/2024

**Rogue2** K0005, E1235 \_\_\_\_\_ **\$2873**  
*Rogue Standard Frame. 275 lb weight capacity.*

**Frame**

Seat Width	
<i>See reference A on page 7</i>	
<input type="checkbox"/> 12"	N/C
<input type="checkbox"/> 13"	N/C
<input type="checkbox"/> 14"	N/C
<input type="checkbox"/> 15"	N/C
<input type="checkbox"/> 16"	N/C
<input type="checkbox"/> 17"	N/C
<input type="checkbox"/> 18"	N/C
<input type="checkbox"/> 19"	N/C
<input type="checkbox"/> 20" (E2201)	\$564

Seat Depth	
<i>See reference B on page 7</i>	
<input type="checkbox"/> 14"	N/C
<input type="checkbox"/> 15"	N/C
<input type="checkbox"/> 16"	N/C
<input type="checkbox"/> 17"	N/C
<input type="checkbox"/> 18"	N/C
<input type="checkbox"/> 19"	N/C
<input type="checkbox"/> 20" (E2203)	\$564

Transit Option	
<input type="checkbox"/> Transit System	\$340
<i>Cannot be retrofitted. Requires Expanded Width Caster Position.</i>	
<input type="checkbox"/> Omit	N/C

Adjustable Depth	
<i>Adds the capability to grow seat in depth</i>	
<input type="checkbox"/> No Rear Adjustment	N/C
<input type="checkbox"/> 1/2"	N/C
<input type="checkbox"/> 1"	\$87
<input type="checkbox"/> 2"	\$169

Sling Position	
<i>In Performance Position (0"), sling is at front of frame. The 1" or 2" sling position extends the frame in front of the sling by that amount. See reference C on page 7.</i>	
<input type="checkbox"/> 0" (Performance)	N/C
<input type="checkbox"/> 1"	N/C
<input type="checkbox"/> 2" (Classic)	N/C

Seat Option	
<input type="checkbox"/> Standard Seat Sling	N/C
<input type="checkbox"/> Aluminum Solid Seat Pan (E2231)	\$272

Front Frame Bend	
<i>Angle between front seat tube and ground. See reference D on page 7.</i>	
<input type="checkbox"/> 90°	\$154
<input type="checkbox"/> 85°	N/C
<input type="checkbox"/> 80°	N/C
<input type="checkbox"/> 75°	N/C
<input type="checkbox"/> 70°	N/C

## Backrest

Back Post / Height Range / Back Height							
Back Post	Price	Height Range	Back Height				
Height Adjustable	N/C	<input type="checkbox"/> Height Adjustable Short Back Post	<input type="checkbox"/> 11" <sup>2</sup>	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	
		<input type="checkbox"/> Height Adjustable Medium Back Post	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"
		<input type="checkbox"/> Height Adjustable Tall Back Post	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	
Fixed Height <sup>1</sup>	N/C	<input type="checkbox"/> Fixed Height Straight Back Post	<input type="checkbox"/> 9" <sup>2</sup>	<input type="checkbox"/> 10" <sup>2</sup>	<input type="checkbox"/> 11" <sup>2</sup>	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"
			<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"

<sup>1</sup> Not available with back upholstery

<sup>2</sup> Not available with Fold Down Push Handle

### Push Handles

- Omit Push Handle N/C
- Standard Push Handle \$97
- Low Profile Push Handle \$154  
*2" shorter than standard push handles*
- Bolt-On Push Handle \$179  
*Handles will be 1" - 2" above back height*
- Fold Down Push Handle \$113

### Stroller Handle

*Stroller Handles are removable and height adjustable*

- Single Stroller Handle \$169
- Double Stroller Handle \$338  
*Not available with Push Handles*

### Rigidizer Bar

- Adjustable Height Rigidizer Bar \$77
- Non-Adjustable Height Rigidizer Bar N/C  
*Not Available with Fixed Height Back Post*

### US Military Patches

- Army  Navy  Marines  Air Force
- Coast Guard  US Purple Heart \$92
- Combo Army/Purple Heart  Combo Navy/Purple Heart
- Combo Marines/Purple Heart  Combo Air Force/Purple Heart
- Combo Coast Guard/Purple Heart \$149

### Seat Back Angle

Standard is 90° (perpendicular to floor). Minus (-°) closes back angle. Plus (+°) opens back angle.

- Standard  +3°  +6°  +9°  -3°  -6°  -9°  -12° N/C

### Back Upholstery Options

Available only with Height Adjustable Back Posts. If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Padded Velcro Adjustable with Standard Cover N/C
- Padded Velcro Adjustable with Air Mesh Breathable Cover \$77
- Tension Adjustable with Standard Cover (E2611) \$359
- Tension Adjustable with Air Mesh Breathable Cover (E2611) \$400
- Omit (\$10)

## Front Seat Height and Casters

### Caster Size / Footrest Taper / Front Seat Height

Select Front Seat Height and Footrest Taper in row corresponding to Caster Size. 6" casters not available with 85° and 90° Front Frame Bend.

Footrest Taper indicates bend of front frame creating footrest width. Taper is measured from outside of seat frame to outside of front tube. 2" Taper not available with 75° Front Frame Bend. Taper is not available with 70° Front Frame Bend. See reference E on page 7.

Front Seat Height is measured from floor to top of seat tube at front of seat sling. Adjustable up or down by 1/2". See reference F on page 7.

Footrest Taper	Price	Caster Size	Front Seat Height														
<input type="checkbox"/> Straight	N/C	3"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"
		4"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	
		5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"		
		6"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"			
<input type="checkbox"/> 1"	N/C	3"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"						
		4"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"							
		5"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"								
		6"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"									
<input type="checkbox"/> 2"	N/C	3"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"						
		4"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"							
		5"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"								
		6"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"									
<input type="checkbox"/> 3.5° V-Front	\$103	3"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"
		4"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	
		5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"		
		6"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"			
<input type="checkbox"/> 7° V-Front	\$103	3"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"
		4"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"	
		5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"		
		6"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"	<input type="checkbox"/> 20.5"			

### Casters

		Caster Size				
		Price	3"	4"	5"	6"
Caster Type	.75" Rollerblade	N/C	<input type="checkbox"/>	<input type="checkbox"/>		
	.75" Lighted Rollerblade	\$72	<input type="checkbox"/>	<input type="checkbox"/>		
	1" Poly	N/C		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1" Poly Aluminum	\$113		<input type="checkbox"/>	<input type="checkbox"/>	
	1" Pneumatic	\$87				<input type="checkbox"/>
	1.5" Poly	\$67		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1.5" Soft Roll Aluminum	\$210	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Caster Position

- See reference L on page 7*
- Optimized \_\_\_\_\_ N/C  
*Caster Placement minimizes the footprint of the chair*
- Expanded Width \_\_\_\_\_ N/C  
*Caster arm is positioned approximately .75" wider than Optimized position*

### Fork Options

- Standard Fork \_\_\_\_\_ N/C
- Single Sided Fork \_\_\_\_\_ \$190
- Frog Legs II Suspension System \_\_\_\_\_ \$461  
*Lightweight shock absorbing fork. Not included in Colored Anodizing Package.*

### Footrest

#### Footplate

- Tubular "Open" Footplate \_\_\_\_\_ N/C  
*Not available with V-Taper*
- Tubular Footplate with ABS Cover \_\_\_\_\_ \$31  
*Not available w/footrest widths greater than 15.5". See Reference G on page 7.  
 Not available with V-Taper.*
- Aluminum Angle Adjustable (K0040) \_\_\_\_\_ \$100
- Aluminum Angle Adjustable with Riser (K0040) \_\_\_\_\_ \$169
- High Mount Aluminum Angle Adjustable (K0040) \_\_\_\_\_ \$169  
*Not available with V-Taper*
- Angle Adjustable Flip Under (K0040, K0037) \_\_\_\_\_ \$231  
*Not available with footrest widths less than 9.5". Not available with V-Taper.  
 See Reference G on page 7.*
- High Mount Angle Adjustable Flip Under (K0040, K0037) \_\_\_\_\_ \$277  
*Not available w/footrest widths greater than 15.5". See Reference G on page 7. Not available with V-Taper.*
- Hybrid Angle Adjustable \_\_\_\_\_ \$92
- Hybrid Angle Adjustable with Wings \_\_\_\_\_ \$133

#### Footrest Protector

- Footrest Protector \_\_\_\_\_ N/C  
*Not available with Angle Adjustable Flip Under, High Mount Angle Adjustable Flip Under or High Mount Aluminum Angle Adjustable Footplates*

#### Seat to Footrest Length

*Measured from front edge of seat sling to top, rear of footrest. Footrest length at least 2.5" shorter than Front Seat Height recommended. See reference H on page 7.*

<input type="checkbox"/> 8"	<input type="checkbox"/> 8.5"	<input type="checkbox"/> 9"	<input type="checkbox"/> 9.5"	<input type="checkbox"/> 10"	<input type="checkbox"/> 10.5"
<input type="checkbox"/> 11"	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"
<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"
<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"

### Rear Seat Height and Wheels

#### Wheel Size / Rear Seat Height

*Select Rear Seat Height in row corresponding to desired Wheel Size. For proper setup, selection is required when omitting rear wheels. Measured from floor to top of seat tube at front of back post. Rear Seat Height is custom to needs. See reference I on page 7.*

Wheel	Rear Seat Height															
20"	<input type="checkbox"/> 12.5" <sup>1</sup>	<input type="checkbox"/> 13"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"
22"	<input type="checkbox"/> 13.5" <sup>1</sup>	<input type="checkbox"/> 14"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"		
24"	<input type="checkbox"/> 14.5" <sup>1</sup>	<input type="checkbox"/> 15"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"				
25"	<input type="checkbox"/> 15" <sup>1</sup>	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"					
26"	<input type="checkbox"/> 15.5" <sup>1</sup>	<input type="checkbox"/> 16"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19"	<input type="checkbox"/> 19.5"	<input type="checkbox"/> 20"						

<sup>1</sup> Available only with Low Profile tires

#### One Arm Drive

- One Arm Drive - Quick Release (E0958) \_\_\_\_\_ \$1532  
*Available only with 20", 22" and 24" spoke wheels and quick release axles. Spoke wheel style provided is exclusive to OAD. Standard handrim will be Aluminum Anodized. Limitations exist based on rear wheel spacing, camber and armrest selections.*
- Left Side Drive \_\_\_\_\_ N/C
- Right Side Drive \_\_\_\_\_ N/C
- Aluminum Drive Handrim \_\_\_\_\_ N/C
- Poly Coated Drive Handrim \_\_\_\_\_ \$205

#### Wheels

*No selection needed with One Arm Drive*

		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	MAXX Spoke	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Mag (Color) <sup>2</sup>	\$120	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Superlight Spoke	\$323			<input type="checkbox"/>	<input type="checkbox"/> <sup>1</sup>	
	Spinergy Spox	\$943		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Spinergy LX	\$1179		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<sup>1</sup> Available only with Superlight and Plastic Coated Handrims

<sup>2</sup> Color selection on page 6

#### Spinergy Options

- Select spoke color and handrim tab for Spinergy wheels*
- Black Spokes  Red Spokes  Blue Spokes  White Spokes \_\_\_\_\_ N/C

## Rear Tire Type

### Full Profile

- Street Pneumatic 1 3/8" (E2211, E2212) \_\_\_\_\_ N/C  
*Not available on 25" wheel size*
- Full Poly 1 3/8" (E2218) \_\_\_\_\_ N/C  
*Not available on Superlight wheel or 25" or 26" wheel sizes*
- Pneumatic with Airless Insert 1 3/8" (E2211, E2213) \_\_\_\_\_ \$97  
*Not available on Superlight wheel or 25" wheel size*
- Schwalbe HS-127 1 3/8" (E2211, E2212) \_\_\_\_\_ \$128  
*Available only on 24" wheel size. Not available on Mag wheel*
- Mountain 2" Black (E2211, E2212) \_\_\_\_\_ \$108  
*Available only on 22", 24" and 25" MAXX wheels and Spinergy wheels*

### Low Profile

- Low Poly 1" (E2218) \_\_\_\_\_ \$77  
*Not available on Superlight wheel, 25" or 26" wheel sizes*
- Sentinel Solid Tire 1" \_\_\_\_\_ \$97  
*Not available on Superlight wheel*
- High Pressure 1" (E2211, E2212) \_\_\_\_\_ \$97  
*Not available on Mag wheel*
- Schwalbe One 1" (E2211, E2212) \_\_\_\_\_ \$128  
*Available only on 24" and 25" wheel sizes. Not available on Mag wheel*
- Schwalbe Marathon Plus Evolution 1" (E2211, E2212) \_\_\_\_\_ \$144  
*Not available on Mag wheel, 20" or 22" wheel sizes*

## Handrims

*No selection needed with One Arm Drive*

### Standard Handrims

- Aluminum Anodized \_\_\_\_\_ N/C
- Superlight Aluminum Anodized \_\_\_\_\_ \$51  
*Available only on 24" and 25" wheel sizes. Available only on MAXX Spokes, Superlight or Spinergy wheels.*
- Plastic Coated \_\_\_\_\_ \$118
- Projection (8 Vertical Projections) \_\_\_\_\_ \$359  
*Available only with full profile tires in combination with 22", 24", and 26" MAXX Spoke wheels, or 22" and 24" Mag wheels*
- Omit Handrim \_\_\_\_\_ N/C

### Ergonomic Handrims

*Not available with 20" or 26" wheels. Tab mount only available with 24" wheels.*

- Ergo Lite Regular Grip \_\_\_\_\_ \$477
- Ergo Standard Regular Grip \_\_\_\_\_ \$456
- Ergo Lite Extra Grip \_\_\_\_\_ \$528
- Ergo Standard Extra Grip \_\_\_\_\_ \$507
- Ergo Lite Omit Grip \_\_\_\_\_ \$405
- Ergo Standard Omit Grip \_\_\_\_\_ \$384

## Handrim Options

### Handrim Mount

*If no selection is made, Standard Mount will be used*

- Close Mount \_\_\_\_\_ \$26

### Micro-Grip

- 3M Micro Grip \_\_\_\_\_ \$67  
*Not available on Plastic Coated or Projection Handrims*

## Axles

- Quick Release \_\_\_\_\_ N/C
- Quad Quick Release \_\_\_\_\_ \$190  
*Not available with Mag or Spinergy wheels*

## Center of Gravity Presets

*Measured from front of back post to center of rear axle. 0" indicates axle will be directly under back post. See reference J on page 7.*

Amputee Positions				
<input type="checkbox"/> -1.5"	<input type="checkbox"/> -1"	<input type="checkbox"/> -.5"	<input type="checkbox"/> 0"	<input type="checkbox"/> .25"
Standard Positions				
<input type="checkbox"/> .5"	<input type="checkbox"/> .75"	<input type="checkbox"/> 1"	<input type="checkbox"/> 1.25"	<input type="checkbox"/> 1.5"
<input type="checkbox"/> 1.75"	<input type="checkbox"/> 2"	<input type="checkbox"/> 2.25"	<input type="checkbox"/> 2.5"	<input type="checkbox"/> 2.75"
<input type="checkbox"/> 3"	<input type="checkbox"/> 3.25"	<input type="checkbox"/> 3.5"	<input type="checkbox"/> 3.75"	<input type="checkbox"/> 4"
<input type="checkbox"/> 4.25"	<input type="checkbox"/> 4.5"	<input type="checkbox"/> 4.75"	<input type="checkbox"/> 5"	

## Rear Wheel Spacing

*1" minimum with 0° camber. 1.5" minimum with T-Arms. The following apply if One Arm Drive is selected: 1.75" minimum with T-Arms. 1.25" minimum with select Rear Seat Heights. See Reference K on page 7.*

- 0.75"  1"  1.25"  1.5"  1.75" \_\_\_\_\_ N/C

## Rear Wheel Camber

*The following apply if One Arm Drive is selected: 6° and 8° not available. 4° not available with T-Arms.*

- 0°  2°  4°  6°  8° \_\_\_\_\_ N/C

## Wheel Locks

### Wheel Lock

- Push to Lock \_\_\_\_\_ N/C
- Push to Lock Flush Mount \_\_\_\_\_ N/C
- Pull to Lock \_\_\_\_\_ N/C
- Under Seat Scissor Lock \_\_\_\_\_ \$36
- Low Profile Scissor Lock \_\_\_\_\_ \$67  
*Not available with Low Poly, Full Poly or Sentinel Solid tires*
- Short Thro Scissor Lock \_\_\_\_\_ \$164
- Omit \_\_\_\_\_ N/C

### Wheel Lock Options

- Extension Handles (E0961) \_\_\_\_\_ \$86  
*Available only with Push to Lock or Pull to Lock Wheel Locks*
- Grade Aids (E0974) \_\_\_\_\_ \$154  
*Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires*

## Armrests

Armrest Type	Length	Price
Height Adjustable T-Arm Tall (11.5"-15") E0973	Desk <input type="checkbox"/>	\$318
	Full <input type="checkbox"/>	
Height Adjustable T-Arm (8.5"-13.5") E0973	Desk <input type="checkbox"/>	\$349
	Full <input type="checkbox"/>	
Height adjustable T-Arm Low (6.5"-10.5") E0973 <sup>1</sup>	Desk <input type="checkbox"/>	\$349
	Full <input type="checkbox"/>	

<sup>1</sup> If requested configuration is unavailable, standard height T-Arm will be used. 14" minimum seat depth with Offset Frame.

<sup>2</sup> Includes direct mount receiver

Armrest Type	Price
Swing Away (8.5"-12.5") E0973 <input type="checkbox"/>	\$244
Swing Away with Side Guard (8.5"-12.5") E0973 <sup>2</sup> <input type="checkbox"/>	\$292
Omit Armrests <input type="checkbox"/>	N/C

## Side Guards

### Direct Mount Side Guards

*Not available with T-Arms or Swing Away with Side Guard Armrests*

- Composite Side Guards - Direct Mount Receiver (K0108) \$149
- Aluminum Straight - Direct Mount Receiver - Black Matte (K0108) \$200
- Aluminum Straight - Direct Mount Receiver - Match Frame Color (K0108) \$251  
*Not available with Hydrographics*
- Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver - Black Matte (K0108) \$282  
*Not available with Adjustable Depth greater than 1/2"*
- Aluminum Adjustable Fender w/Foam Pad - Direct Mount Receiver - Match Frame Color (K0108) \$333  
*Not available with Adjustable Depth greater than 1/2". Not available with Hydrographics.*
- Carbon Fiber Adjustable Fender - Direct Mount Receiver (K0108) \$497  
*Not available with Adjustable Depth greater than 1/2"*

### Clamp-On Mount Side Guards

*Not available with T-Arms or Swing Away with Side Guard Armrests*

- Composite Side Guards - Clamp-On Receiver (K0108) \$179
- Aluminum Straight - Clamp On Receiver - Black Matte (K0108) \$231
- Aluminum Straight - Clamp-On Receiver - Match Frame Color (K0108) \$282  
*Not available with Hydrographics*
- Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Black Matte (K0108) \$313
- Aluminum Adjustable Fender w/Foam Pad - Clamp-On Receiver - Match Frame Color (K0108) \$364  
*Not available with Hydrographics*
- Carbon Fiber Adjustable Fender - Clamp-On Receiver (K0108) \$528

## Accessories

### Anti-Tips

- Rear Anti-Tips (E0971) \$169
- User Activated - Pair \$349
- User Activated - Left \$174
- User Activated - Right \$174
- Omit N/C

### Pelvic Positioning Belt

- 1 1/2" Auto Buckle Pelvic Positioning Belt (E0978) \$87
- 1 1/2" Padded Auto Buckle Pelvic Positioning Belt (E0978) \$174
- Small Bodypoint Evoflex (16-24") (E0978) \$193
- Medium Bodypoint Evoflex (19-29") (E0978) \$193
- Large Bodypoint Evoflex (25-37") (E0978) \$200

### Belt Mounting Kit

- Ki Mobility Belt Mounting Clamp \$36
- Bodypoint Belt Mounting Bracket Kit \$113

### Calf Strap

- Velcro Adjustable (K0038) \$43
- Padded Velcro Adjustable (K0038) \$70
- Bodypoint Padded AeroMesh (K0038) \$121

### Other Accessories

- Tool Kit \$51
- Removable Underseat Pouch \$51
- Ki Mobility Backpack \$67
- Luggage Carrier \$87
- Black Neoprene Impact Guards \$97
- Under Chair Carrier and Removable Bag \$87  
*Not available with Seat Widths less than 14" or scissor locks*
- Cane and Crutch Holder \$200

## Color

Options	Price	Colors																			
		Gloss												Matte							
		Black	Blue Skies	Bubblegum Pink	Candy Blue	Candy Purple	Candy Red	Copper	Granny Smith	MRGallic Orange	Silver	Teal	Yellow	Charcoal Gray	Emerald Green	Iridescent Burgundy	Iridescent Navy	Military Green	Retro Red	Royal Blue	Shadow Black
Frame	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Backrest Release Bar	\$51	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	N/C	<input type="checkbox"/> Black Anodized																			
MAXX Mag Color		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Colored Anodizing Package

Select desired color. If no selection is made, black will be used.

	Color	Black <sup>1</sup>	Blue	Bright Green	Gold	Orange	Pewter	Pink	Purple	Red
Retail		N/C	\$26	\$26	\$26	\$26	\$26	\$26	\$26	\$26
Option	Back Plates <sup>2</sup>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Camber Tube <sup>3</sup>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Towers	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Side Guard Receiver	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Caster Fork <sup>3</sup>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

<sup>1</sup> Black is no charge

<sup>2</sup> Not available with Transit Option

<sup>3</sup> Available only with Standard or Single Sided

**References**

**A. Seat Width** See Diagram 1

Measured from outside of frame tube on one side to outside of frame tube on other side.

**B. Seat Depth** See Diagram 3

Measured from front of back posts to front edge of seat sling. Seat sling starts at beginning of bend at front of frame.

**C. Sling Position** See Diagram 2

In Performance Position (0°), sling is at front of frame to keep chair shorter and more maneuverable. The 1° or 2° seat sling position extends the frame in front of the sling by that amount. Adding more frame can improve stability and provide support to aid in transfers.

**D. Front Frame Bend** See Diagram 3

Angle between front seat tube and ground.

**E. Footrest Taper** See Diagram 1

Indicates bend of front frame creating footrest. Taper is measured from outside of seat frame to outside of front tube.

**F. Front Seat Height** See Diagram 1

Measured from floor to top of seat tube at front of seat sling.

**G. Footrest Width** See Diagram 1

Measured from inside of front tube to inside of front tube on other side. Footrest width is listed below for each seat width and taper.

Seat Width	Straight	1° Y Taper	2° Y Taper	3.5° V Taper	7° V Taper
Outside Measurement	2" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width	4" Narrower than Seat Width	6" Narrower than Seat Width
12"	9.5"	7.5"		7.5"	
13"	10.5"	8.5"		8.5"	
14"	11.5"	9.5"		9.5"	
15"	12.5"	10.5"		10.5"	
16"	13.5"	11.5"	9.5"	11.5"	9.5"
17"	14.5"	12.5"	10.5"	12.5"	10.5"
18"	15.5"	13.5"	11.5"	13.5"	11.5"
19"	16.5"	14.5"	12.5"	14.5"	12.5"
20"	17.5"	15.5"	13.5"	15.5"	13.5"

**H. Seat to Footrest Length** See Diagram 3

Measured from front edge of seat sling to top rear of footrest. Footrest length of at least 2.5" shorter than front seat height recommended.

**I. Rear Seat Height** See Diagram 3

Measured from floor to top of seat tube at front of back post. Rear seat height is custom to needs.

**J. Center of Gravity Preset** See Diagram 3

Measured from front of back post to center of rear axle. 0° indicates axle will be directly under back post.

**K. Rear Wheel Spacing** See Diagram 1

Measured from outside of seat back to inside of rear tire. Adjustable out .5" from setting.

**L. Caster Position** See Diagram 4

The Optimized Position keeps the casters tighter to the frame for enhanced maneuverability. Expanded Position pushes the casters out .75" past the Optimized Position for increased stability.

Diagram 1




Diagram 2




Diagram 3





Diagram 4





Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom G - General Use Cushion

Weight Capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		\$138	\$138	\$138	\$138	\$138	\$138	\$138	\$229
HCPC		E2601	E2601	E2601	E2601	E2601	E2601	E2601	E2602
Depth	14"	<input type="checkbox"/> XG1414							
	15"	<input type="checkbox"/> XG1415	<input type="checkbox"/> XG1515	<input type="checkbox"/> XG1615					
	16"	<input type="checkbox"/> XG1416	<input type="checkbox"/> XG1516	<input type="checkbox"/> XG1616	<input type="checkbox"/> XG1716	<input type="checkbox"/> XG1816	<input type="checkbox"/> XG1916	<input type="checkbox"/> XG2016	
	17"	<input type="checkbox"/> XG1417	<input type="checkbox"/> XG1517	<input type="checkbox"/> XG1617	<input type="checkbox"/> XG1717	<input type="checkbox"/> XG1817	<input type="checkbox"/> XG1917	<input type="checkbox"/> XG2017	
	18"	<input type="checkbox"/> XG1418	<input type="checkbox"/> XG1518	<input type="checkbox"/> XG1618	<input type="checkbox"/> XG1718	<input type="checkbox"/> XG1818	<input type="checkbox"/> XG1918	<input type="checkbox"/> XG2018	<input type="checkbox"/> XG2218
	19"		<input type="checkbox"/> XG1519	<input type="checkbox"/> XG1619	<input type="checkbox"/> XG1719	<input type="checkbox"/> XG1819	<input type="checkbox"/> XG1919	<input type="checkbox"/> XG2019	<input type="checkbox"/> XG2219
	20"		<input type="checkbox"/> XG1520	<input type="checkbox"/> XG1620	<input type="checkbox"/> XG1720	<input type="checkbox"/> XG1820	<input type="checkbox"/> XG1920	<input type="checkbox"/> XG2020	<input type="checkbox"/> XG2220
	21"			<input type="checkbox"/> XG1621		<input type="checkbox"/> XG1821	<input type="checkbox"/> XG1921	<input type="checkbox"/> XG2021	<input type="checkbox"/> XG2221
22"			<input type="checkbox"/> XG1622		<input type="checkbox"/> XG1822	<input type="checkbox"/> XG1922	<input type="checkbox"/> XG2022	<input type="checkbox"/> XG2222	

## Axiom G - General Use Cushion with Extra Outer Cover

Weight Capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		\$224	\$224	\$224	\$224	\$224	\$224	\$224	\$315
HCPC		E2601	E2601	E2601	E2601	E2601	E2601	E2601	E2602
Depth	14"	<input type="checkbox"/> XG1414C2							
	15"	<input type="checkbox"/> XG1415C2	<input type="checkbox"/> XG1515C2	<input type="checkbox"/> XG1615C2					
	16"	<input type="checkbox"/> XG1416C2	<input type="checkbox"/> XG1516C2	<input type="checkbox"/> XG1616C2	<input type="checkbox"/> XG1716C2	<input type="checkbox"/> XG1816C2	<input type="checkbox"/> XG1916C2	<input type="checkbox"/> XG2016C2	
	17"	<input type="checkbox"/> XG1417C2	<input type="checkbox"/> XG1517C2	<input type="checkbox"/> XG1617C2	<input type="checkbox"/> XG1717C2	<input type="checkbox"/> XG1817C2	<input type="checkbox"/> XG1917C2	<input type="checkbox"/> XG2017C2	
	18"	<input type="checkbox"/> XG1418C2	<input type="checkbox"/> XG1518C2	<input type="checkbox"/> XG1618C2	<input type="checkbox"/> XG1718C2	<input type="checkbox"/> XG1818C2	<input type="checkbox"/> XG1918C2	<input type="checkbox"/> XG2018C2	<input type="checkbox"/> XG2218C2
	19"		<input type="checkbox"/> XG1519C2	<input type="checkbox"/> XG1619C2	<input type="checkbox"/> XG1719C2	<input type="checkbox"/> XG1819C2	<input type="checkbox"/> XG1919C2	<input type="checkbox"/> XG2019C2	<input type="checkbox"/> XG2219C2
	20"		<input type="checkbox"/> XG1520C2	<input type="checkbox"/> XG1620C2	<input type="checkbox"/> XG1720C2	<input type="checkbox"/> XG1820C2	<input type="checkbox"/> XG1920C2	<input type="checkbox"/> XG2020C2	<input type="checkbox"/> XG2220C2
	21"			<input type="checkbox"/> XG1621C2		<input type="checkbox"/> XG1821C2	<input type="checkbox"/> XG1921C2	<input type="checkbox"/> XG2021C2	<input type="checkbox"/> XG2221C2
22"			<input type="checkbox"/> XG1622C2		<input type="checkbox"/> XG1822C2	<input type="checkbox"/> XG1922C2	<input type="checkbox"/> XG2022C2	<input type="checkbox"/> XG2222C2	

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0236.10 Pricing and specifications subject to change without notice.





Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom S - Skin Protection Cushion

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		\$235	\$235	\$235	\$235	\$235	\$235	\$235	\$352
HCPC		E2603	E2603	E2603	E2603	E2603	E2603	E2603	E2604
Depth	14"	<input type="checkbox"/> XS1414							
	15"	<input type="checkbox"/> XS1415	<input type="checkbox"/> XS1515	<input type="checkbox"/> XS1615					
	16"	<input type="checkbox"/> XS1416	<input type="checkbox"/> XS1516	<input type="checkbox"/> XS1616	<input type="checkbox"/> XS1716	<input type="checkbox"/> XS1816	<input type="checkbox"/> XS1916	<input type="checkbox"/> XS2016	
	17"	<input type="checkbox"/> XS1417	<input type="checkbox"/> XS1517	<input type="checkbox"/> XS1617	<input type="checkbox"/> XS1717	<input type="checkbox"/> XS1817	<input type="checkbox"/> XS1917	<input type="checkbox"/> XS2017	
	18"	<input type="checkbox"/> XS1418	<input type="checkbox"/> XS1518	<input type="checkbox"/> XS1618	<input type="checkbox"/> XS1718	<input type="checkbox"/> XS1818	<input type="checkbox"/> XS1918	<input type="checkbox"/> XS2018	<input type="checkbox"/> XS2218
	19"		<input type="checkbox"/> XS1519	<input type="checkbox"/> XS1619	<input type="checkbox"/> XS1719	<input type="checkbox"/> XS1819	<input type="checkbox"/> XS1919	<input type="checkbox"/> XS2019	<input type="checkbox"/> XS2219
	20"		<input type="checkbox"/> XS1520	<input type="checkbox"/> XS1620	<input type="checkbox"/> XS1720	<input type="checkbox"/> XS1820	<input type="checkbox"/> XS1920	<input type="checkbox"/> XS2020	<input type="checkbox"/> XS2220
	21"			<input type="checkbox"/> XS1621		<input type="checkbox"/> XS1821	<input type="checkbox"/> XS1921	<input type="checkbox"/> XS2021	<input type="checkbox"/> XS2221
	22"			<input type="checkbox"/> XS1622		<input type="checkbox"/> XS1822	<input type="checkbox"/> XS1922	<input type="checkbox"/> XS2022	<input type="checkbox"/> XS2222

## Axiom S - Skin Protection Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		\$333	\$333	\$333	\$333	\$333	\$333	\$333	\$450
HCPC		E2603	E2603	E2603	E2603	E2603	E2603	E2603	E2604
Depth	14"	<input type="checkbox"/> XS1414C2							
	15"	<input type="checkbox"/> XS1415C2	<input type="checkbox"/> XS1515C2	<input type="checkbox"/> XS1615C2					
	16"	<input type="checkbox"/> XS1416C2	<input type="checkbox"/> XS1516C2	<input type="checkbox"/> XS1616C2	<input type="checkbox"/> XS1716C2	<input type="checkbox"/> XS1816C2	<input type="checkbox"/> XS1916C2	<input type="checkbox"/> XS2016C2	
	17"	<input type="checkbox"/> XS1417C2	<input type="checkbox"/> XS1517C2	<input type="checkbox"/> XS1617C2	<input type="checkbox"/> XS1717C2	<input type="checkbox"/> XS1817C2	<input type="checkbox"/> XS1917C2	<input type="checkbox"/> XS2017C2	
	18"	<input type="checkbox"/> XS1418C2	<input type="checkbox"/> XS1518C2	<input type="checkbox"/> XS1618C2	<input type="checkbox"/> XS1718C2	<input type="checkbox"/> XS1818C2	<input type="checkbox"/> XS1918C2	<input type="checkbox"/> XS2018C2	<input type="checkbox"/> XS2218C2
	19"		<input type="checkbox"/> XS1519C2	<input type="checkbox"/> XS1619C2	<input type="checkbox"/> XS1719C2	<input type="checkbox"/> XS1819C2	<input type="checkbox"/> XS1919C2	<input type="checkbox"/> XS2019C2	<input type="checkbox"/> XS2219C2
	20"		<input type="checkbox"/> XS1520C2	<input type="checkbox"/> XS1620C2	<input type="checkbox"/> XS1720C2	<input type="checkbox"/> XS1820C2	<input type="checkbox"/> XS1920C2	<input type="checkbox"/> XS2020C2	<input type="checkbox"/> XS2220C2
	21"			<input type="checkbox"/> XS1621C2		<input type="checkbox"/> XS1821C2	<input type="checkbox"/> XS1921C2	<input type="checkbox"/> XS2021C2	<input type="checkbox"/> XS2221C2
	22"			<input type="checkbox"/> XS1622C2		<input type="checkbox"/> XS1822C2	<input type="checkbox"/> XS1922C2	<input type="checkbox"/> XS2022C2	<input type="checkbox"/> XS2222C2

## Axiom P - Positioning Cushion

Weight Capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		\$407	\$407	\$407	\$407	\$407	\$407	\$407	\$528
HCPC		E2605	E2605	E2605	E2605	E2605	E2605	E2605	E2606
Depth	14"	<input type="checkbox"/> XP1414							
	15"	<input type="checkbox"/> XP1415	<input type="checkbox"/> XP1515	<input type="checkbox"/> XP1615					
	16"	<input type="checkbox"/> XP1416	<input type="checkbox"/> XP1516	<input type="checkbox"/> XP1616	<input type="checkbox"/> XP1716	<input type="checkbox"/> XP1816	<input type="checkbox"/> XP1916	<input type="checkbox"/> XP2016	
	17"	<input type="checkbox"/> XP1417	<input type="checkbox"/> XP1517	<input type="checkbox"/> XP1617	<input type="checkbox"/> XP1717	<input type="checkbox"/> XP1817	<input type="checkbox"/> XP1917	<input type="checkbox"/> XP2017	
	18"	<input type="checkbox"/> XP1418	<input type="checkbox"/> XP1518	<input type="checkbox"/> XP1618	<input type="checkbox"/> XP1718	<input type="checkbox"/> XP1818	<input type="checkbox"/> XP1918	<input type="checkbox"/> XP2018	<input type="checkbox"/> XP2218
	19"		<input type="checkbox"/> XP1519	<input type="checkbox"/> XP1619	<input type="checkbox"/> XP1719	<input type="checkbox"/> XP1819	<input type="checkbox"/> XP1919	<input type="checkbox"/> XP2019	<input type="checkbox"/> XP2219
	20"		<input type="checkbox"/> XP1520	<input type="checkbox"/> XP1620	<input type="checkbox"/> XP1720	<input type="checkbox"/> XP1820	<input type="checkbox"/> XP1920	<input type="checkbox"/> XP2020	<input type="checkbox"/> XP2220
	21"			<input type="checkbox"/> XP1621		<input type="checkbox"/> XP1821	<input type="checkbox"/> XP1921	<input type="checkbox"/> XP2021	<input type="checkbox"/> XP2221
	22"			<input type="checkbox"/> XP1622		<input type="checkbox"/> XP1822	<input type="checkbox"/> XP1922	<input type="checkbox"/> XP2022	<input type="checkbox"/> XP2222



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom P - Positioning Cushion with Extra Outer Cover

Weight Capacity: 350 lbs

		Width							
		14"	15"	16"	17"	18"	19"	20"	22"
Retail		\$517	\$517	\$517	\$517	\$517	\$517	\$517	\$638
HCPC		E2605	E2605	E2605	E2605	E2605	E2605	E2605	E2606
Depth	14"	<input type="checkbox"/> XP1414C2							
	15"	<input type="checkbox"/> XP1415C2	<input type="checkbox"/> XP1515C2	<input type="checkbox"/> XP1615C2					
	16"	<input type="checkbox"/> XP1416C2	<input type="checkbox"/> XP1516C2	<input type="checkbox"/> XP1616C2	<input type="checkbox"/> XP1716C2	<input type="checkbox"/> XP1816C2	<input type="checkbox"/> XP1916C2	<input type="checkbox"/> XP2016C2	
	17"	<input type="checkbox"/> XP1417C2	<input type="checkbox"/> XP1517C2	<input type="checkbox"/> XP1617C2	<input type="checkbox"/> XP1717C2	<input type="checkbox"/> XP1817C2	<input type="checkbox"/> XP1917C2	<input type="checkbox"/> XP2017C2	
	18"	<input type="checkbox"/> XP1418C2	<input type="checkbox"/> XP1518C2	<input type="checkbox"/> XP1618C2	<input type="checkbox"/> XP1718C2	<input type="checkbox"/> XP1818C2	<input type="checkbox"/> XP1918C2	<input type="checkbox"/> XP2018C2	<input type="checkbox"/> XP2218C2
	19"		<input type="checkbox"/> XP1519C2	<input type="checkbox"/> XP1619C2	<input type="checkbox"/> XP1719C2	<input type="checkbox"/> XP1819C2	<input type="checkbox"/> XP1919C2	<input type="checkbox"/> XP2019C2	<input type="checkbox"/> XP2219C2
	20"		<input type="checkbox"/> XP1520C2	<input type="checkbox"/> XP1620C2	<input type="checkbox"/> XP1720C2	<input type="checkbox"/> XP1820C2	<input type="checkbox"/> XP1920C2	<input type="checkbox"/> XP2020C2	<input type="checkbox"/> XP2220C2
	21"			<input type="checkbox"/> XP1621C2		<input type="checkbox"/> XP1821C2	<input type="checkbox"/> XP1921C2	<input type="checkbox"/> XP2021C2	<input type="checkbox"/> XP2221C2
	22"			<input type="checkbox"/> XP1622C2		<input type="checkbox"/> XP1822C2	<input type="checkbox"/> XP1922C2	<input type="checkbox"/> XP2022C2	<input type="checkbox"/> XP2222C2

## Axiom SP Visco - Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$493	\$493	\$493	\$493	\$493	\$493	\$493	\$603	\$603
HCPC		E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2608	E2608
Depth	14"	<input type="checkbox"/> XSP1414V								
	15"	<input type="checkbox"/> XSP1415V	<input type="checkbox"/> XSP1515V	<input type="checkbox"/> XSP1615V						
	16"	<input type="checkbox"/> XSP1416V	<input type="checkbox"/> XSP1516V	<input type="checkbox"/> XSP1616V	<input type="checkbox"/> XSP1716V	<input type="checkbox"/> XSP1816V	<input type="checkbox"/> XSP1916V	<input type="checkbox"/> XSP2016V		<input type="checkbox"/> XSP2216V
	17"	<input type="checkbox"/> XSP1417V	<input type="checkbox"/> XSP1517V	<input type="checkbox"/> XSP1617V	<input type="checkbox"/> XSP1717V	<input type="checkbox"/> XSP1817V	<input type="checkbox"/> XSP1917V	<input type="checkbox"/> XSP2017V		<input type="checkbox"/> XSP2217V
	18"	<input type="checkbox"/> XSP1418V	<input type="checkbox"/> XSP1518V	<input type="checkbox"/> XSP1618V	<input type="checkbox"/> XSP1718V	<input type="checkbox"/> XSP1818V	<input type="checkbox"/> XSP1918V	<input type="checkbox"/> XSP2018V	<input type="checkbox"/> XSP2118V	<input type="checkbox"/> XSP2218V
	19"		<input type="checkbox"/> XSP1519V	<input type="checkbox"/> XSP1619V	<input type="checkbox"/> XSP1719V	<input type="checkbox"/> XSP1819V	<input type="checkbox"/> XSP1919V	<input type="checkbox"/> XSP2019V	<input type="checkbox"/> XSP2119V	<input type="checkbox"/> XSP2219V
	20"		<input type="checkbox"/> XSP1520V	<input type="checkbox"/> XSP1620V	<input type="checkbox"/> XSP1720V	<input type="checkbox"/> XSP1820V	<input type="checkbox"/> XSP1920V	<input type="checkbox"/> XSP2020V	<input type="checkbox"/> XSP2120V	<input type="checkbox"/> XSP2220V
	21"			<input type="checkbox"/> XSP1621V		<input type="checkbox"/> XSP1821V	<input type="checkbox"/> XSP1921V	<input type="checkbox"/> XSP2021V	<input type="checkbox"/> XSP2121V	<input type="checkbox"/> XSP2221V
	22"			<input type="checkbox"/> XSP1622V		<input type="checkbox"/> XSP1822V	<input type="checkbox"/> XSP1922V	<input type="checkbox"/> XSP2022V	<input type="checkbox"/> XSP2122V	<input type="checkbox"/> XSP2222V

## Axiom SP Visco - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$614	\$614	\$614	\$614	\$614	\$614	\$614	\$724	\$724
HCPC		E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2608	E2608
Depth	14"	<input type="checkbox"/> XSP1414VC2								
	15"	<input type="checkbox"/> XSP1415VC2	<input type="checkbox"/> XSP1515VC2	<input type="checkbox"/> XSP1615VC2						
	16"	<input type="checkbox"/> XSP1416VC2	<input type="checkbox"/> XSP1516VC2	<input type="checkbox"/> XSP1616VC2	<input type="checkbox"/> XSP1716VC2	<input type="checkbox"/> XSP1816VC2	<input type="checkbox"/> XSP1916VC2	<input type="checkbox"/> XSP2016VC2		<input type="checkbox"/> XSP2216VC2
	17"	<input type="checkbox"/> XSP1417VC2	<input type="checkbox"/> XSP1517VC2	<input type="checkbox"/> XSP1617VC2	<input type="checkbox"/> XSP1717VC2	<input type="checkbox"/> XSP1817VC2	<input type="checkbox"/> XSP1917VC2	<input type="checkbox"/> XSP2017VC2		<input type="checkbox"/> XSP2217VC2
	18"	<input type="checkbox"/> XSP1418VC2	<input type="checkbox"/> XSP1518VC2	<input type="checkbox"/> XSP1618VC2	<input type="checkbox"/> XSP1718VC2	<input type="checkbox"/> XSP1818VC2	<input type="checkbox"/> XSP1918VC2	<input type="checkbox"/> XSP2018VC2	<input type="checkbox"/> XSP2118VC2	<input type="checkbox"/> XSP2218VC2
	19"		<input type="checkbox"/> XSP1519VC2	<input type="checkbox"/> XSP1619VC2	<input type="checkbox"/> XSP1719VC2	<input type="checkbox"/> XSP1819VC2	<input type="checkbox"/> XSP1919VC2	<input type="checkbox"/> XSP2019VC2	<input type="checkbox"/> XSP2119VC2	<input type="checkbox"/> XSP2219VC2
	20"		<input type="checkbox"/> XSP1520VC2	<input type="checkbox"/> XSP1620VC2	<input type="checkbox"/> XSP1720VC2	<input type="checkbox"/> XSP1820VC2	<input type="checkbox"/> XSP1920VC2	<input type="checkbox"/> XSP2020VC2	<input type="checkbox"/> XSP2120VC2	<input type="checkbox"/> XSP2220VC2
	21"			<input type="checkbox"/> XSP1621VC2		<input type="checkbox"/> XSP1821VC2	<input type="checkbox"/> XSP1921VC2	<input type="checkbox"/> XSP2021VC2	<input type="checkbox"/> XSP2121VC2	<input type="checkbox"/> XSP2221VC2
	22"			<input type="checkbox"/> XSP1622VC2		<input type="checkbox"/> XSP1822VC2	<input type="checkbox"/> XSP1922VC2	<input type="checkbox"/> XSP2022VC2	<input type="checkbox"/> XSP2122VC2	<input type="checkbox"/> XSP2222VC2



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom SP - Fluid

Weight Capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$504	\$504	\$504	\$504	\$504	\$504	\$504	\$617	\$617
HCPC		E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2608	E2608
Depth	14"	<input type="checkbox"/> XSP1414F								
	15"	<input type="checkbox"/> XSP1415F	<input type="checkbox"/> XSP1515F	<input type="checkbox"/> XSP1615F						
	16"	<input type="checkbox"/> XSP1416F	<input type="checkbox"/> XSP1516F	<input type="checkbox"/> XSP1616F	<input type="checkbox"/> XSP1716F	<input type="checkbox"/> XSP1816F	<input type="checkbox"/> XSP1916F	<input type="checkbox"/> XSP2016F		<input type="checkbox"/> XSP2216F
	17"	<input type="checkbox"/> XSP1417F	<input type="checkbox"/> XSP1517F	<input type="checkbox"/> XSP1617F	<input type="checkbox"/> XSP1717F	<input type="checkbox"/> XSP1817F	<input type="checkbox"/> XSP1917F	<input type="checkbox"/> XSP2017F		<input type="checkbox"/> XSP2217F
	18"	<input type="checkbox"/> XSP1418F	<input type="checkbox"/> XSP1518F	<input type="checkbox"/> XSP1618F	<input type="checkbox"/> XSP1718F	<input type="checkbox"/> XSP1818F	<input type="checkbox"/> XSP1918F	<input type="checkbox"/> XSP2018F	<input type="checkbox"/> XSP2118F	<input type="checkbox"/> XSP2218F
	19"		<input type="checkbox"/> XSP1519F	<input type="checkbox"/> XSP1619F	<input type="checkbox"/> XSP1719F	<input type="checkbox"/> XSP1819F	<input type="checkbox"/> XSP1919F	<input type="checkbox"/> XSP2019F	<input type="checkbox"/> XSP2119F	<input type="checkbox"/> XSP2219F
	20"		<input type="checkbox"/> XSP1520F	<input type="checkbox"/> XSP1620F	<input type="checkbox"/> XSP1720F	<input type="checkbox"/> XSP1820F	<input type="checkbox"/> XSP1920F	<input type="checkbox"/> XSP2020F	<input type="checkbox"/> XSP2120F	<input type="checkbox"/> XSP2220F
	21"			<input type="checkbox"/> XSP1621F		<input type="checkbox"/> XSP1821F	<input type="checkbox"/> XSP1921F	<input type="checkbox"/> XSP2021F	<input type="checkbox"/> XSP2121F	<input type="checkbox"/> XSP2221F
	22"			<input type="checkbox"/> XSP1622F		<input type="checkbox"/> XSP1822F	<input type="checkbox"/> XSP1922F	<input type="checkbox"/> XSP2022F	<input type="checkbox"/> XSP2122F	<input type="checkbox"/> XSP2222F

## Axiom SP Fluid - Skin Protection & Positioning Cushion with Extra Outer Cover

Weight Capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$625	\$625	\$625	\$625	\$625	\$625	\$625	\$738	\$738
HCPC		E2607	E2607	E2607	E2607	E2607	E2607	E2607	E2608	E2608
Depth	14"	<input type="checkbox"/> XSP1414FC2								
	15"	<input type="checkbox"/> XSP1415FC2	<input type="checkbox"/> XSP1515FC2	<input type="checkbox"/> XSP1615FC2						
	16"	<input type="checkbox"/> XSP1416FC2	<input type="checkbox"/> XSP1516FC2	<input type="checkbox"/> XSP1616FC2	<input type="checkbox"/> XSP1716FC2	<input type="checkbox"/> XSP1816FC2	<input type="checkbox"/> XSP1916FC2	<input type="checkbox"/> XSP2016FC2		<input type="checkbox"/> XSP2216FC2
	17"	<input type="checkbox"/> XSP1417FC2	<input type="checkbox"/> XSP1517FC2	<input type="checkbox"/> XSP1617FC2	<input type="checkbox"/> XSP1717FC2	<input type="checkbox"/> XSP1817FC2	<input type="checkbox"/> XSP1917FC2	<input type="checkbox"/> XSP2017FC2		<input type="checkbox"/> XSP2217FC2
	18"	<input type="checkbox"/> XSP1418FC2	<input type="checkbox"/> XSP1518FC2	<input type="checkbox"/> XSP1618FC2	<input type="checkbox"/> XSP1718FC2	<input type="checkbox"/> XSP1818FC2	<input type="checkbox"/> XSP1918FC2	<input type="checkbox"/> XSP2018FC2	<input type="checkbox"/> XSP2118FC2	<input type="checkbox"/> XSP2218FC2
	19"		<input type="checkbox"/> XSP1519FC2	<input type="checkbox"/> XSP1619FC2	<input type="checkbox"/> XSP1719FC2	<input type="checkbox"/> XSP1819FC2	<input type="checkbox"/> XSP1919FC2	<input type="checkbox"/> XSP2019FC2	<input type="checkbox"/> XSP2119FC2	<input type="checkbox"/> XSP2219FC2
	20"		<input type="checkbox"/> XSP1520FC2	<input type="checkbox"/> XSP1620FC2	<input type="checkbox"/> XSP1720FC2	<input type="checkbox"/> XSP1820FC2	<input type="checkbox"/> XSP1920FC2	<input type="checkbox"/> XSP2020FC2	<input type="checkbox"/> XSP2120FC2	<input type="checkbox"/> XSP2220FC2
	21"			<input type="checkbox"/> XSP1621FC2		<input type="checkbox"/> XSP1821FC2	<input type="checkbox"/> XSP1921FC2	<input type="checkbox"/> XSP2021FC2	<input type="checkbox"/> XSP2121FC2	<input type="checkbox"/> XSP2221FC2
	22"			<input type="checkbox"/> XSP1622FC2		<input type="checkbox"/> XSP1822FC2	<input type="checkbox"/> XSP1922FC2	<input type="checkbox"/> XSP2022FC2	<input type="checkbox"/> XSP2122FC2	<input type="checkbox"/> XSP2222FC2

## Axiom AS - Adjustable Skin Protection Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$641	\$641	\$641	\$641	\$641	\$641	\$641	\$743	\$743
HCPC		E2622	E2622	E2622	E2622	E2622	E2622	E2622	E2622	E2602
Height	14"	<input type="checkbox"/> XAS1414	<input type="checkbox"/> XAS1514	<input type="checkbox"/> XAS1614						
	15"	<input type="checkbox"/> XAS1415	<input type="checkbox"/> XAS1515	<input type="checkbox"/> XAS1615	<input type="checkbox"/> XAS1715					
	16"	<input type="checkbox"/> XAS1416	<input type="checkbox"/> XAS1516	<input type="checkbox"/> XAS1616	<input type="checkbox"/> XAS1716	<input type="checkbox"/> XAS1816	<input type="checkbox"/> XAS1916	<input type="checkbox"/> XAS2016		<input type="checkbox"/> XAS2216
	17"	<input type="checkbox"/> XAS1417	<input type="checkbox"/> XAS1517	<input type="checkbox"/> XAS1617	<input type="checkbox"/> XAS1717	<input type="checkbox"/> XAS1817	<input type="checkbox"/> XAS1917	<input type="checkbox"/> XAS2017		<input type="checkbox"/> XAS2217
	18"	<input type="checkbox"/> XAS1418	<input type="checkbox"/> XAS1518	<input type="checkbox"/> XAS1618	<input type="checkbox"/> XAS1718	<input type="checkbox"/> XAS1818	<input type="checkbox"/> XAS1918	<input type="checkbox"/> XAS2018	<input type="checkbox"/> XAS2118	<input type="checkbox"/> XAS2218
	19"	<input type="checkbox"/> XAS1419	<input type="checkbox"/> XAS1519	<input type="checkbox"/> XAS1619	<input type="checkbox"/> XAS1719	<input type="checkbox"/> XAS1819	<input type="checkbox"/> XAS1919	<input type="checkbox"/> XAS2019	<input type="checkbox"/> XAS2119	<input type="checkbox"/> XAS2219
	20"	<input type="checkbox"/> XAS1420	<input type="checkbox"/> XAS1520	<input type="checkbox"/> XAS1620	<input type="checkbox"/> XAS1720	<input type="checkbox"/> XAS1820	<input type="checkbox"/> XAS1920	<input type="checkbox"/> XAS2020	<input type="checkbox"/> XAS2120	<input type="checkbox"/> XAS2220
	21"		<input type="checkbox"/> XAS1521	<input type="checkbox"/> XAS1621	<input type="checkbox"/> XAS1721	<input type="checkbox"/> XAS1821	<input type="checkbox"/> XAS1921	<input type="checkbox"/> XAS2021	<input type="checkbox"/> XAS2121	<input type="checkbox"/> XAS2221
	22"		<input type="checkbox"/> XAS1522	<input type="checkbox"/> XAS1622	<input type="checkbox"/> XAS1722	<input type="checkbox"/> XAS1822	<input type="checkbox"/> XAS1922	<input type="checkbox"/> XAS2022	<input type="checkbox"/> XAS2122	<input type="checkbox"/> XAS2222



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Axiom AS - Adjustable Skin Protection Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$762	\$762	\$762	\$762	\$762	\$762	\$762	\$864	\$864
HCPC		E2622	E2622	E2622	E2622	E2622	E2622	E2622	E2622	E2602
Height	14"	<input type="checkbox"/> XAS1414C2	<input type="checkbox"/> XAS1514C2	<input type="checkbox"/> XAS1614C2						
	15"	<input type="checkbox"/> XAS1415C2	<input type="checkbox"/> XAS1515C2	<input type="checkbox"/> XAS1615C2	<input type="checkbox"/> XAS1715C2					
	16"	<input type="checkbox"/> XAS1416C2	<input type="checkbox"/> XAS1516C2	<input type="checkbox"/> XAS1616C2	<input type="checkbox"/> XAS1716C2	<input type="checkbox"/> XAS1816C2	<input type="checkbox"/> XAS1916C2	<input type="checkbox"/> XAS2016C2		<input type="checkbox"/> XAS2216C2
	17"	<input type="checkbox"/> XAS1417C2	<input type="checkbox"/> XAS1517C2	<input type="checkbox"/> XAS1617C2	<input type="checkbox"/> XAS1717C2	<input type="checkbox"/> XAS1817C2	<input type="checkbox"/> XAS1917C2	<input type="checkbox"/> XAS2017C2		<input type="checkbox"/> XAS2217C2
	18"	<input type="checkbox"/> XAS1418C2	<input type="checkbox"/> XAS1518C2	<input type="checkbox"/> XAS1618C2	<input type="checkbox"/> XAS1718C2	<input type="checkbox"/> XAS1818C2	<input type="checkbox"/> XAS1918C2	<input type="checkbox"/> XAS2018C2	<input type="checkbox"/> XAS2118C2	<input type="checkbox"/> XAS2218C2
	19"	<input type="checkbox"/> XAS1419C2	<input type="checkbox"/> XAS1519C2	<input type="checkbox"/> XAS1619C2	<input type="checkbox"/> XAS1719C2	<input type="checkbox"/> XAS1819C2	<input type="checkbox"/> XAS1919C2	<input type="checkbox"/> XAS2019C2	<input type="checkbox"/> XAS2119C2	<input type="checkbox"/> XAS2219C2
	20"	<input type="checkbox"/> XAS1420C2	<input type="checkbox"/> XAS1520C2	<input type="checkbox"/> XAS1620C2	<input type="checkbox"/> XAS1720C2	<input type="checkbox"/> XAS1820C2	<input type="checkbox"/> XAS1920C2	<input type="checkbox"/> XAS2020C2	<input type="checkbox"/> XAS2120C2	<input type="checkbox"/> XAS2220C2
	21"		<input type="checkbox"/> XAS1521C2	<input type="checkbox"/> XAS1621C2	<input type="checkbox"/> XAS1721C2	<input type="checkbox"/> XAS1821C2	<input type="checkbox"/> XAS1921C2	<input type="checkbox"/> XAS2021C2	<input type="checkbox"/> XAS2121C2	<input type="checkbox"/> XAS2221C2
	22"		<input type="checkbox"/> XAS1522C2	<input type="checkbox"/> XAS1622C2	<input type="checkbox"/> XAS1722C2	<input type="checkbox"/> XAS1822C2	<input type="checkbox"/> XAS1922C2	<input type="checkbox"/> XAS2022C2	<input type="checkbox"/> XAS2122C2	<input type="checkbox"/> XAS2222C2

## Axiom ASP - Adjustable Skin Protection & Positioning Cushion

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$661	\$661	\$661	\$661	\$661	\$661	\$661	\$764	\$764
HCPC		E2624	E2624	E2624	E2624	E2624	E2624	E2624	E2624	E2602
Height	14"	<input type="checkbox"/> XASP1414	<input type="checkbox"/> XASP1514	<input type="checkbox"/> XASP1614						
	15"	<input type="checkbox"/> XASP1415	<input type="checkbox"/> XASP1515	<input type="checkbox"/> XASP1615	<input type="checkbox"/> XASP1715	<input type="checkbox"/> XASP1815				
	16"	<input type="checkbox"/> XASP1416	<input type="checkbox"/> XASP1516	<input type="checkbox"/> XASP1616	<input type="checkbox"/> XASP1716	<input type="checkbox"/> XASP1816	<input type="checkbox"/> XASP1916	<input type="checkbox"/> XASP2016		<input type="checkbox"/> XASP2216
	17"	<input type="checkbox"/> XASP1417	<input type="checkbox"/> XASP1517	<input type="checkbox"/> XASP1617	<input type="checkbox"/> XASP1717	<input type="checkbox"/> XASP1817	<input type="checkbox"/> XASP1917	<input type="checkbox"/> XASP2017		<input type="checkbox"/> XASP2217
	18"	<input type="checkbox"/> XASP1418	<input type="checkbox"/> XASP1518	<input type="checkbox"/> XASP1618	<input type="checkbox"/> XASP1718	<input type="checkbox"/> XASP1818	<input type="checkbox"/> XASP1918	<input type="checkbox"/> XASP2018	<input type="checkbox"/> XASP2118	<input type="checkbox"/> XASP2218
	19"	<input type="checkbox"/> XASP1419	<input type="checkbox"/> XASP1519	<input type="checkbox"/> XASP1619	<input type="checkbox"/> XASP1719	<input type="checkbox"/> XASP1819	<input type="checkbox"/> XASP1919	<input type="checkbox"/> XASP2019	<input type="checkbox"/> XASP2119	<input type="checkbox"/> XASP2219
	20"	<input type="checkbox"/> XASP1420	<input type="checkbox"/> XASP1520	<input type="checkbox"/> XASP1620	<input type="checkbox"/> XASP1720	<input type="checkbox"/> XASP1820	<input type="checkbox"/> XASP1920	<input type="checkbox"/> XASP2020	<input type="checkbox"/> XASP2120	<input type="checkbox"/> XASP2220
	21"		<input type="checkbox"/> XASP1521	<input type="checkbox"/> XASP1621	<input type="checkbox"/> XASP1721	<input type="checkbox"/> XASP1821	<input type="checkbox"/> XASP1921	<input type="checkbox"/> XASP2021	<input type="checkbox"/> XASP2121	<input type="checkbox"/> XASP2221
	22"		<input type="checkbox"/> XASP1522	<input type="checkbox"/> XASP1622	<input type="checkbox"/> XASP1722	<input type="checkbox"/> XASP1822	<input type="checkbox"/> XASP1922	<input type="checkbox"/> XASP2022	<input type="checkbox"/> XASP2122	<input type="checkbox"/> XASP2222

## Axiom ASP - Adjustable Skin Protection & Positioning Cushion with Extra Outer Cover

Weight capacity: 350 lbs

		Width								
		14"	15"	16"	17"	18"	19"	20"	21"	22"
Retail		\$782	\$782	\$782	\$782	\$782	\$782	\$782	\$885	\$885
HCPC		E2624	E2624	E2624	E2624	E2624	E2624	E2624	E2624	E2602
Height	14"	<input type="checkbox"/> XASP1414C2	<input type="checkbox"/> XASP1514C2	<input type="checkbox"/> XASP1614C2						
	15"	<input type="checkbox"/> XASP1415C2	<input type="checkbox"/> XASP1515C2	<input type="checkbox"/> XASP1615C2	<input type="checkbox"/> XASP1715C2	<input type="checkbox"/> XASP1815C2				
	16"	<input type="checkbox"/> XASP1416C2	<input type="checkbox"/> XASP1516C2	<input type="checkbox"/> XASP1616C2	<input type="checkbox"/> XASP1716C2	<input type="checkbox"/> XASP1816C2	<input type="checkbox"/> XASP1916C2	<input type="checkbox"/> XASP2016C2		<input type="checkbox"/> XASP2216C2
	17"	<input type="checkbox"/> XASP1417C2	<input type="checkbox"/> XASP1517C2	<input type="checkbox"/> XASP1617C2	<input type="checkbox"/> XASP1717C2	<input type="checkbox"/> XASP1817C2	<input type="checkbox"/> XASP1917C2	<input type="checkbox"/> XASP2017C2		<input type="checkbox"/> XASP2217C2
	18"	<input type="checkbox"/> XASP1418C2	<input type="checkbox"/> XASP1518C2	<input type="checkbox"/> XASP1618C2	<input type="checkbox"/> XASP1718C2	<input type="checkbox"/> XASP1818C2	<input type="checkbox"/> XASP1918C2	<input type="checkbox"/> XASP2018C2	<input type="checkbox"/> XASP2118C2	<input type="checkbox"/> XASP2218C2
	19"	<input type="checkbox"/> XASP1419C2	<input type="checkbox"/> XASP1519C2	<input type="checkbox"/> XASP1619C2	<input type="checkbox"/> XASP1719C2	<input type="checkbox"/> XASP1819C2	<input type="checkbox"/> XASP1919C2	<input type="checkbox"/> XASP2019C2	<input type="checkbox"/> XASP2119C2	<input type="checkbox"/> XASP2219C2
	20"	<input type="checkbox"/> XASP1420C2	<input type="checkbox"/> XASP1520C2	<input type="checkbox"/> XASP1620C2	<input type="checkbox"/> XASP1720C2	<input type="checkbox"/> XASP1820C2	<input type="checkbox"/> XASP1920C2	<input type="checkbox"/> XASP2020C2	<input type="checkbox"/> XASP2120C2	<input type="checkbox"/> XASP2220C2
	21"		<input type="checkbox"/> XASP1521C2	<input type="checkbox"/> XASP1621C2	<input type="checkbox"/> XASP1721C2	<input type="checkbox"/> XASP1821C2	<input type="checkbox"/> XASP1921C2	<input type="checkbox"/> XASP2021C2	<input type="checkbox"/> XASP2121C2	<input type="checkbox"/> XASP2221C2
	22"		<input type="checkbox"/> XASP1522C2	<input type="checkbox"/> XASP1622C2	<input type="checkbox"/> XASP1722C2	<input type="checkbox"/> XASP1822C2	<input type="checkbox"/> XASP1922C2	<input type="checkbox"/> XASP2022C2	<input type="checkbox"/> XASP2122C2	<input type="checkbox"/> XASP2222C2



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



### AS & ASP Modifications

*Extended lead times will apply. Makes cushion non-refundable*

#### Rail Cut

XASNR

Bilateral Angled Seat Rail Relief Cut N/C

#### Notch for Back Posts

XASNB

2" x 2" Bilateral Rear Notch for Back Post Clearance N/C

## Axiom Air LTE - General Use Air Cushion

*Weight capacity: 315 lbs*

		Width			
		16"		18"	
		Cushion w/ Standard Cover	Cushion w/ Extra Standard Cover	Cushion w/ Standard Cover	Cushion w/ Extra Standard Cover
	Retail	\$160	\$244	\$160	\$244
	HCPC	E2601	E2601	E2601	E2601
Depth	16"	<input type="checkbox"/> XLTE1616	<input type="checkbox"/> XLTE1616C2	<input type="checkbox"/> XLTE1816	<input type="checkbox"/> XLTE1816C2
	18"	<input type="checkbox"/> XLTE1618	<input type="checkbox"/> XLTE1618C2	<input type="checkbox"/> XLTE1818	<input type="checkbox"/> XLTE1818C2

## Seating Accessories

Axiom Solid Seat Insert \$70



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Star StarLock Cushions

### Star StarLock Cushion (3" Tall Cells) and Extra Standard / Incontinent Cover

		Width														
		13.5"			15.5"			17.25"			19.25"			21"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent		
Retail	\$589	\$71	\$82	\$589	\$71	\$82	\$589	\$71	\$82	\$589	\$71	\$82	\$589	\$71	\$82	
HCPCS	E2624			E2624			E2624			E2624			E2624			
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	21"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

### Star StarLock Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover

		Width														
		13.5"			15.5"			17.25"			19.25"			21"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent		
Retail	\$589	\$71	\$82	\$589	\$71	\$82	\$589	\$71	\$82	\$589	\$71	\$82	\$589	\$71	\$82	
HCPCS	E2624			E2624			E2624			E2624			E2624			
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	17.25"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	21"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

### Star StarLock Cushion (5" Tall Cells) and Extra Standard / Incontinent Cover

		Width														
		14"			16"			18"			20"			22"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent		
Retail	\$604	\$71	\$82	\$604	\$71	\$82	\$604	\$71	\$82	\$604	\$71	\$82	\$704	\$93	\$106	
HCPCS	E2624			E2624			E2624			E2624			E2625			
Depth	14"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	18"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	20"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	22"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0903.2 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Star Stabil-Air Cushions

### Star Stabil-Air Cushion (3" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$563	\$71	\$82	\$563	\$71	\$82	\$563	\$71	\$82	\$563	\$71	\$82	\$563	\$71	\$82	\$692	\$93	\$106
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>															
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	21"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Star Stabil-Air Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$563	\$71	\$82	\$563	\$71	\$82	\$563	\$71	\$82	\$563	\$71	\$82	\$563	\$71	\$82	\$692	\$93	\$106
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>															
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	21"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Star Standard Air Cushions

### Star Standard Air Cushion (2" Tall Cells) and Extra Standard / Incontinent Cover

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent
Retail		\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$595	\$93	\$106
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>															
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	21"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



**Star Standard Air Cushion (3" Tall Cells) and Extra Standard / Incontinent Cover**

		Width														
		13.5"			15.5"			17.25"			19.25"			21"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	
Retail		\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82
HCPCS		E2622			E2622			E2622			E2622			E2622		
Depth	13.5"															
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	21"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**Star Standard Air Cushion (4" Tall Cells) and Extra Standard / Incontinent Cover**

		Width																	
		13.5"			15.5"			17.25"			19.25"			21"			22.75"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent		
Retail		\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$595	\$93	\$106
HCPCS		E2622			E2622			E2622			E2622			E2622			E2623		
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	19.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	21"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	

**Star Standard Air Cushions Dual Valve**


**Star Standard Air Cushion Dual Valve (4" Tall Cells) and Extra Standard / Incontinent Cover**

		Width														
		13.5"			15.5"			17.25"			19.25"			21"		
		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover		Cushion	Cover	
		Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	Standard	Incontinent	
Retail		\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82	\$487	\$71	\$82
HCPCS		E2622			E2622			E2622			E2622			E2622		
Depth	13.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>												
	15.5"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	17.25"				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	19.25"							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
	21"										<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12" -14" (30.5cm -35.5cm)
14" (35.5cm)	12" -13" (30.5cm -33.0cm)	14" -16" (35.5cm -40.5cm)
16" (40.5cm)	14" -15" (35.5cm -38.0cm)	16" -18" (40.5cm -45.5cm)
18" (45.5cm)	16" - 17" (40.5cm - 43.0cm)	18" - 20" (45.5cm - 50.5cm)
20" (51.0cm)	18" -19" (46.0cm -48.5cm)	20" -22" (50.5cm -55.5cm)
22" (56.0cm)	20" -21" (51.0cm -53.5cm)	22" -24" (56.0cm -61.0cm)

Hardware Description	Code
Improved DEX Hardware, removable adjustable hardware (300 lb max user weight)	<b>DEX</b>
New BTX Hardware, removable adjustable hardware not recommended for, or installed on Ultralight Rigid Frame Wheelchairs (400 lb max user weight)	<b>BTX</b>
New Fixed Hardware, non-removeable hardware up to (300 lb max user weight)	<b>FXT</b>



## Axiom Back

Hardware is included with back. Low Back selections are not available with headrests, modifications, or posture supports. Weight capacity: 300 lbs.

### Axiom Posterior Back (PB)

Includes Lumbar Pad

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		\$625	\$625	\$625	\$625	\$625	\$941	
HCPC		E2613	E2613	E2613	E2613	E2613	E2614	
Height	Low	8"	<input type="checkbox"/> AXPB1208	<input type="checkbox"/> AXPB1408	<input type="checkbox"/> AXPB1608	<input type="checkbox"/> AXPB1808	<input type="checkbox"/> AXPB2008	
		10"		<input type="checkbox"/> AXPB1410	<input type="checkbox"/> AXPB1610	<input type="checkbox"/> AXPB1810	<input type="checkbox"/> AXPB2010	<input type="checkbox"/> AXPB2210
		12"	<input type="checkbox"/> AXPB1212	<input type="checkbox"/> AXPB1412	<input type="checkbox"/> AXPB1612	<input type="checkbox"/> AXPB1812	<input type="checkbox"/> AXPB2012	
	Mid	14"		<input type="checkbox"/> AXPB1414	<input type="checkbox"/> AXPB1614	<input type="checkbox"/> AXPB1814	<input type="checkbox"/> AXPB2014	<input type="checkbox"/> AXPB2214
		16"	<input type="checkbox"/> AXPB1216	<input type="checkbox"/> AXPB1416	<input type="checkbox"/> AXPB1616	<input type="checkbox"/> AXPB1816	<input type="checkbox"/> AXPB2016	
		18"		<input type="checkbox"/> AXPB1418	<input type="checkbox"/> AXPB1618	<input type="checkbox"/> AXPB1818	<input type="checkbox"/> AXPB2018	<input type="checkbox"/> AXPB2218

### Back Hardware

DEX Hardware  BTX Hardware  FXT Hardware \_\_\_\_\_ N/C

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN0473.10 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



### Lateral Support for AXPB Backs

Pair	Description	Item Code	HCPC	Pair
<input type="checkbox"/>	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PBLF3x4PP	E0956	\$313
<input type="checkbox"/>	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PBLF4x5PP	E0956	\$313
<input type="checkbox"/>	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLS3X4PP	E0956, E1028	\$666
<input type="checkbox"/>	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLS4X5PP	E0956, E1028	\$666
<input type="checkbox"/>	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLAS3X4PP	E0956, E1028	\$810
<input type="checkbox"/>	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PBLAS3X4PP	E0956, E1028	\$810

### Axiom Posterior Lateral Back (PL)

*Includes Lumbar Pad*

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		\$722	\$722	\$722	\$722	\$722	\$1002	
HCPC		E2615	E2615	E2615	E2615	E2615	E2616	
Height	Low	8"	<input type="checkbox"/> AXPL1208	<input type="checkbox"/> AXPL1408	<input type="checkbox"/> AXPL1608	<input type="checkbox"/> AXPL1808	<input type="checkbox"/> AXPL2008	
		10"		<input type="checkbox"/> AXPL1410	<input type="checkbox"/> AXPL1610	<input type="checkbox"/> AXPL1810	<input type="checkbox"/> AXPL2010	<input type="checkbox"/> AXPL2210
		12"	<input type="checkbox"/> AXPL1212	<input type="checkbox"/> AXPL1412	<input type="checkbox"/> AXPL1612	<input type="checkbox"/> AXPL1812	<input type="checkbox"/> AXPL2012	
	Mid	14"		<input type="checkbox"/> AXPL1414	<input type="checkbox"/> AXPL1614	<input type="checkbox"/> AXPL1814	<input type="checkbox"/> AXPL2014	<input type="checkbox"/> AXPL2214
		16"	<input type="checkbox"/> AXPL1216	<input type="checkbox"/> AXPL1416	<input type="checkbox"/> AXPL1616	<input type="checkbox"/> AXPL1816	<input type="checkbox"/> AXPL2016	
		18"		<input type="checkbox"/> AXPL1418	<input type="checkbox"/> AXPL1618	<input type="checkbox"/> AXPL1818	<input type="checkbox"/> AXPL2018	<input type="checkbox"/> AXPL2218
	Tail	20"	<input type="checkbox"/> AXPL1220	<input type="checkbox"/> AXPL1420	<input type="checkbox"/> AXPL1620	<input type="checkbox"/> AXPL1820	<input type="checkbox"/> AXPL2020	
		22"		<input type="checkbox"/> AXPL1422	<input type="checkbox"/> AXPL1622	<input type="checkbox"/> AXPL1822	<input type="checkbox"/> AXPL2022	<input type="checkbox"/> AXPL2222

### Back Hardware

DEX Hardware  BTX Hardware  FXT Hardware \_\_\_\_\_ N/C

### Axiom Lateral Back (LB)

*Does not include Lumbar Pad*

		Width						
		12"	14"	16"	18"	20"	22"	
Retail		\$872	\$872	\$872	\$872	\$872	\$1014	
HCPC		E2620	E2620	E2620	E2620	E2620	E2621	
Height	Low	8"	<input type="checkbox"/> AXLB1208	<input type="checkbox"/> AXLB1408	<input type="checkbox"/> AXLB1608	<input type="checkbox"/> AXLB1808	<input type="checkbox"/> AXLB2008	
		10"		<input type="checkbox"/> AXLB1410	<input type="checkbox"/> AXLB1610	<input type="checkbox"/> AXLB1810	<input type="checkbox"/> AXLB2010	<input type="checkbox"/> AXLB2210
		12"	<input type="checkbox"/> AXLB1212	<input type="checkbox"/> AXLB1412	<input type="checkbox"/> AXLB1612	<input type="checkbox"/> AXLB1812	<input type="checkbox"/> AXLB2012	
	Mid	14"		<input type="checkbox"/> AXLB1414	<input type="checkbox"/> AXLB1614	<input type="checkbox"/> AXLB1814	<input type="checkbox"/> AXLB2014	<input type="checkbox"/> AXLB2214
		16"	<input type="checkbox"/> AXLB1216	<input type="checkbox"/> AXLB1416	<input type="checkbox"/> AXLB1616	<input type="checkbox"/> AXLB1816	<input type="checkbox"/> AXLB2016	
		18"		<input type="checkbox"/> AXLB1418	<input type="checkbox"/> AXLB1618	<input type="checkbox"/> AXLB1818	<input type="checkbox"/> AXLB2018	<input type="checkbox"/> AXLB2218
	Tail	20"	<input type="checkbox"/> AXLB1220	<input type="checkbox"/> AXLB1420	<input type="checkbox"/> AXLB1620	<input type="checkbox"/> AXLB1820	<input type="checkbox"/> AXLB2020	
		22"		<input type="checkbox"/> AXLB1422	<input type="checkbox"/> AXLB1622	<input type="checkbox"/> AXLB1822	<input type="checkbox"/> AXLB2022	<input type="checkbox"/> AXLB2222

### Back Hardware

DEX Hardware  BTX Hardware  FXT Hardware \_\_\_\_\_ N/C



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



### Lateral Support for AXPL & AXLB Backs

Pair	Description	Item Code	HCPC	Pair
<input type="checkbox"/>	3 x 4 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PLLF3x4PP	E0956	\$313
<input type="checkbox"/>	4 x 5 Lateral Pads with Fixed Mounting Hardware - Pair (E0956)	PLLF4x5PP	E0956	\$313
<input type="checkbox"/>	3 x 4 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLS3x4PP	E0956, E1028	\$666
<input type="checkbox"/>	4 x 5 Lateral Pads with Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLS4x5PP	E0956, E1028	\$666
<input type="checkbox"/>	3 x 4 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLAS3x4PP	E0956, E1028	\$810
<input type="checkbox"/>	4 x 5 Lateral Pads with Adjustable Swing Away Mounting Hardware - Pair (E0956,E1028)	PLLAS4x5PP	E0956, E1028	\$810

### Axiom Deep Lateral Back (DL)

*Does not include Lumbar Pad*

		Width						
		12"	14"	16"	18"	20"	22"	
	Retail	\$872	\$872	\$872	\$872	\$872	\$1014	
Height	Mid	14"	<input type="checkbox"/> AXDL1214	<input type="checkbox"/> AXDL1414	<input type="checkbox"/> AXDL1614	<input type="checkbox"/> AXDL1814	<input type="checkbox"/> AXDL2014	<input type="checkbox"/> AXDL2214
		16"		<input type="checkbox"/> AXDL1416	<input type="checkbox"/> AXDL1616	<input type="checkbox"/> AXDL1816	<input type="checkbox"/> AXDL2016	
		18"	<input type="checkbox"/> AXDL1218	<input type="checkbox"/> AXDL1418	<input type="checkbox"/> AXDL1618	<input type="checkbox"/> AXDL1818	<input type="checkbox"/> AXDL2018	<input type="checkbox"/> AXDL2218
	Tall	20"	<input type="checkbox"/> AXDL1220	<input type="checkbox"/> AXDL1420	<input type="checkbox"/> AXDL1620	<input type="checkbox"/> AXDL1820	<input type="checkbox"/> AXDL2020	
		22"		<input type="checkbox"/> AXDL1422	<input type="checkbox"/> AXDL1622	<input type="checkbox"/> AXDL1822	<input type="checkbox"/> AXDL2022	<input type="checkbox"/> AXDL2222

### Back Hardware

DEX Hardware  BTX Hardware  FXT Hardware \_\_\_\_\_ N/C

### Headrest

Description	Std. Straight Removable Adjustable Hardware		Multi-Axis Removable Adjustable Hardware		Posalinc Removable Adjustable Hardware		HCPCS
	Item Code	Retail	Item Code	Retail	Item Code	Retail	
Comfort Soft Headrest 8"	<input type="checkbox"/> CS08SS	\$347	<input type="checkbox"/> CS08MA	\$369	<input type="checkbox"/> CS08PL	\$519	E0955 / E1028
Comfort Soft Headrest 10"	<input type="checkbox"/> CS10SS	\$347	<input type="checkbox"/> CS10MA	\$369	<input type="checkbox"/> CS10PL	\$519	E0955 / E1028
Comfort Soft Headrest 14"	<input type="checkbox"/> CS14SS	\$347	<input type="checkbox"/> CS14MA	\$369	<input type="checkbox"/> CS14PL	\$519	E0955 / E1028
Contour Positioning Pad, 8" X 6.25"	<input type="checkbox"/> CP08SS	\$347	<input type="checkbox"/> CP08MA	\$369	<input type="checkbox"/> CP08PL	\$519	E0955 / E1028
Contour Positioning Pad, 9" X 6.5"	<input type="checkbox"/> CP09SS	\$347	<input type="checkbox"/> CP09MA	\$369	<input type="checkbox"/> CP09PL	\$519	E0955 / E1028
Contour Positioning Pad, 10" X 7.5"	<input type="checkbox"/> CP10SS	\$347	<input type="checkbox"/> CP10MA	\$369	<input type="checkbox"/> CP10PL	\$519	E0955 / E1028
Tri-soft Headrest 10"	<input type="checkbox"/> TS10SS	\$347	<input type="checkbox"/> TS10MA	\$369	<input type="checkbox"/> TS10PL	\$519	E0955 / E1028
Tri-soft Headrest 14"	<input type="checkbox"/> TS14SS	\$347	<input type="checkbox"/> TS14MA	\$369	<input type="checkbox"/> TS14PL	\$519	E0955 / E1028

Description	Item Code	Retail	HCPCS
Universal Headrest Mounting Plate	<input type="checkbox"/> UMP	\$51	K0108

*Not required for Axiom Headrests*



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Anterior Thoracic Supports

*Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.*

Anterior Thoracic Supports		
Description	Item Code	Retail
Dynaform Chest Support: 2-Piece with Hook&Loop, 3" High Pad, Small	<input type="checkbox"/> 14941	\$125
Dynaform Chest Support: 2-Piece with Hook&Loop, 4" High Pad, Medium	<input type="checkbox"/> 14942	\$125
Dynaform Chest Support: 2-Piece with Hook&Loop, 5" High Pad, Large	<input type="checkbox"/> 14943	\$125
Dynaform Chest Support: 1-Piece with Buckle, 3" Height x 20" Width, Small	<input type="checkbox"/> 14924	\$110
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 20" Width, Medium	<input type="checkbox"/> 14927	\$110
Dynaform Chest Support: 1-Piece with Buckle, 4" Height x 22" Width, Large	<input type="checkbox"/> 14928	\$110
Dynaform Anterior Trunk Support: Standard, Dynamic, X-Small	<input type="checkbox"/> 14851	\$130
Dynaform Anterior Trunk Support: Standard, Dynamic, Small	<input type="checkbox"/> 14852	\$130
Dynaform Anterior Trunk Support: Standard, Dynamic, Medium	<input type="checkbox"/> 14853	\$130
Dynaform Anterior Trunk Support: Standard, Dynamic, Large	<input type="checkbox"/> 14854	\$130
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, X-Small	<input type="checkbox"/> 14856	\$130
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Small	<input type="checkbox"/> 14857	\$130
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Medium	<input type="checkbox"/> 14858	\$130
Dynaform Anterior Trunk Support: Slim Cut, Dynamic, Large	<input type="checkbox"/> 14859	\$130
AirLogic Anterior Trunk Support: Standard, Structure, Small	<input type="checkbox"/> 14282	\$130
AirLogic Anterior Trunk Support: Standard, Structure, Medium	<input type="checkbox"/> 14283	\$175
AirLogic Anterior Trunk Support: Standard, Structure, Large	<input type="checkbox"/> 14284	\$175
AirLogic Anterior Trunk Support: Standard, Structure, X-Large	<input type="checkbox"/> 14285	\$175
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Small	<input type="checkbox"/> 14257	\$175
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Medium	<input type="checkbox"/> 14258	\$175
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, Large	<input type="checkbox"/> 14259	\$175
AirLogic Anterior Trunk Support: Slim Cut, Dynamic, X-Large	<input type="checkbox"/> 14260	\$175

## Accessories

- Adjustable Shoulder Guides (K0108) \$97  
*Includes Anterior Thoracic Mounting Modification. Not available with Low Backs.*
- Extra Lumbar Pad (K0108) \$51  
*Available only with Axiom Lateral Back and Axiom Deep Lateral back*
- Privacy Flap (K0108) \$51



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com



## Axiom Back Modifications

### Axiom Back Modifications

*Extended lead times will apply for all modifications. Makes back non-refundable.*

- Anterior Thoracic Mounting Modification N/C  
*Provides holes for Anterior Thoracic Supports and Shoulder Guides*

#### Gloss Finish

*If no selection is made, standard Shadow Black will be used with no additional lead time*

- Black Back Shell N/C
- Blue Skies Back Shell N/C
- Bubblegum Pink Back Shell N/C
- Candy Blue Back Shell N/C
- Candy Purple Back Shell N/C
- Candy Red Back Shell N/C
- Copper Back Shell N/C
- Granny Smith Back Shell N/C
- Metallic Orange Back Shell N/C
- Yellow Back Shell N/C
- Silver Back Shell N/C
- Teal Back Shell N/C

#### Matte Finish

*If no selection is made, standard Shadow Black will be used with no additional lead time*

- Military Green Back Shell N/C
- Charcoal Gray Back Shell N/C
- Emerald Green Back Shell N/C
- Iridescent Burgundy Back Shell N/C
- Iridescent Navy Back Shell N/C
- Retro Red Back Shell N/C
- Royal Blue Back Shell N/C
- Smooth White Back Shell N/C