

permobil

ROHO PRODIGY Mattress Overlay System

A non-powered, air filled, adjustable, zoned, reactive, full support surface; PRODIGY is intended to conform to an individual's shape to provide skin/soft tissue and deep tissue injury protection. Flexible, interconnected DRY FLOATATION® air cells divided into three sections adjust to match the individual's body contour by adjusting the three sections individually to address specific, changing needs of the user.

Designed for individuals who are at risk for breakdown and/or have Stage 1 or Stage 2 ischemic ulcers (pressure sores).

PRODUCT FEATURES

- Promotes a comfortable environment to assist with pain management.
- Accommodates other positioning devices beneath the overlay.
- Hydrostatic distribution of load across the entire contact area.
- Provides 3" (7.5 cm) of immersion.
- Durable polyvinyl material can be easily cleaned and may be patched for extended use.

Size of Complete Overlay: 81.75" Length x 36" Wide x 3.5" High

Weight Limit of User: 300 pounds (136 Kg)

Warranty: 12 month limited warranty



*Contact your medical equipment
supplier for more information or
to purchase.*

U.S. Medicare Code:
E0197 Group 1 Support Surfaces



*Prevention Through
Innovation.*

DISTRIBUTED BY:

5 MINDS MOBILITY GROUP

ggarland@5mindsmobility.com

orders@5mindsmobility.com

555 Washington Ave., Suite 310

St. Louis, MO 63101

Ph: 866-709-1332 Fax: (877)-426-5697