













## S-Ergo 115-TP

### *Ultra Lightweight Ergonomic Manual Transporter*

- S-Style Ergonomic Seat
- High strength, weighs only 22 lbs. (w/o footrests)
- Silver ¼" Removable Anti-bacterial cushion
- Folding backrest/folding seat for transporting in vehicle or travel
- Fixed armrests w/ wider concaved armpads
- Swing In & Away footrests allowing maximum safety while entering or exiting the chair.
- "Tube-in Center" foot-plate, assures better side leg support
- Extending push to lock manual brakes
- Pouch for carrying small items attached to upholstery behind backrest.
- 14" Rear polyurethane, high profile, flat free wheels
- 7"x1" Polyurethane front caster
- 16x17 or 18x17 Seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black Breathable Mesh Bottom
- Weight Capacity of 250 lbs

#### *Features for S-Ergo Series*

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slipage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/concaved armpad, providing comfort & support

	items	F/R wheels											Max. User Weight
	Width												
S-Ergo 115	~ 16"x17"	7"/14"	16"	17"	22 ½"	12"	8"	35 ¾"	18 ½"	16"	36"	22 lbs.	250 lbs.
S-Ergo 115	~ 18"x17"	7"/14"	18"	17"	23 ½"	12"	8"	35 ¾"	18 ½"	16"	36"	22 lbs.	250 lbs.

*\*Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam, UV Resistant Canopy, Rigidifying Push-Bar, Reflecting Spoke Guard, Foldable/Adjustable Headrest, Frog Leg Suspension, Magnesium Wheel, Vulcan Grip Wheel*



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335