



ErgoFlight S-2512 ***Ultra Lightweight Ergonomic Wheelchair***

- 14.5 lb Aluminum Frame. This Manual Wheelchair Competes Against Titanium Weight!
- Patented S-Style Ergonomic Seat Frame
- 6061 T-6 Aircraft-grade Aluminum Only 19.8 lbs!
- Companion Brakes Located on Handles for Added Safety and Comfort.
- Built in Silver Aegis Anti-bacterial Seat and Back Cushion
- Folding Backrest/Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/Ergonomic Arm Pads
- Removable Swing-in & Out Footrests
- 20" Rear Polyurethane, High Tread, Flat Free Wheels
- 6" x 1" Polyurethane Front Casters
- 18" x 17", 16" x 17" Seat Width available
- Available Quick Release Rear Wheels
- Frame Color: Pearl Silver
- Weight Capacity: 220 lbs
- Limited Lifetime Warranty on Frame

Taking Aluminum to Titanium Heights

Titanium Weight at Aluminum price!

***U.K. ISO 7176/19 Dynamic
Crash Test Passed***

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	Wt. Cap. (lb.)
S-2512F18SS	18	17	8	18.5	17.5	36	27	15	38	19	220
S-2512F16SS	16	17	8	18.5	17.5	36	25	15	38	19	220
S-2512Q18SS	18	17	8	18.5	17.5	36	27	15	38	19	220
S-2512Q16SS	16	17	8	18.5	17.5	36	25	15	38	19	220

****Available Options ~ Quick Release wheels, Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam cushions, UV Resistant Canopy, Rigidifying Push Bar, Foldable/Adjustable Headrest***

karma KARMAN HEALTHCARE Inc.

19255 San Jose Avenue, City of Industry, CA 91748-1418

Order. 1-800-805-2762 / Tel.626-581-2235 / Fax.626-581-2335