



KARMAN HEALTHCARE
 19255 San Jose Ave.
 City Of Industry, CA 91748
Tel: 626-581-2235 **Fax:** 626-581-2335
Toll Free: 800-805-2762



QR Reader Smartphone Info



S-2512 Ergo Flight SERIES ORDER FORM

Submitting for QUOTE ORDER

Mark For: _____
 Date: _____
 Dealer Acct#: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: (____) _____ FAX: (____) _____
 Confirmation Email: _____

Purchase Order #: _____
 Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: _____

The HCPCS CODES are not intended to be, nor should be considered billing or legal advice. They are provided based on code verification through PDAC or our interpretations of billing. The use of HCPCS CODES does not ensure coverage or payment. For coverage information, please verify the appropriate payer's coverage policy

S-2512 Series:

S-2512 Ergonomic Wheelchair std features: Swing away footrest, folding backrest, cloth side guard, frame weight 14.5 lbs, weight w/o footrest 19.8 lbs, weight capacity 220 lbs, Anti bacterial seat and back upholstery, companion brakes.

BASE FRAME				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	S-2512F16SS	K0004	16x17 Pearl Silver frame	\$1,227.00
<input type="checkbox"/>	S-2512F18SS	K0004	18x17 Pearl Silver frame	\$1,227.00
<input type="checkbox"/>	S-2512Q16SS	K0004	16x17 Pearl Silver frame w/Quick release axles	\$1,437.00
<input type="checkbox"/>	S-2512Q18SS	K0004	18x17 Pearl Silver frame w/Quick release axles	\$1,437.00

SEAT CUSHIONS				
<input type="checkbox"/>	CU-ERGO-18	E2601	18x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-ERGO-16	E2601	16x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-FT-18	E2601	18x16 seat foam cushion	\$72.00
<input type="checkbox"/>	CU-FT-16	E2601	18x16 seat foam cushion	\$72.00
<input type="checkbox"/>	CU-GFR-18	E2601	18x16 seat gel cushion	\$98.00
<input type="checkbox"/>	CU-GFR-16	E2601	18x16 seat gel cushion	\$98.00

BACK CUSHIONS				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	CU-ERGO-18-back	E2611	18" wide memory foam 2" back cushion	\$116.00
<input type="checkbox"/>	CU-ERGO-16-back	E2611	16" wide memory foam 2" back cushion	\$116.00
<input type="checkbox"/>	BKF-1817	E2611	18" wide back foam cushion	\$98.00
<input type="checkbox"/>	BKF-1617	E2611	16" wide back foam cushion	\$98.00

OPTIONS & ACCESSORIES				
<input type="checkbox"/>	SB22	E0978	Clamp style	\$19.00
<input type="checkbox"/>	SB88	E0978	Wrap around seat belt 1 piece	\$19.00
<input type="checkbox"/>	SB99-48	E0978	Auto Style	\$24.00
<input type="checkbox"/>	HR-FLD-115	E0955 + E1028	16-18" width Foldable Headrest with hardware	\$238.00
<input type="checkbox"/>	AT-2512	E0971	Anti Tippers-pair	\$56.00
<input type="checkbox"/>	HL-2512	E0951	Heel Loop	std
<input type="checkbox"/>	LS16	K0038	16" width leg strap-calf support	\$11.00
<input type="checkbox"/>	LS18	K0038	18" width leg strap-calf support	\$11.00
<input type="checkbox"/>	BKR-EXT-S-105-16		16" width back rest extension	\$96.00
<input type="checkbox"/>	BKR-EXT-S-105-18		18" width back rest extension	\$96.00
<input type="checkbox"/>	CANO-115		Canopy for S-100 series. UV protectant	\$248.00
<input type="checkbox"/>	PB-115-16		Foldable rigidifying push bar for 16" width	\$88.00
<input type="checkbox"/>	PB-115-18		Foldable rigidifying push bar for 16" width	\$88.00
<input type="checkbox"/>	CP2	E1399	Universal carry pouch-small	\$11.00
<input type="checkbox"/>	CP4	E1399	Universal carry pouch-large	\$19.00