

S-Ergo 106


Lightweight Ergonomic Manual

Wheelchair W/ User Angle Adjustable Back



Features for S-Ergonomic Series

- Ergonomic ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slippage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/Concaved arm pad, providing comfort & support



- Medicare code K0004
- S-Style Ergonomic Seat
- High strength, weighs only 27 lbs. (w/ footrests)
- Silver ¼" Anti-Bacterial cushion, able to be removed to wash
- Has 3 easy angle adjustments: "Working" Angle (95°), "Moving" Angle (100°) & "Resting" Angle (105°)
- "Rigidifying" push-bar & folds to strengthen back support, comes standard on this model!
- Folding backrest/folding seat for transporting vehicle or travel
- Extending push to lock manual breaks
- Fixed armrests w/wider concaved armpads
- Fixed footrests
- "Tube-in Center" foot-plate, assures side leg support
- Pouch for carrying small items attached to upholstery behind backrest.
- 24" Rear polyurethane, high tread, flat free wheels
- 7x1" Front Polyurethane casters
- 16x17, 18x17 Seat widths available
- Anti-Tippers standard on this model!
- Frame Color: Pearl Silver Only
- Upholstery: Black breathable mesh bottom/ AEGIS top

	seat width	seat depth	armrest height	seat height	back height	overall height	overall open width	folded width	overall length	weight w/o riggings	weight capacity
Model No.	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(in.)	(lb.)	(lb.)
S-106F16	16	17	8	19	17	36	23	12	39	27	250
S-106F18	18	17	8	19	17	36	25	12	39	27	250

**Available Options ~ Seat belt, Memory Foam Cushion, Adjustable UV Resistant Canopy, Rear Wheel Reflector/Spoke Guard, Magnesium Wheels, Vulcan Grip Wheels, Foldable/ Adjustable Headrest, Frog Leg Suspension*



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335