




Aqua Creek Products

The Total Solution for Pool and Spa Access



The Titan 600 Lift

Item# F-600TL

The Titan 600 Lift shares a lot of the same great features as the Revolution Lift but is designed to accommodate patrons up to 600 lbs. No other ADA compliant lift on the market has this weight capacity. Now independently verified by Wyle Labs to be ADA compliant.

This lift provides the ability to service multiple bodies of water as well as both above-ground spas and in-ground pools and spas.

FEATURES:

- 600 Pound (272 kg) weight capacity
- Independently verified to ADA standards
- Rechargeable battery with wall mount charger
- Clears 30" high walls, up to 10" thick
- 360° powered rotation in either direction
- Adjustable footrest & lap belt
- Dual flip-up arm rests
- Water resistant push-button remote control
- 5 year structural warranty
- 5 year pro-rated electronics warranty



U.S. Patent No. D507,769 S

OPTIONS:

- **Anchor, standard** (additional) - **Item #F-612TLA**
- **Upgrade Pack** - **Item #F-700RLAP**
Includes battery, headrest, chest strap & cover
- **Headrest** - **Item #F-422PLH-2**
- **Chest strap** - **Item #F-423CS**
- **Lockable Cover** - **Item #F-720RLC**
- **Transport Cart** - **Item #F-614TTC**
- **Pull-out leg rest** - **Item #F-10SAR**
- **Cycle Attachment** - **Item #F-019CA**
- **Solar charger** - **Item #F-045SCH**



Transport Cart

PROUDLY SOLD BY:

The Titan 600 Lift

ADA Compliant