



CHERRY RED



PARADISE BLUE



SOLAR FLARE ORANGE  
*Limited Edition*



ENVY GREEN  
*Limited Edition*



SATIN SILVER  
*Limited Edition*



SUNBURST YELLOW  
*Limited Edition*

**Seat Color**



Black Vinyl



*Available in 6 New  
Vibrant Colors with  
a Larger Foot Plate and  
New Shroud Design!*

**Disassembles  
Quickly &  
Easily!**

LiteRider<sup>®</sup> GOLDEN  
ENVY  
2018

The all new LiteRider Envy Power Wheelchair is both stylish and durable with a tubular design that is unique to the Golden LiteRider family! The LiteRider Envy, featuring 2-22 amp hour batteries for an operating range of up to 15.5 miles, is Golden's Retail Cash Power Wheelchair! The LiteRider Envy has a larger footplate and an impressive turning radius of 28.5", and like the other LiteRiders in the family, it disassembles quickly and easily for transport. The heaviest piece is only 37 lbs., making storage convenient in the trunk of any vehicle. When you are looking for a compact power wheel chair that is fun to drive and easy to transport – look no further than the Golden LiteRider Envy!



<b>Model Number</b>	<b>GP162</b>
Medicare Code	Not Coded
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	3.5 mph
Operating Range <sup>1</sup>	15.5 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Type Batteries	2-22 AH
Colors	Cherry Red, Paradise Blue, Solar Flare Orange, Envy Green, Satin Silver, Sunburst Yellow
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG VR-2 50 Amp Remote
<b>Dimensions:</b>	
Length	36"
Width	22.5"
Height (ground to top of back of seat)	33" - 34"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	20.5" - 21.5"
Deck to Top of Seat	N/A
Number of Seat Height Positions	2
Size of Increments	1"
<b>Weight of Unit: (assembled)</b>	
Front Half	37 lbs.
Rear Half	36 lbs.
Seat w/Arms	23 lbs.
Batteries	27 lbs./pack
<b>Tires:</b>	
Drive Wheels	Flat Free
Casters	9"
Anti-Tip Wheels	6"
<b>Standard Seat:</b>	
Back Height - no headrest	16"
Back Height with headrest	N/A
Width x Depth	17" x 16"
Color	Black Vinyl
<b>Optional Seats:</b>	
Back Height - no headrest	16"
Back Height with headrest	N/A
Width x Depth	20" x 18"
Color	Black Vinyl
<b>Seat Options:</b>	
Pneumatic Rise Seat Option	N/A
Electric Rise Seat Option	N/A
<b>Footrest:</b>	
Height adjustable	Yes
Angle adjustable	Yes
<b>Warranty:</b>	
Frame	Five Years
Drive Train	13 Month Warranty
Electronics <sup>2</sup>	13 Month Warranty
Batteries	12 Month Warranty

<sup>1</sup> Battery range tested at 200 lbs. and will vary due to rider weight, drive surface and terrain.

<sup>2</sup> Electronics warranty excludes batteries.

## Disassembles easily for transport!



1. Remove seat



2. Pull out battery pack



3. Grasp the frame and locking handle; lift the handle to unlock, and lift the frame up and off the drive train.



Perfect for transport!

## Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Pack N' Go
- Seat Belt (Optional)

Some accessories require a mounting bracket or mounting clips.