



Harmar

MOBILITY MADE EASY

SL500HD

Summit High Capacity Stairway Lift

The breakthrough Harmar SL500HD has the heaviest carrying capacity (500lbs) in the industry as well as a 23" wide, high-backed contoured seat that helps ensure a comfortable, secure ride. Also standard is reinforced footrest for added stability and heavy-duty armrests that help provide more secure transfers.



Specifications:

Capacity:	500Lbs
Seat Size:	23" W x 18.5" Deep x 16" Tall seat back
Seat Rotation:	90 degrees at both to and bottom landings
Power Source:	115 VAC household current
Warranty:	3-Year Parts, 5-Year Drivetrain