


Our MOST POPULAR SAFETY BED with safety rails ranging up to 27" above the mattress!

SleepSafe® II (Medium Bed)



SPRING TO SAFETY RAIL HEIGHT

Fixed Foundation with Box Spring	23 3/4", 20" or 16"
Fixed Foundation with Bunkie Board	33" or 27 3/4"
Manual Articulation	27", 23" or 19 1/4"
Manual HI-LO	LO 27 1/4" to HI 15 1/2" or LO 22" to HI 10 1/4"
Manual HI-LO Combo	LO 27 1/4" to HI 15 1/2" or LO 22" to HI 10 1/4"
Manual Electric Combo	LO 27 1/4" to HI 15 1/2" or LO 22" to HI 10 1/4"
Electric Articulation	27", 23" or 19 1/4"
Electric HI-LO Combo	LO 27 1/4" to HI 17" or LO 22" to HI 11 3/4"

SleepSafe® II (Medium Bed) Measurements and Specifications:

SAFETY RAIL HEIGHT

Fixed Foundation with Box Spring	17 3/4", 14" or 10"
Fixed Foundation with Bunkie Board	27" or 21 3/4"
Manual Articulation	21", 17" or 13 1/4"
Manual HI-LO	LO 21 1/4" to HI 9 1/2" or LO 16" to HI 4 1/4"
Manual HI-LO Combo	LO 21 1/4" to HI 9 1/2" or LO 16" to HI 4 1/4"
Manual Electric Combo	LO 21 1/4" to HI 9 1/2" or LO 16" to HI 4 1/4"
Electric Articulation	21", 17" or 13 1/4"
Electric HI-LO Combo	LO 21 1/4" to HI 11" or LO 16" to HI 5 3/4"

OVERALL OUTSIDE DIMENSIONS

Twin	42 1/4" W x 80 1/2" L x 50 3/4" H
Full	56" W x 80" L x 50 3/4" H
Queen	65" W x 80" L x 50 3/4" H

MATTRESS SIZE (Without pads / With pads / Length)

Twin	38" W / 36" W / 77" L
Full	52" W / 50" W / 77" L
Queen	61" W / 60" W / 77" L

MATTRESS TO FLOOR HEIGHT

Fixed Foundation with Box Spring	28", 32" or 35 1/2"
Fixed Foundation with Bunkie Board	18 3/4" or 24"
Manual Articulation	25", 28 3/4" or 32 1/2"
Manual HI-LO	LO 24 1/2" to HI 36 1/4" or LO 29 3/4" to HI 41 1/2"
Manual HI-LO Combo	LO 24 1/2" to HI 36 1/4" or LO 29 3/4" to HI 41 1/2"
Manual Electric Combo	LO 24 1/2" to HI 36 1/4" or LO 29 3/4" to HI 41 1/2"
Electric Articulation	25", 28 3/4" or 32 1/2"
Electric HI-LO Combo	LO 24 1/2" to HI 34 3/4" or LO 29 3/4" to HI 40"

FLOOR TO RAIL HEIGHT

45 3/4"

UNDER BED FLOOR CLEARANCE

8"

TRANSFER HEIGHT

With rails attached / detached 32" / 30 1/2"

AVAILABLE OPTIONS

Padding over or around windows (6 colors available);

Solid clear or mesh windows; IV pole; Tubing Channel Access

Cut Out; Headboard and/or footboard windows

Wood Choices

SleepSafe® Beds can be built with wood or color finishes:



Color Choices



All bed finishes are highly durable surfaces available in wood and colors. Colors can be used alone or in combination with others.

Multi-Color...



...Mix & Match!

SleepSafe® Beds are designed to address the issues of falls and entrapment for those with special needs. Models include SleepSafe® -Low Bed, SleepSafe® II -Medium Bed, & SleepSafe® -Tall Bed, each offering more safety rail to mattress height.



SleepSafeBed.com / 866-852-2337

FDA Registered / Made in USA / Call for DME information

SleepSafe® Beds are built to order, offering twin, full or queen size safety beds in fixed, articulating or HiLo foundations in electric, and/or manual configurations. padding; wood or color selections, and other helpful options for those with special needs.