

IMPORTANT! When ordering your bed, please indicate with or without padding. If padded, see options below.

SleepSafe® Bed Padding Options



Padding Color Choices



Materials:

- 70 Denier Nylon Fabric Cover with a Fluid Resistant Healthcare Coating
- 3/4" High Density Foam Padding

Padding features:

- Bacteria resistant
- Malodor Protection
- Fluid Resistant
- Padding extends to the top of the headboard and footboard.
- Padding extends up and over the safety rails and attaches to the outside.
- Padding covers the corner posts.
- Fungal resistant
- Flame Retardant
- Tear Resistant
- Non-allergenic
- Stain Resistant
- Anti-Static

Padding questions?


- Please contact customer service for any specific information about padding.

Padding attachment:

- Trim-Lok secures padding to the Safety Rails.
- Velcro® secures padding to the inside panels.

Padding maintenance:

- For cleaning, use warm water and a mild cleanser. If the padding is over the windows, it can be removed and cleaned. If the padding is around the windows, please do not remove. It is best to clean this padding in its place.



Important notes:

Padding is installed at the factory and will cover ALL 4 sides. Padding only 1 or 2 sides is NOT an option. Padding removal will cause a gap between the mattress and the bed's interior.



SleepSafeBed.com / 866-852-2337

Made in USA / Call for DME information / FDA Registered