


Our ORIGINAL SAFETY BED with safety rails ranging up to 17¼" above the mattress!

SleepSafe® (Low Bed)



SleepSafe® (Low Bed) Measurements and Specifications:

SAFETY RAIL HEIGHT (Top of rail to top of mattress)

| | |
|------------------------------------|------------------------------|
| Fixed Foundation with Box Spring | 12" or 8¾" |
| Fixed Foundation with Bunkie Board | 17¼" or 14" |
| Manual Articulation | 15" or 11¾" |
| Manual HI-LO | LO 11½" to EVEN with rail |
| Manual HI-LO Combo | LO 11½" to EVEN with rail |
| Manual Electric Combo | LO 11½" to EVEN with rail |
| Electric Articulation | 15" or 11¾" |
| Electric HI-LO Combo | LO 11½" to HI 1¼" below rail |

MATTRESS TO FLOOR HEIGHT

| | |
|------------------------------------|--------------------|
| Fixed Foundation with Box Spring | 23¾" or 27" |
| Fixed Foundation with Bunkie Board | 18½" or 21¾" |
| Manual Articulation | 20¾" or 24" |
| Manual HI-LO | LO 24¼" to HI 36" |
| Manual HI-LO Combo | LO 24¼" to HI 36" |
| Manual Electric Combo | LO 24¼" to HI 36" |
| Electric Articulation | 20¾" or 24" |
| Electric HI-LO Combo | LO 24¼" to HI 34½" |

OVERALL OUTSIDE DIMENSIONS

| | |
|-------|--------------------------|
| Twin | 42¼" W x 80½" L x 40¾" H |
| Full | 56" W x 80" L x 40¾" H |
| Queen | 65" W x 80" L x 40¾" H |

FLOOR TO RAIL HEIGHT

| | |
|--|------|
| | 35¾" |
|--|------|

UNDER BED FLOOR CLEARANCE

| | |
|--|----|
| | 8" |
|--|----|

TRANSFER HEIGHT

| | |
|--------------------------------|-------------|
| With rails attached / detached | 22¼" / 20¾" |
|--------------------------------|-------------|

AVAILABLE OPTIONS

Padding over or around windows (6 colors available);
Solid clear or mesh windows; IV pole; Tubing Channel Access
Cut Out; Headboard and/or footboard windows

MATTRESS SIZE (Without pads / With pads / Length)

| | |
|-------|-----------------------|
| Twin | 38" W / 36" W / 77" L |
| Full | 52" W / 50" W / 77" L |
| Queen | 61" W / 60" W / 77" L |

SleepSafe® Beds are designed to address the issues of falls and entrapment for those with special needs.

Models include SleepSafe® - Low Bed, SleepSafe® II - Medium Bed, & SleepSafe® - Tall Bed, each offering more safety rail to mattress height.



SleepSafeBed.com / 866-852-2337

FDA Registered / Made in USA / Call for DME information

SleepSafe® Beds are built to order, offering twin, full or queen size safety beds in fixed, articulating or HiLo foundations, in electric, and/or manual configurations; padding; wood or color selections, and other helpful options for those with special needs.

Wood Choices

SleepSafe® Beds can be built with wood or color finishes:



Color Choices



All bed finishes are highly durable surfaces available in wood and colors. Colors can be used alone or in combination with others.

Multi-Color...

