

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (with riggings)	Product Weight	Carton Shipping Weight	Weight Limit
16"	21"	12.5"	36"	17.5"-19.5"	16"	16"	16"	42"	35.8 lbs.	39.8 lbs.	250 lbs.
18"	28"	12.5"	36"	17.5"-19.5"	16"	18"	16"	42"	37 lbs.	41 lbs.	250 lbs.
20"	30"	12"	36"	17.5"-19.5"	16"	20"	16"	42"	39 lbs.	43 lbs.	250 lbs.