

## Scout 4

HCPCS: K0800-1, K0800



FIGURE A



FIGURE B




FIGURE C




FIGURE D




FIGURE E

- Quick and easy disassembly (Figure A)
- Color-through color panels hide scuffs and scratches
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Delta tiller
- Easy free wheel operation
- Available with 12AH or 20AH batteries
- 20AH batteries increase scooter's maximum range

### PRODUCTS

### Description

|              |                                       |
|--------------|---------------------------------------|
| SFSCOUT4     | Scout 4-Wheel Scooter, Red/Blue, 1/ea |
| SFSCOUT4-EXT | Scout 4-Wheel 20AH Battery, Red/Blue  |

### SPECIFICATIONS

### Description

|                       |   |
|-----------------------|---|
| Base Weight           | 58 lbs                                  |
| Batteries (pair)      | 12V x 12AH, 12V x 20AH (EXT Model Only) |
| Battery Weight (pair) | 20 lbs                                  |
| Brakes                | Electromagnetic                         |
| Charger               | 2A Offboard                             |
| Climbing Angle        | 6°                                      |
| Controller            | CTE 45A                                 |
| Dimensions            | 42.5" (L) x 21.5" (W)                   |
| Floor to Seat Height  | 17" - 19"                               |
| Front Wheels          | 8" x 2" Flat-free                       |
| Ground Clearance      | 2.5"                                    |
| HCPCS                 | K0800                                   |
| Heaviest Piece Weight | 39 lbs                                  |
| Maximum Range         | 9 miles (12AH)   15 miles (20AH)        |
| Motor                 | 24V x 270W x 4700 rpm                   |
| Rear Wheels           | 8" x 2" Flat-free                       |

| SPECIFICATIONS   | Description         |
|--|---------------------|
| Seat Dimensions  | 16.5" (W) x 16" (D) |
| Seat Weight  | 16 lbs              |
| Top Speed  | 4.25 mph            |
| Total Weight (excl. hand control)                        | 94 lbs              |
| Turning Radius   | 53.75"              |
| Warranty on Batteries                                    | 12 Months           |
| Warranty on Electronic Controller/Drive Train Components | 24 Months           |
| Warranty on Frame  | Lifetime            |
| Weight Capacity  | 300 lbs             |
| Wheel Base Dimensions (L x W)                            | 31.9" x 16.9"       |

| ACCESSORIES | Description                        |
|-------------|------------------------------------|
| SF8050      | Sun Shade for Scooter,1/ea, EA     |
| SF8020      | Rear Basket Set, LType, 1/ea, EA   |
| SF8030      | Crutch Holder Set, SType, 1/ea, EA |
| SF8010      | Oxygen Holder Set, L Type,1/ea, EA |

| PARTS         | Description                                  |
|---------------|--|
| LRM402109     | 12AH Battery,POV,1/ea, EA                    |
| LRM412103-11  | Speed Control Knob, EA                       |
| BC350103      | Connecting Rod Set,Bobcat 4,1set/pkg, PKG    |
| C08-069-01000 | Tiller Lock Set, Bobcat X3 / X4, 1/ea, EA    |
| S31026        | Anti Tip Wheel,Phoenix/Bobcat,1/ea, EA       |
| S35003-P      | PlasticBasketBobcat,1ea, EA                  |
| S31052        | Armrest Knob POV 1/ea, EA                    |
| BC20ACB       | Bobcat Breaker 20A with Cap,1ea, EA          |
| BC35051       | BasketHook Set,Phoenix,Scout,Ventura,1/e, EA |
| BC953014      | Arm Pad Only,Scout, Bobcat, 1/ea, EA         |
| BC35040L      | Bobcat Complete Arm Rest,Left,1ea, EA        |
| BC35040R      | Bobcat Complete Arm Rest,Right,1ea, EA       |
| BC35045       | Bobcat Seat,Black,1ea, EA                    |
| BC31021       | Tiller Bearing Set,Bobcat 4,1set/pkg, PKG    |
| LRM412123     | 21AH Battery,POV, 1/ea, EA                   |
| LRM302111     | 2ah Battery Charger, POV,1/ea, EA            |
| BC31006       | Bobcat Ignition Switch,1ea, EA               |
| BC31001       | Scout,Bobcat, Speed Dial, 1/ea, EA           |
| BCKEY         | Key, Bobcat, Scout, 1ea, EA                  |
| BC350112      | Bobcat Transaxle Only,1ea, EA                |