



Rise




Regal
PR-751TY




Saddle



Orion



Havana



VINO



Winston
PR-410



Royal
PR-752



Space Saver
PR-906

Signature Series

The Signature Series offers four models, each with its own unique, luxurious style. Choose from the Regal with dual open arm storage; the Winston with overstuffed, generous proportions; the Royal with its plush, three-pillow back; or the Space Saver that requires only 6" from the wall for a full recline. The Signature Series offers a style like no other.



Rise.

Roll.

Ride.



GOLDEN
technologies

YOUR LIFE IN MOTION



SIGNATURE SERIES



	Winston	Regal	Royal	Space Saver
Model Number	PR-410	PR-751TY	PR-752	PR-906
Size	Large	Medium	Medium	Medium
Weight Capacity	375 Lb.	375 Lb.	375 Lb.	250 Lb.
Recline Positions	3	3	3	3
Standard Fabrics	4	4	4	4
Fabric Options	ALL	ALL	ALL	ALL
Chaise Seat	YES	YES	YES	YES
Back Type	Biscuit	3 Pillow	3 Pillow	Seam
Head Rest Cover	YES	YES	YES	YES
Arm Covers	NO	NO	YES	NO
Battery Back Up	YES	YES	YES	YES
Overall Width	43"	36"	32.5"	35.5"
Overall Height	46"	44"	44"	43"
Width Between Arms	21.5"	22"	21"	18"
Floor to Top of Seat	21"	20"	19.5"	22"
Seat Depth	21"	21.5"	21"	22"
Seat to Top of Back	29"	28"	28"	27"
Seat to Top of Arm	6.5"	8"	6"	8"
Floor to Back Edge of Seat Raised	29"	28"	27.5"	30"
Floor to Front Edge of Seat Raised	24.5"	24"	23.5"	24"
Seat Angle Raised	20°	18°	15°	24°
Floor to Top of Foot Rest Reclined	23.5"	24"	23.5"	21.5"
Floor to Top of Back Reclined	36"	33"	33"	37"
Back Angle Closed	110°	110°	115°	120°
Back Angle Reclined	158°	140°	145°	140°
Foot Rest Angle Closed	90°	90°	90°	90°
Foot Rest Angle Opened	180°	180°	180°	180°
Distance Required from Wall Reclined	19"	18"	19"	6"
Extension Reclined from Back of Seat to End of Foot Rest	45"	42"	42"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES	1 Yr.
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.	1 Yr.
Pro Rated Years 4 Thru 7	YES	YES	YES	NO
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.	1 Yr.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.	1 Yr.
Removable Back Standard	YES	YES	YES	YES
Smart Tek System Standard	YES	YES	YES	YES
Customizable	NO	NO	NO	NO
Left Hand Control Option Available	YES	NO	YES	YES
Extra Pocket Option Available	YES	YES	YES	YES
Heat and Massage Option Available	YES	YES	YES	YES
Foot Rest Extension Option Available	NO	NO	YES	NO
Straight Lift Option Available	YES	YES	YES	NO

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68
Listing # E112417
FNAE 80-214A



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com

SHT-SIGN-111513