

Size Chart

Special Tomato Soft-Touch Liners	Size 1	Size 2	Size 3	Size 4	Size 5
Standard Seat Cushion Width	8" (20 cm)	10" (25 cm)	12" (30 cm)	14" (35.5 cm)	16" (40.6 cm)
Standard Seat Cushion Depth	8" (20 cm)	10" (25 cm)	12" (30 cm)	14" (35.5 cm)	16" (40.6 cm)
Hip-Flex® Seat Cushion Width	8" (20 cm)	10" (25 cm)	12" (30 cm)	14" (35.5 cm)	16" (40.6 cm)
Hip-Flex® Seat Cushion Depth	8" (20 cm)	10" (25 cm)	12" (30 cm)	14" (35.5 cm)	16" (40.6 cm)
Back Cushion Width	7.5" (18.75 cm)	8.5" (21.25 cm)	10" (25 cm)	13" (33cm)	14" (35.5 cm)
Back Cushion Height	10" (25 cm)	13.5" (33.75 cm)	14.5" (36.25 cm)	15.5" (39.3 cm)	16" (40.6 cm)
Overall Dimensions of Standard Seat Cushion (DxW)	10.5"x11" (27x28 cm)	13"x12" (33x31 cm)	15"x14" (38x36 cm)	17"x17" (43x43 cm)	19.5"x19.5" (49.5x49.5 cm)
Overall Height of Standard Seat Cushion	3.5" (8.75 cm)	4" (10 cm)	4.5" (11.25 cm)	4.5" (11.25 cm)	4.5" (11.25 cm)
Overall Dimensions of Hip-Flex® Seat Cushion (DxW)	10.5"x11" (27x28 cm)	13"x12" (33x31 cm)	15"x14" (38x36 cm)	17"x17" (43x43 cm)	19.5"x19.5" (49.5x49.5 cm)
Overall Height of Hip-Flex® Seat Cushion	3.5" (8.75 cm)	4" (10 cm)	4.5" (11.25 cm)	4.5" (11.25 cm)	4.5" (11.25 cm)
Overall Dimensions of Back Cushion (H x W)	10"x10.5" (26x27 cm)	13.5"x12.25" (33.75x31.25 cm)	14.5"x13.25" (36.25x33.5 cm)	15.5"x14.5" (39.3x36.8 cm)	16.5"x16.5" (41.9x41.9 cm)
Weight of Standard Seat Cushion	1.5 lbs. (0.68 kg)	2 lbs. (0.91kg)	2.5 lbs. (1.13 kg)	3 lbs. (1.36 kg)	5 lbs. (2.27 kg)
Weight of Hip-Flex® Seat Cushion	1.5 lbs. (0.68 kg)	2 lbs. (0.91 kg)	2.5 lbs. (1.13 kg)	3 lbs. (1.36 kg)	5 lbs. (2.27 kg)
Weight of Back Cushion	1.5 lbs. (0.68kg)	2 lbs. (0.91 kg)	2.5 lbs. (1.13 kg)	3 lbs. (1.36 kg)	5 lbs. (2.27 kg)