

Sling Specifications

	SMALL	MEDIUM	LARGE	BARIATRIC
	A - circumference around the shoulders (inches)			
Universal Split Leg	28-45	37-57	47-67	57-80
Hammock	28-43	35-51	45-60	N/A
	B - circumference around the chest (inches)			
Pediatric	18-28	N/A	N/A	N/A
Sani/Toileting	18-40	26-46	28-50	34-73
Sit-to-Stand	26-42	30-46	34-53	46-64
	Body Weight (lbs)			
Universal Split Leg	85-140	140-200	200-400	400-600
R-Series	N/A	100-180	160-280	200-400
*R140/R141 300-600lbs				

Above dimensions are guidelines and may vary by user.

