

SMOOV

POWER ASSIST FOR ACTIVE WHEELCHAIRS

NEW



SMOOV one.

The perfect solution designed for active wheelchair users.

The SMOOV one is an innovative electric drive that you can easily dock on to and off of your wheelchair whenever you need some extra power.

This modern and sporty power unit is designed to significantly increase the driver's mobility and range.

The ergonomic control unit is easy to operate and can be mounted on either the left or right side.



Speed up to
6 miles/hr *



Range up to
12 miles *



Slope up to
16 % *



Weight
16 lbs

*varies based on the weight of the user, temperature, and terrain

Find out more about SMOOV at
www.smoov.com

Live life – without limits.

The SMOOV one is the perfect power assist designed for active wheelchair users. Thanks to its discreet drive and control unit brackets, your wheelchair can still be used manually: simply dock and undock.

With a range of up to 12 miles,* this unique drive unit gets you to your destination in style. You can easily charge the battery from any outlet.

With a maximum speed of 6 miles/h,* the SMOOV is your adaptable everyday companion. You can even drive over curbs and edges thanks to the flexible wheel and the 360° swivel fork.

The SMOOV can take you to places that once seemed beyond reach. It can be relied on to handle some of the most difficult terrains. Even hills and slopes of up to 16%* are no problem with the SMOOV.



App connected

Download now.



SMOOV

Alber USA

1005 International Dr
Oakdale, PA 15071
Phone (724) 695-7822
Fax (724) 695-3710
info@alber-usa.com

