



- Hazelnut
- Sterling
- Copper
- Calypso
- Shiraz



Are you short on space, but searching for extreme comfort? You don't have to look any further than our sizing collection of space saving power lift recliners. With a seamed back and forward-gliding recline motion, this design is offered in small, medium, and large sizes and requires only 4", 5" or 6" inches from any wall. The Golden Space Saver will fit you perfectly in almost any sized room!

PR931





PR931

S I G N A T U R E S E R I E S

	Space Saver	Space Saver	Space Saver
Model Number	PR931S	PR931M	PR931L
Size	Small	Medium	Large
Weight Capacity	350 LBS	350 LBS	350 LBS
Recline Positions	Space Saver	Space Saver	Space Saver
Standard Fabrics	5	5	5
Fabric Options	ALL	ALL	ALL
Chaise Seat	YES	YES	YES
Coil Springs	NO	NO	NO
Back Type	Seam	Seam	Seam
Head Rest Cover	NO	NO	NO
Arm Covers	NO	NO	NO
Battery Back Up	YES	YES	YES
Overall Width	30.5"	32.5"	35"
Overall Height	39.5"	41.5"	43.5"
Width Between Arms	17"	20"	22.5"
Floor to Top of Seat	18"	19"	20"
Seat Depth	20"	21"	22.5"
Seat to Top of Back	26.5"	28.5"	30.5"
Seat to Top of Arm	5.5"	5.5"	6"
Seat Angle Raised	20°	20°	20°
Floor to Top of Foot Rest Reclined	22"	22"	22"
Floor to Top of Back Reclined	34.5"	35.5"	37"
Back Angle Closed	105°	105°	109°
Back Angle Reclined	132°	135°	135°
Foot Rest Angle Closed	90°	90°	90°
Foot Rest Angle Opened	180°	180°	180°
Distance Required from Wall Reclined	4"	5"	6"
Extension reclined from Back of Seat to end of Foot Rest	41"	42.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.
Removable Back Option Available	YES	YES	YES
Smart Tek System Standard	YES	YES	YES
Customizable	NO	NO	NO
Left Hand Control Option Available	NO	NO	NO
Extra Pocket Option Available	NO	NO	NO
Heat & Massage Option Available	YES	YES	YES
Gel Memory Foam Option Available	NO	NO	NO
Power Pillow Option Available	NO	NO	NO
AutoFlex Option Available	NO	NO	NO
Foot Rest Extension Option Available	NO	NO	NO
Straight Lift Option Available	NO	NO	NO
User Height (suggested)	5'0" - 5'3"	5'4" - 5'10"	5'11" - 6'2"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68
Listing # E112417
FNAE 80-214A



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com

SHT-PR931-092616