




Item #	Product HCPCS	UOM
SFEX2317FS-12	KO800	1/ea
SFEX2317FS-21	KO800	1/ea



accessories

SF8010	Oxygen Tank Holder	1/ea
SF8020	Rear Basket	1/ea
SF8030	Cane/Crutch Holder	1/ea



White



Black



Silver



Blue



Pink

Spitfire EX2 3-Wheel Travel Scooter with 12AH Battery

- 6 sets of interchangeable color panels in red, blue, silver, black, white and pink
- Swivel mount cup holder included
- Extra bright LED headlight
- Charge ports located on the tiller and battery pack for easy charging
- Height adjustable swivel seat
- Adjustable padded armrests with reflectors
- Large plastic basket with removable liner
- Disassembles into 5 easy to carry pieces
- Delta tiller for easy throttle control

specifications

Weight Capacity	300 lbs.	Batteries	12AH 21AH
Max Speed	5 mph	Battery Charger	2A Offboard
Maximum Range	9 Miles 15 Miles	Brakes	Electromagnetic
Max. Climbing Angle	6°	Freewheel Mode	Yes
Turning Radius	54"	Front Wheels	7.5" x 2"
Ground Clearance	3"	Rear Wheels	8.5" x 3"
Overall Dimensions	42" x 20" x 36"	Base Weight	60 lbs.
Seat Width	17"	Battery Weight	22 lbs. (12AH) 28 lbs. (21AH)
Seat Depth	16"	Seat Weight	17 lbs.
Back Height	16"	Total Weight	99 lbs. (12AH) 105 lbs. (21AH)
Seat-to-Floor Height	19" – 21"	Warranty on Frame	Limited Lifetime
Motor	24V x 270W x 4700 rpm	Warranty on Electronics	24 Months
Controller	45A	Warranty on Batteries	12 Months