



ORDER FORM

# Squiggles + Stander

Prone, upright, and supine standing with up to 60° abduction.

Customer	
Account	<input type="text"/>
Account #	<input type="text"/> Date <input type="text"/>
<input type="checkbox"/> Quote	<input type="checkbox"/> Order PO # <input type="text"/>
Buyer	<input type="text"/>
Mark For	<input type="text"/>

Shipping Address	
Name	<input type="text"/>
Address	<input type="text"/>
City	<input type="text"/> State <input type="text"/>
ZIP	<input type="text"/> Phone <input type="text"/>
Email / Fax	<input type="text"/>

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## Squiggles+ Stander

Weight Capacity: 48.4 lbs

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	HCPCS	RETAIL
<b>Packages</b>			\$
Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, and Small Sandals with Straps.			
<input type="checkbox"/> PSPSTB	Squiggles Plus Stander Package - Blue	E0641	\$6,766
<input type="checkbox"/> PSPSTP	Squiggles Plus Stander Package - Purple	E0641	\$6,766
<input type="checkbox"/> PSPSTR	Squiggles Plus Stander Package - Red	E0641	\$6,766
<input type="checkbox"/> PSPSTRV	Squiggles Plus Stander Package - Red Vinyl	E0641	\$6,766
<b>Standard Model</b>			\$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplates and tray attachment.

<input type="checkbox"/> 184-600	Standard Model with Support Shell	\$3,173
<input type="checkbox"/> 184-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	\$1,552

### Positioning Pads \$

Selection Required. Includes chest and hip pad covers and knee support covers.

<input type="checkbox"/> 184-1800-02	Squiggles Plus Stander Positioning Pads - Blue	\$705
<input type="checkbox"/> 184-1800-04	Squiggles Plus Stander Positioning Pads - Purple	\$705
<input type="checkbox"/> 184-1800-09	Squiggles Plus Stander Positioning Pads - Red	\$705
<input type="checkbox"/> 184-1800-09VIN	Squiggles Plus Stander Positioning Pads - Red Vinyl	\$705

### Accessories \$

<input type="checkbox"/> 118-619	Tray - Gray	\$420
<input type="checkbox"/> 188-1854-05	Padded Tray Insert - Black	\$340
<input type="checkbox"/> 117-769	Grab Rail / Toy Bar for Tray	\$305
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)	\$273
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)	\$298
<input type="checkbox"/> 117-787	CITY <input type="text"/> Sandal Raiser 1" (Each) Requires Sandal Selection.	\$95

### Headrest \$

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

<input type="checkbox"/> 117-697	<b>Flat Headrest Support Hardware</b>	\$236
<input type="checkbox"/> 184-1812-02	Flat Headrest Cushions - Blue	\$125
<input type="checkbox"/> 184-1812-04	Flat Headrest Cushions - Purple	\$125
<input type="checkbox"/> 184-1812-09	Flat Headrest Cushions - Red	\$125
<input type="checkbox"/> 184-1812-09VIN	Flat Headrest Cushions - Red Vinyl (Pair)	\$125
<input type="checkbox"/> 118-667	<b>Lateral Bracket Assembly Hardware (Pair)</b>	\$217
<input type="checkbox"/> 184-1814-02	Flat Headrest and Lateral Cushions - Blue	\$218
<input type="checkbox"/> 184-1814-04	Flat Headrest and Lateral Cushions - Purple	\$218
<input type="checkbox"/> 184-1814-09	Flat Headrest and Lateral Cushions - Red	\$218
<input type="checkbox"/> 184-1814-09VIN	Flat Headrest and Lateral Cushions - Red Vinyl (Pair)	\$218

## Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Distance Between Midline of Knees	5.5" - 8.4"	140 mm - 215 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.5" - 31.5"	480 mm - 800 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 mm L x 390 mm W x 390 mm H
Pivot Base Weight	14.3 lbs	6.5 kgs
Pivot Base Dimensions (Unfolded)	37.4" L x 23.2" W x 20.3" H	950 mm L x 590 mm W x 515 mm H
Pivot Base Dimensions (Folded)	26" L x 23.2" W x 20.3" H	660 mm L x 590 mm W x 515 mm H
Pivot Base Angle Range	Prone 90-70° / Supine 90-160°	Prone 90-70° / Supine 90-160°

### Sunrise Medical (US) LLC

2842 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel.: 1-800-333-1000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-129855 Rev.F ©01/2024



800.333.4000

[www.SunriseMedical.com](http://www.SunriseMedical.com)