



ORDER FORM

# Squiggles Seat

Ultimate combination of postural seating and everyday comfort for young children while still maintaining a fun, colorful design.

Customer	
Account	<input type="text"/>
Account #	<input type="text"/> Date <input type="text"/>
<input type="checkbox"/> Quote	<input type="checkbox"/> Order PO # <input type="text"/>
Buyer	<input type="text"/>
Mark For	<input type="text"/>

Shipping Address	
Name	<input type="text"/>
Address	<input type="text"/>
City	<input type="text"/> State <input type="text"/>
ZIP	<input type="text"/> Phone <input type="text"/>
Email / Fax	<input type="text"/>

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## Squiggles Seat

Weight Capacity: 48.4 lbs

<input checked="" type="checkbox"/> OPTION	DESCRIPTION	HCPCS	RETAIL
<b>Seat Shell</b>			\$
Selection Required.			
<input type="checkbox"/> 120-610S	Squiggles Seat Shell		
	120-610-1 Seat Base with Depth Adjustment	E2294	\$1,097
	120-610-2 Seat Back - Fixed	E2293	\$1,097
	120-610-3 Seat Back Reclining Mechanism	E1014	\$479
	120-610-4 Seat Base Hip Guides / Support	E0956	\$133
<b>Seat Cover</b>			\$
Selection Required. Includes seat base and backrest cushion covers.			
<input type="checkbox"/> 120-710-03	Squiggles Seat Covers - Blue		N/C
<input type="checkbox"/> 120-710-01	Squiggles Seat Covers - Green		N/C
<input type="checkbox"/> 120-710-02	Squiggles Seat Covers - Orange		N/C
<input type="checkbox"/> 120-710-04	Squiggles Seat Covers - Pink		N/C
<b>Base</b>			\$
<input type="checkbox"/> 120-800	Size 1 Hi Low Chassis (Foot Pedal Operated.)	E1399	\$1,964
	122-710 Hi Low Chassis Push Handle		\$286
<input type="checkbox"/> 102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$2,251
	120-818 Adapter to Mount Squiggles Seat on Size 2 Hi Low Chassis		\$219
<input type="checkbox"/> 194-600	Mobility Base Interface Plate Allows mounting Squiggles Seat onto Sunrise Wheelchairs	K0108	\$516
<b>Foot Positioning</b>			\$
Sandals and Footplate Raisers require Squiggles Seat Adjustable Footrest			
<input type="checkbox"/> 120-639	Squiggles Seat Adjustable Footrest	K0040	\$243
<input type="checkbox"/> 152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$273
<input type="checkbox"/> 152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$298
<input type="checkbox"/> 120-816	Qty <input type="text" value=""/> Footplate Raiser 1" (Each)	E0951, E0952	\$92
<b>Harness</b>			\$
<input type="checkbox"/> 120-844-03	Pelvic Harness Small - Blue	E0978	\$165
<input type="checkbox"/> 120-844-01	Pelvic Harness Small - Green	E0978	\$165
<input type="checkbox"/> 120-844-02	Pelvic Harness Small - Orange	E0978	\$165
<input type="checkbox"/> 120-844-04	Pelvic Harness Small - Pink	E0978	\$165
<input type="checkbox"/> 117-848-03	Pelvic Harness Medium - Blue	E0978	\$165
<input type="checkbox"/> 117-848-01	Pelvic Harness Medium - Green	E0978	\$165
<input type="checkbox"/> 117-848-04	Pelvic Harness Medium - Pink	E0978	\$165
<input type="checkbox"/> 120-789	Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$92
<input type="checkbox"/> 155-N001-L600	Butterfly Harness Small	E0960	\$406
<input type="checkbox"/> 155-N001-L601	Butterfly Harness Medium	E0960	\$406

Squiggles Seat

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Lateral Support		\$	
Cover selection required with Lateral Support Hardware			
<input type="checkbox"/> 120-855	<b>Flip-Away Lateral Support Hardware (Pair)</b>		\$345
<input type="checkbox"/> 120-838-03	Lateral Support Covers - Blue (Pair)	E0956	\$158
<input type="checkbox"/> 120-838-01	Lateral Support Covers - Green (Pair)	E0956	\$158
<input type="checkbox"/> 120-838-02	Lateral Support Covers - Orange (Pair)	E0956	\$158
<input type="checkbox"/> 120-838-04	Lateral Support Covers - Pink (Pair)	E0956	\$158
<input type="checkbox"/> 120-656	<b>Rigid Lateral Support Hardware (Pair)</b>	E1028	\$210
<input type="checkbox"/> 120-838-03	Lateral Support Covers - Blue (Pair)	E0956	\$158
<input type="checkbox"/> 120-838-01	Lateral Support Covers - Green (Pair)	E0956	\$158
<input type="checkbox"/> 120-838-02	Lateral Support Covers - Orange (Pair)	E0956	\$158
<input type="checkbox"/> 120-838-04	Lateral Support Covers - Pink (Pair)	E0956	\$158
<input type="checkbox"/> 120-954	<b>Laterals w/ Protraction Pads - Black (Pair)</b>	E0956	\$443
<input type="checkbox"/> 137-808-03	Chest Support - Blue (Requires Laterals)	E0960	\$193
<input type="checkbox"/> 137-808-01	Chest Support - Green (Requires Laterals)	E0960	\$193
<input type="checkbox"/> 137-808-02	Chest Support - Orange (Requires Laterals)	E0960	\$193
<input type="checkbox"/> 137-808-04	Chest Support - Pink (Requires Laterals)	E0960	\$193
Tray		\$	
<input type="checkbox"/> 120-672	Tray - Gray	E0950	\$410
<input type="checkbox"/> 120-688	Extra High Tray +2" Adjustment	E0950	\$519
<input type="checkbox"/> 188-1854-05	Padded Tray Insert - Black	K0108	\$386
<input type="checkbox"/> 117-769	Grab Rail / Toy Bar for Tray	K0108	\$305
Pommel		\$	
<input type="checkbox"/> 120-641	Pommel	E0957	\$116

## Squiggles Stander

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Headrest			\$
Cushion selection required with hardware.			
Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.			
<input type="checkbox"/> 120-859	<b>Contoured Headrest Support Hardware</b>	E1028	\$385
<input type="checkbox"/> 137-800-03	Contoured Headrest Cushion - Blue	E0955	\$176
<input type="checkbox"/> 137-800-01	Contoured Headrest Cushion - Green	E0955	\$176
<input type="checkbox"/> 137-800-02	Contoured Headrest Cushion - Orange	E0955	\$176
<input type="checkbox"/> 137-800-04	Contoured Headrest Cushion - Pink	E0955	\$176
<input type="checkbox"/> 117-697	<b>Flat Headrest Support Hardware</b>	E1028	\$236
<input type="checkbox"/> 117-850-03	Flat Headrest Cushion - Blue	E0955	\$121
<input type="checkbox"/> 117-850-01	Flat Headrest Cushion - Green	E0955	\$121
<input type="checkbox"/> 117-850-04	Flat Headrest Cushion - Pink	E0955	\$121
<input type="checkbox"/> 118-667	<b>Lateral Bracket Assembly Hardware (Pair)</b> Requires Flat Headrest Support Hardware and Cushion	E1028	\$217
<input type="checkbox"/> 137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	E0956	\$197
<input type="checkbox"/> 137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	E0956	\$197
<input type="checkbox"/> 137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	E0956	\$197
<input type="checkbox"/> 117-745	Removable Head Support Bracket for 1/2" Box Section	E1028	\$73
Accessories			\$
<input type="checkbox"/> 120-779-03	1" Back Spacer Pad and Flip-Away Lateral Supports - Blue	E0956, E1297	\$832
<input type="checkbox"/> 120-779-01	1" Back Spacer Pad and Flip-Away Lateral Supports - Green	E0956, E1297	\$832
<input type="checkbox"/> 120-779-04	1" Back Spacer Pad and Flip-Away Lateral Supports - Pink	E0956, E1297	\$832
<input type="checkbox"/> 120-702	<b>Shoulder Support Hardware</b>	K0108	\$281
<input type="checkbox"/> 120-828-03	Shoulder Support Cover - Blue	K0108	\$160
<input type="checkbox"/> 120-828-01	Shoulder Support Cover - Green	K0108	\$160
<input type="checkbox"/> 120-828-04	Shoulder Support Cover - Pink	K0108	\$160

## Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (Approximate)	1 - 5 years	1 - 5 years
Stature	29.5" - 43.3"	750 mm - 1100 mm
Maximum User Weight	48.4 lbs	22 kg
Seat Depth	7.5" - 11.4"	190 mm - 290 mm
Seat Width	6.3" - 10.2"	160 mm - 260 mm
Hip Width	6.3" - 8.7"	160 mm - 220 mm
Backrest Height	12.6" - 15"	320 mm - 380 mm
Backrest Height with Shoulder Support	16.5" - 18.9"	420 mm - 480 mm
Minimum Chest Width	6.3"	160 mm
Maximum Knee Width	3.5"	90 mm
Minimum Distance Between Pommel and Side of Seat	7.5"	190 mm
Minimum Seat to Footplate	4.8"	122.5 mm
Seat to Floor	11.4" - 26.4"	290 mm - 670 mm
Backrest Angle (Prone / Recline)	10° / 25°	10° / 25°
Footplate Angle (Plantar / Dorsiflexion)	10° / 10°	10° / 10°
Tray Size	18.1" x 22"	460 mm x 560 mm
Seat Weight	13.2 lbs	6 kgs
Tray Height	4.5" - 6.7"	115 mm x 170 mm
Squiggles Seat Hi-Lo Chassis Foot Operated (Excluding Push Bar)	Imperial	Metric
Weight	19.8 lbs	9 kgs
Length	25.2"	640 mm
Width	21.3"	540 mm
Minimum Dimensions for Storage	12.4"	315 mm

### Sunrise Medical (US) LLC

2842 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel: 1-800-333-4000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-129856 Rev. H ©01/2024



800.333.4000

[www.SunriseMedical.com](http://www.SunriseMedical.com)