

PRODUCT

SUMMIT STAIRLIFT EVALUATION / ORDER FORM

Quote Order

Date: _____

Sales Representative: _____

Dealer Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Fax: _____

Email: _____

Shipping Method: Pick up Common Carrier

Account Number: _____

Projected Close Date: _____

Purchase Order Number: _____

Ship To Address: _____

Contact: _____

Phone: _____

Opportunity Reference / Project Name: _____

PARTS & OPTIONS

MODEL NUMBER

SL350OD (Outdoor)

SL350AC (AC Power)

RAIL LENGTH

Standard Rail Length: 15' 6"

SLRU20 20' track upgrade

SLRU25 25' track upgrade

CUSTOM CUT RAIL

Custom cut part numbers include track

SLR16C Custom cut track less than 16'

SLR20C 16' - 20' custom cut track

SLR25C 20' - 25' custom cut track

NOTE: If track length is greater than 16', maximum lift capacity is reduced to 250 lb

OPTIONS

SL310 Key lock on lift only


SL315 Key lock on lift and call/send

SL320 Flip-up arms

SL325 21" extra wide seat

SL330 Grocery basket (mounts to side of seat)

HOW TO MEASURE



MEASUREMENTS

_____ **A TO B** Inches from top landing to bottom floor

_____ **A TO H** Inches from top landing to nearest obstruction: (i.e. a door)

RISE _____ inches **TREAD** _____ inches *NOTE: Min angle is 35°*

STAIRWELL WIDTH _____ inches

RAIL LOCATION

Left hand

Right hand

