

SunLite

Bright Light Therapy Lamp


Product #: P801 - Silver
P802 - Brown

Product UPC: 0 23601 03801 2
0 23601 03802 9

Retail Specs: 14.75" x 9" x 2.12"

Case Pack: 4 units

MSRP: \$99.99



Bright light therapy is an effective, non-invasive treatment option for:

- The Winter Blues
- Irregular sleep patterns
- Jet Lag
- Low energy

SunLite produces 10,000 LUX of glare-free white light, the maximum level recommended by experts for therapeutic use.



SunLite is the bright light therapy lamp that provides therapeutic white light at a comfortable distance. Beautiful styling and two contemporary colour themes make SunLite the bright choice for home or office.

- 10,000 LUX at 10"
- 99.3% UV-filtered for safety
- Quick change bulbs for "no tool" maintenance
 - Lightweight at 6.6 lbs for portability
- From the leaders in Bright Light Therapy lamps

SunLite

Bright Light Therapy Lamp

Product #: P801 - Silver
P802 - Brown

Product UPC: 0 23601 03801 2
0 23601 03802 9

Retail Specs: 14.75" x 9" x 2.12"

Case Pack: 4 units

MSRP: \$99.99



Product Specifications

UL & C-UL Listed

Enclosure Size: 12.4" x 8.9" x 2" / 31.5 cm x 22.5 cm x 5 cm

Weight: 6.6 lb / 3.0 kg

Boxed Weight: 7.2 lb / 3.3 kg

Electronic Ballasts: Instant start, no flicker

Electrical: 120VAC, 50/60 Hz, 1.5A

Lens Material: High-impact, glare-free polycarbonate

Diffusing Filter: 99.3% UV filter

Light Intensity: 10,000 LUX at 10" / 2,500 LUX at 20"

Light Tubes: 2 x 36 watt compact fluorescent

4000 Kelvin colour temperature, CRI - 85

Warranty: One-year limited