

# 600

## Breezy 600

**B**reezy 600 is the lightest weight feature-packed chair in its class, making it the ideal solution for foot propellers and those needing environmental access. Ten different seat height adjustments allow this chair to fit every client.

### BREEZY 600

High-Strength Lightweight



### OPTIONS



**One-Arm Drive** positive lock strut system increases propelling efficiency and is easily removed for folding. Fits Breezy 503, 510 and Breezy 600 models.



Standard, hemi and ultra-hemi configurations (shown above).

# Breezy Features & Options

- High Strength/  
Lightweight
- Lightweight



Description	Breezy 600	Breezy 510	Breezy 500	Breezy 503
<b>• FRAME</b>				
Width	14" - 20"	14" - 20"	16", 18"	16" - 18"
Depth	16", 18"	16", 18"	16"	16"
Front Seat Ht	13" - 20.5"	13" - 20.5"	17.5" & 18.5"	17.5" & 18.5"
Backrest Ht	16" - 19"	17" - 19"	17" - 19"	17" Fixed
Reclining Backrest (180°)		Opt. (w/Adj. Pivot Pt.)		
Seat Sling	Padded & Non-Padded Nylon	Padded & Non-Padded Nylon	Padded & Non-Padded Nylon	Padded & Non-Padded Nylon
Frame Colors	Forest Green, Silver Texture & Blue	Silver Texture, Blue & Chrome	Silver Texture	Silver Texture
<b>• ARMREST</b>				
Desk Length	Padded	Padded	Padded	Padded
Full Length	Padded	Padded	Padded	Padded
Fixed Non-Removable				
Dual Post Removable		Recliner Only		
Flip Back	Standard	Standard	Standard	Standard
Height Adjustable	Option	Option	Option	Option
Tall Fixed Ht (+ 2.25")				
<b>• HANGERS</b>				
Fixed Non-Removable				
60 Degree Swing-Away	Option	Option	Option	Option
70 Degree Swing-Away	Option	Option	Option	Option
Elevating Legrest	Option	Option	Option	Option
Articulating Legrest	Option	Option	Option	Option
Footplate Style	Composite, Aluminum & Angle Adjustable	Composite & Aluminum	Composite	Composite
Heel Loops	Standard	Standard	Standard	Standard
<b>• REAR AXLE</b>				
Fixed Position Axle				
Forward Offset Axle	Standard	Standard	Standard	Standard
Amputee Axle		Option	Option	Option
One Arm Drive	Option	Option	Option	Option
Quick-Release Axle	Option	Option	Option	Option
<b>• REAR WHEEL</b>				
Size	20", 22" & 24"	20", 22" & 24"	24"	24"
Tire Type	Low Profile, Full Profile, Pneumatic, Pneumatic/ Airless	Low Profile, Full Profile, Pneumatic, Pneumatic/ Airless	Low Profile	Low Profile
Handrim Options	Composite, Aluminum Painted, Aluminum Anodized, Plastic Coated & Projection	Composite, Aluminum Painted, Aluminum Anodized, Plastic Coated & Projection	Composite, Aluminum Painted & Plastic Coated	Composite, Aluminum Painted, Aluminum Anodized, Plastic Coated
<b>• WHEEL LOCKS</b>	Push-to-Lock, Pull-to-Lock	Push-to-Lock, Pull-to-Lock	Push-to-Lock, Pull-to-Lock	Push-to-Lock, Pull-to-Lock
<b>• WEIGHT CAPACITY</b>	250 lbs.	250 lbs.	250 lbs.	250 lbs.
<b>• WEIGHT OF CHAIR (LBS.)</b>	28	32	<34	34.5
<b>• ACCESSORIES</b>				
6" Wheel Lock Ext. Handles	Option	Option	Option	Option
Anti-Fold Device				
Anti-Tippers	Option	Option	Option	Option
Calf Pad Cover	Standard	Standard	Standard	Standard
Chart Pocket				
Crutch Holder	Option	Option	Option	Option
Grade Aids	Option	Option		
Telescoping IV Pole				
Oxygen Tank Holder	Option	Option	Option	Option
Oxygen Tank Holder/IV Pole				
Positioning Belt	Option	Option	Option	Option
Wheelchair Tray	Option	Option	Option	Option