



# Mighty Lite™





# Mighty Lite™

A lightweight, portable transport solution for parents and kids, Mighty Lite offers the most adjustability of any product in its class.

The fully folding frame makes the Mighty Lite exceptionally portable and easy to use.

Four optional tie-down points on every frame ensures the Mighty Lite is Transit Option-ready.

Standard features include 8" front casters and 10" rear wheels, and easy-to-clean, quick-drying upholstery available in two versatile colors.

To accommodate the availability requirements of all our customers, Mighty Lite offers both Quick-Ship and fully Configured ordering options.

## FEATURES & BENEFITS

- ↘ Economical
- ↘ Portable
- ↘ Versatile
- ↘ Lightweight
- ↘ Foldable
- ↘ Adjustable
- ↘ Transit Option ready

## SPECIFICATIONS


Frame color: Black  
 Upholstery: Navy or Raspberry  
 Adjustable Seat angle: 15° and 22°  
 Adjustable Back angle: 85°, 90°, 95°  
 Footrest position range: from 4" up to 26" in 3/4" increments

## DIMENSIONS

Frame Size	12	14	16	18
Seat/Back Width**	11.75"	13.75"	15.5"	17.5"
Seat Depth*: Minimum	11.5"	13"	13.5"	15"
Seat Depth*: Maximum	13.5"	15"	15.5"	17"
Top of Backrest to Seat	22.25"	21.25"	23.75"	23.5"
Headrest Extension	7"	7"	7"	7"



## UPHOLSTERY COLORS



\*\* Measured between armrest tubes

\* Seat depth measurements for Seat Angle of 22 degrees, and Back Angle of 90 degrees

NOTE 1: Weight for 14" wide with 10" rear wheels and fully Transit-approved is 27.5 lbs.

NOTE 2: Specifications are subject to change without notice. Please refer to owner's manual and on-chair labeling.



# Mighty Lite Accessories

Choose from a full range of accessories and positioning options to customize your Mighty Lite.

Please note, Accessories are not available with Quick-Ship option.



## FULL TORSO SUPPORT VEST

Completely adjustable, this vest provides the most support with maximum comfort. Helps prevent forward slumping and maintains midline seated position.

## UPPER EXTREMITY SUPPORT TRAY

Provides trunk stability and is the perfect place to have a meal.

## ANGLE-ADJUSTABLE FOOTPLATES

Near-infinite range of angles to support foot positions

## FOOT & ANKLE POSITIONERS

Soft support for the feet. Available in four sizes —adjustable straps provide custom fit

**PADDED HEADWINGS**  
Provides additional lateral support for the head.

**MEDIAL THIGH SUPPORT**  
Allows for better positioning of thighs.

- Our breathable ComforFit™ material is standard (black only)
- Navy or Raspberry pack-cloth is available as a no-charge option.



## CANOPY

Keep cool and dry with this removable canopy. Height adjustable, the canopy provides cover over the headrest and includes a flap-covered observation window.



## HEADREST EXTENSION

Fully adjustable in both width and height.

- Black nylon fabric is standard.
- Navy or Raspberry Cordora is available as a no-charge option.



## TRANSIT OPTION

Crash-tested assurance for tying-down chair during transportation.



To order, call toll-free 888.593.5278

MORE Accessories →



## Accessories



### ADJUSTABLE LATERAL SUPPORT & SCOLI STRAP

Fully adjustable strap provides increased trunk support and can be used to accommodate scoliosis.

- Our breathable ComforFit material is standard (black only)
- Navy or Raspberry pack-cloth fabric is available as a no-charge option.

### H-HARNESS & PADDED COVERS

Imitation lamb's wool increases comfort, while the H-shape of the harness provides additional support.

### 3-POINT POSITIONING BELT

Increases pelvic stabilization and aids in positioning.

### UTILITY BAG

Keep essential items close at hand with this roomy, zippered bag. Easily removed from your Mighty Lite chair.

### FOAM WEDGE

Modify positioning and increase stability with this wedge insert.

### WHAT IS COMFORFIT MATERIAL?

It is a unique open-cell fabric that allows air to circulate through and moisture to escape.

The dynamic property of ComforFit provides gentle, lightly cushioned support with a bit of stretch to move with the client.

Machine washable.



To order Mighty Lite and Mighty Lite Accessories, please contact us at

888.593.5278 | f 888.400.6195

2010 E Spruce Circle | Olathe, KS 66062

artgrouprehab.com