

Patient-focused Prevention and Treatment.



Synergy® Air Elite Low Air Loss Therapy
A flexible choice for low air loss and alternating pressure.

A faded background image shows a caregiver wearing white gloves gently holding the hand of a patient lying in a hospital bed. The patient's head is resting on a pillow, and the caregiver's hands are positioned to provide support and comfort.

Hill-Rom

Enhancing Outcomes for Patients and Their Caregivers.™

The Flexible Choice for Wound Prevention and Healing.

Synergy® Air Elite Low Air Loss Therapy

With both Low Air Loss and alternating pressure treatment, the Synergy® Air Elite Low Air Loss Therapy allows you, as a caregiver, to create customized therapy and prevention programs for your patients. This mattress replacement helps enable the cost conscious hospital to implement therapy and prevention programs in accordance with Joint Commission standards, without the expense of purchasing an entire new fleet of beds. The Synergy® Air Elite surface helps enhance the level of care you provide your patients while helping improve the patient safety standards in your facility.



At a glance

- Alternating pressure and low air loss mattress replacement system:
 - Ideally suited for patients in a variety of care settings
 - Polyurethane fitted cover reduces effects of shear and friction
 - Utilizes existing bed frames, reducing rental costs and storage issues
- Prevention/Treatment Stage I, II and uncomplicated Stage III / IV wounds
- Fully inflates in less than two minutes
- Air flow addresses maceration related problems and promotes a dry healing environment
- Available for rental or purchase



Synergy® Air Elite system is available for rent on several Hill-Rom® beds such as the Advance® frame and can be used in conjunction with other flat deck beds.*

Additional Features and Benefits

Performance:

- Three therapy choices:
 - Low air loss
 - Alternating pressure/pressure relief
 - Pressure relief
- Flexibility to customize therapy for prevention and healing
- Individual air cushions designed for maximum pressure distribution, reduction in shear
- Multiple function choices powered by blower technology allow mattress to perform rhythmic inflation and deflation to redistribute and reduce patient interface pressure
- Vented air cushions provide effective low air loss therapy
- Waterproof and vapor permeable ticking

Safety:

- Auto firm function provides uniform firmness to accommodate caregiver procedures
- Safety features:
 - Power failure triggers auto alarm
 - Main power switch can be shut off in an emergency
 - CPR hose disconnects for rapid deflation
 - 2.5" foam base provides support in the event of a power failure

Comfort and Convenience:

- Individual air cushions:
 - Connected to side panel with connector for easy detachment and cleaning
- Double insulated blower provides quiet operation

**Confirm compatibility with frame manufacturer.*

Partnering for Wound Prevention and Treatment.



The average cost per hospitalization for patients who develop a Stage III & IV pressure ulcers has been reported to be \$43,180.*

Specifications

Control Unit

Inflation time	under 2 minutes
Functions.....	alternating pressure, pressure relief, max inflate, low air loss
Weight.....	13.2 lbs
Alternating pressure	5-95 minutes
Max flow rate.....	1200 lpm
Rated voltage.....	110-120V
Max current.....	5 AMPS
Rated frequency.....	60 Hz
Protection type.....	Class II, Type BF
Mode of operation	continuous

Mattress

Material:

Latex-free comforter.....	polyurethane-coated taffeta nylon
Cell.....	polyurethane-laminated taffeta nylon
Base.....	high denier polyester

Dimensions (H" x W" x L'):

Inflated.....	10" x 36" x 85"
	10" x 39" x 85"
	10" x 48" x 85"
Suggested Medicare code	E0277
Product weight	33 lbs (Standard) 40 lbs (Bariatric)
Weight range.....	36" Width 600 lbs 39" Width 1,000 lbs 48" Width 1,000 lbs

*Center for Medicare & Medicaid Services Office of Public Affairs, April 14, 2008, Fact Sheet for: CMS PROPOSES ADDITIONS TO LIST OF HOSPITAL-ACQUIRED CONDITIONS FOR FISCAL YEAR 2009.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2008 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

153727 Rev. 2 12/30/08



Enhancing Outcomes for Patients and Their Caregivers.™

USA **800-638-2546**
Canada **800-445-3730**

www.hill-rom.com