

# TRU COMFORT™

## FEATURES

- Deep contoured high density molded foam seat base
- Viscoelastic total top layer
- Ergonomically contoured back cushion
- Soft vinyl incontinence cover

## BENEFITS

- Provides lateral and anti-thrust support for clients in need of mild to moderate positioning
- Designed for clients who need mild skin protection
- Helps promote good posture
- Easily cleaned



## TRU-COMFORT™

Seat Cushion:  
E2601/E2602  
General Use

Positioning  
Back Cushion:  
E2620/E2621  
Planar Back  
with Lateral  
Supports

# TRU COMFORT™ plus

## FEATURES

- Deep contoured high density molded foam seat base
- Viscoelastic total top layer
- Ergonomically contoured back cushion
- Moisture wicking four-way stretch incontinence cover

## BENEFITS

- Provides lateral and anti-thrust support for clients in need of mild to moderate positioning
- Designed for clients who need mild to moderate skin protection
- Helps promote good posture
- Reduces moisture build-up
- Promotes better circulation



## TRU-COMFORT™ PLUS

Seat Cushion:  
E2603/E2604  
Skin Protection

Positioning  
Back Cushion:  
E2620/E2621  
Planar Back with  
Lateral Supports

# SHAPE

## FEATURES

- Ergonomically contoured
- High density base foam with a soft viscoelastic upper layer
- One hand quick release
- Solid locking durable hardware
- Highly adjustable in depth, width, and back angle
- Available in standard, low and tall back heights

## BENEFITS

- Helps promote good posture
- Offers conforming comfort without allowing bottom-out
- Allows the client to maintain an independent active lifestyle
- Rigid feel on a manual wheelchair for better performance
- Allows the medical professional to adjust for unique positioning needs of the client
- Provides support for a wide range of body types



## Synergy® SHAPE






## SHAPE

Positioning  
Back Cushion:  
E2613/E2614  
Posterior

# SYNERGY®

cushions and backs



Economically Sound. Patient First.

(US) 866-800-2002 • (Canada) 888-570-1113  
Visit us on the web: [www.quantumrehab.com](http://www.quantumrehab.com)

**QUANTUM®**  
MOVING FORWARD

A Division of Pride Mobility Products® Corp.


# SIMPLICITY®

## FEATURES

- Deep contoured high density molded foam base
- Coccyx and sacral relief
- Moisture wicking four-way stretch incontinence cover

## BENEFITS

- Provides lateral support for clients in need of mild positioning
- Designed for clients who have mild skin protection needs
- Reduces moisture build-up



E2601/E2602  
General Use


# SOLUTION® 1

## FEATURES

- Deep contoured high density molded foam base
- Coccyx and sacral relief
- Silicone soft gel top layer with air-flow design
- Moisture wicking four-way stretch incontinence cover

## BENEFITS

- Provides lateral support for clients in need of moderate positioning
- Designed for clients who need mild to moderate skin protection
- Reduces shear, moisture build-up, and heat build-up
- Promotes better circulation



E2603/E2604  
Skin Protection


# SOLUTION®

## FEATURES

- Lateral and pre-ischial shelf contoured into high density molded foam base
- Twin-cell aqueous gel insert under at-risk bony areas
- Moisture wicking four-way stretch cover
- Inner liner for incontinence and to reduce shear
- Viscoelastic total top layer

## BENEFITS

- Provides lateral and anti-thrust support for clients in need of moderate positioning
- Designed for clients who need mild to moderate skin protection
- Reduces shear and moisture build-up
- Promotes better circulation



E2607/E2608  
Skin Protection & Positioning


# STRUCTURE™ 2

## FEATURES

- Pre-ischial shelf contoured into a waterproof polyethylene foam base
- Twin-cell lightweight fluid insert under at-risk bony areas built into overlay
- Moisture wicking four-way stretch cover
- Solid seat insert included
- Extra twin-cell lightweight fluid insert included for extra emersion and positioning
- Multi-piece adjustable positioning kit available

## BENEFITS

- Provides anti-thrust support for clients in need of moderate positioning
- Provides excellent stability
- Suitable for clients with changing postural needs over time
- Designed for clients who need moderate to high skin protection
- Reduces shear and moisture build-up
- Promotes better circulation



K0734/K0735  
Adjustable Skin Protection

# SPECTRUM®

## GEL FEATURES

- Lateral and pre-ischial shelf contoured into high density molded foam base
- Twin-cell light-weight fluid insert under at-risk bony areas
- Moisture wicking four-way stretch cover
- Inner liner for incontinence and to reduce shear
- Viscoelastic total top layer
- Solid seat insert included

## GEL BENEFITS

- Provides lateral and anti-thrust support for clients in need of moderate positioning
- Provides excellent stability
- Designed for clients who need moderate to high skin protection
- Reduces shear and moisture build-up
- Promotes better circulation

E2605/E2606 – Positioning


## VISCO FEATURES

- Lateral and pre-ischial shelf contoured into high density molded foam base
- Multi-layer viscoelastic foam insert under at-risk bony areas
- Moisture wicking four-way stretch cover
- Inner liner for incontinence and to reduce shear
- Viscoelastic total top layer
- Solid seat insert included

## VISCO BENEFITS

- Provides lateral and anti-thrust support for clients in need of moderate positioning
- Provides excellent stability
- Designed for clients who need moderate to high skin protection
- Reduces shear and moisture build-up
- Promotes better circulation

E2607/E2608 – Skin Protection & Positioning



# SPECTRUM® AIR COMING SOON

## FLAT FEATURES

- Multi-cell air bladder
- Unique air-flow properties of the Sil-Air silicone foam base helps keep skin cool and dry
- Pre-ischial shelf contoured into Sil-Air silicone foam base
- Waterproof, antibacterial Sil-Air base
- Precision air volume adjustment
- Compact hand pump fits in integrated carrying pouch
- Solid seat insert included

## FLAT BENEFITS

- Provides excellent stability and support
- Designed for clients who need moderate to high skin protection
- Provides anti-thrust support for clients in need of mild to moderate positioning
- Greatly reduces bacteria and odor
- Prevents bottom-out or high pressure from over-inflation
- Convenient on-the-go air pressure adjustment

K0734/K0735  
Adjustable Skin Protection


## CONTOURED FEATURES

- Multi-cell air bladder
- Unique air-flow properties of the Sil-Air silicone foam base helps keep skin cool and dry
- Lateral and pre-ischial shelf contoured into Sil-Air silicone foam base
- Waterproof, antibacterial Sil-Air base
- Precision air volume adjustment
- Compact hand pump fits in integrated carrying pouch
- Solid seat insert included

## CONTOURED BENEFITS

- Provides excellent stability and support
- Designed for clients who need moderate to high skin protection
- Provides lateral and anti-thrust support for clients in need of moderate positioning
- Greatly reduces bacteria and odor
- Prevents bottom-out or high pressure from over-inflation
- Convenient on-the-go air pressure adjustment

K0736/K0737  
Adjustable Skin Protection and Positioning



FLAT BASE

CONTOURED BASE