

Titan LTE Portable Powerchair

HCPCS: K0821, K0821




FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F




FIGURE G

- Disassembles for transport & storage
- Flip-back armrests with height, width and angle adjustments
- Interchangeable color panels in red & blue
- 300 lbs. weight capacity
- Flat-free, non-marking tires
- Bluetooth electronics for easy troubleshooting and diagnosing

PRODUCTS	Description
TITANLTE-18FS	Titan LTE 18" Folding Seat, Red & Blue
SPECIFICATIONS	Description
Anti-Tip Wheels	Standard
Base Weight	70 lbs. 92 lbs.
Batteries (pair)	(2) 12V x 20AH
Battery Weight (pair)	22 lbs.
Brakes	Electromagnetic
Casters	8" x 2"
Charger	3A Offboard
Climbing Angle	6°
Controller	Dynamic Linx 40A
Drive Wheels	9" x 3"
Freewheel Mode	Yes
Ground Clearance	2.5"
HCPCS	K0821
Heaviest Piece	39 lbs.
Max Speed	4.2 mph
Maximum Range	8 miles
Motor	200W x 3900 rpm
Product Code/SKU	TITANLTE-18FS
Product Type	Power Wheelchair
Seat Depth	18"
Seat To Floor Height	20" - 22"
Seat Weight	30 lbs.
Seat Width	18"
Steering Type	Directional Joystick

SPECIFICATIONS	Description
Turning Radius	27"
Warranty	Limited Lifetime Warranty on Frame, 24 Month Limited Warranty on Electronics, 12 Month Limited Warranty on Batteries
ACCESSORIES	Description
SF8030	Crutch Holder Set, SType, 1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8020	Rear Basket Set, LType, 1/ea, EA
PARTS	Description
LINX-JOYSTICK	Joystick,LiNX Electronics, 1/ea, EA
LDM400001-07	Controller Tube Knob, EA
LDM400001-04	Controller Tube - Square, EA
LDM400001-01	JoystickMountingScrewsACMPower,1pr, PR
LDR410011L	Left ArmrestAssembly,Titan,Intrepid,1/ea, EA
DST-013	SeatW/SeatPlate&W/OArmrests,ScoutDST,1e, EA
LDR410011R	RightArmrestAssembly,Titan,Intrepid,1/ea, EA
LDR410110	Armrest Pad - Titan,Intrepid,1/ea, EA
C07-069-00201	Seat Post, ScoutDST, 1/ea, EA
DST-016	Seat Plate, ScoutDST, 1/ea, EA
SEC9031-B	Bearings & Spacer Set for SEC9031, 1/pkg, PKG
DST-010	ForkStmBearingSet,ScoutDST/TITANLTE,1/pk, PKG
VENT-86	Shroud Knob, Ventura, 1/ea, EA
LRM412123	21AH Battery,POV, 1/ea, EA
LRM302107	20 AMP Circuit Breaker, EA
LRM302111	2ah Battery Charger, POV,1/ea, EA