





TRANSITIONS® Modular Entry Ramp Technical Specifications

MODEL	TMER 1	TMER 1.5	TMER 2	TMER 3	TMER 4	TMER 5	TMER 6
MADE OF ALUMINUM	✓	✓	✓	✓	✓	✓	✓
SLIP-RESISTANT SURFACE	✓	✓	✓	✓	✓	✓	✓
MODULAR DESIGN	✓	✓	✓	✓	✓	✓	✓
FLAP DIMS (L x W)	1.781" x 34"	1.781" x 34"	1.781" x 34"	1.781" x 34"	1.781" x 34"	1.781" x 34"	1.781" x 34"
USABLE SIZE (H x L x W)	1"x6.5"x34"	1.5"x9.5"x34"	2"x12"x34"	3"x17"x34"	4"x22.25"x34"	5"x27.5"x34"	6"x32.75"x34"
WEIGHT*	3 lbs.	5 lbs.	5.5 lbs.	8 lbs.	11 lbs.	15 lbs.	20 lbs.
WEIGHT CAPACITY (DO NOT EXCEED)	700 lbs.	700 lbs.	700 lbs.	700 lbs.	700 lbs.	700 lbs.	700 lbs.
FOR WHEELCHAIRS & SCOOTERS	✓	✓	✓	✓	✓	✓	✓
USAGE GUIDELINES (RISE)**	.75" - 1.25"	1.25" - 1.75"	1.75" - 2.25"	2.75" - 3.25"	3.75" - 4.25"	4.75" - 5.25"	5.75" - 6.25"
WARRANTY	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited	Lifetime Limited

* Weights and dimensions are approximate and subject to change.


** Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.






TRANSITIONS® Modular Entry Ramp Technical Specifications


TMER 2



TMER 3



TMER 4





TRANSITIONS® Modular Entry Ramp Technical Specifications

