

# Supernal<sup>®</sup> 5

## SLEEP SYSTEM

The Supernal<sup>®</sup> 5 Adjustable Bed System offers you the **full functionality of a hospital bed** in an elegant design that blends seamlessly with your bedroom decor while offering full-body support and wellness-focused positioning.

Standard features like **massage to reduce body tension** and various bed positioning make comfort and independence simple. At the press of a button, you'll enjoy **ease, style, and functionality** — all in one bed built to support your well-being.



Experience a bed that adapts to you — delivering comfort, independence, and restorative sleep so you wake refreshed and ready for the day.

***Dreams, Made Easy<sup>™</sup>***



# Comfort Designed with Purpose

The Supernal® 5 combines the care you need with the luxury you want helping you stay active while aging with more comfort and independence.

## ✓ Effortless adjustment:

Head, foot, height (Hi-Low) to help you find that perfect sleep positioning

## ✓ Available in: Full, Queen, Dual King sizes

## ✓ Advanced Positioning:

Help manage health conditions both at home and in the facility

- Tilt into Trendelenburg to support circulation
- Reverse Trendelenburg to make breathing easier or help if you suffer from conditions such as acid reflux
- Cardiac Chair position to improve circulation and go from lying flat to sitting upright with less strain

## ✓ Customized height: Super low height of 12":

(for bed transfers) up to 24" (ideal for caregivers to provide care) making it easier to get in and out of bed

## ✓ Sturdy design: 400 and 500 lb weight capacities

## ✓ Soothing rolling casters: let you move or lock the bed when needed

## ✓ Soothing massage: helps you unwind and reduce body tension

## ✓ Emergency power-down: function for added peace of mind

## ✓ Backlit wireless remote: makes night time adjustments easy



## Optional Add-Ons

### ✓ Ascent or Soft Touch Mattresses: Offer superior pressure relief and support

### ✓ Support Rails: Added assistance for getting in and out of bed

### ✓ Headboard: Optional laminate headboard to go against the wall

### ✓ Bamboo Rail Covers: Soft, stylish protection for extra comfort

### ✓ Battery Backup: Keeps you powered even during outages

GF Health Products, Inc. is an  
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 678.291.3190

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.transfermaster.com](http://www.transfermaster.com)



## Comfortable In Every Setting

Created with your comfort in mind, the Supernal® 5 bed brings you the benefits of an adjustable hospital bed with advanced positioning in a beautifully discreet style that blends seamlessly into your home.

## Designed for the Way You Live

Backed by more than 100 years of Graham-Field's Made in USA expertise, our sleep solutions are designed for the way you live — delivering **Dreams, Made Easy™** for every stage of active aging.

# Supernal® 5

SLEEP SYSTEM



For more information visit us online at  
[www.transfermaster.com](http://www.transfermaster.com)

**Transfer Master**  
by graham-field