

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Riggings	Product Weight	Carton Shipping Weight	Weight Limit
TR37E	18 "	8 "	39 "	19 "	18 "	17 "	16 "	33 "	26 lbs.	33 lbs.	250 lbs.
TR39E	20 "	8 "	39 "	19 "	18 "	19 "	16 "	33 "	26 lbs.	33 lbs.	250 lbs.