



Invacare® Top End® Eliminator™ OSR Racing Wheelchair

Product Description:

Designed for speed, the Invacare® Top End® Eliminator™ OSR racing wheelchair offers exceptional stability and aerodynamics. The Eliminator™ OSR racing wheelchair features a relaxed head tube angle for great cornering and maximum traction. Combined with the stiff aerodynamic T-frame design, precision track compensator system, integral wrap around fenders, aluminum side panels and rubber coated pushrims plus a choice of three cage styles to maximize personal performance, it is easy to see why it the choice of champions around the world.

About New Smart Forms:

Top End is excited to introduce a new way to order Top End products. Our new Smart Forms provide an easy and electronic way to fill out order forms, save copies, and submit them via e-mail. Please see below for the features and useful information. Please note that because of the many upgrades available on our Invacare® Top End® Eliminator™ OSR Wheelchair, the base price does not include wheels. The standard selections are set as the default for your convenience and the overall price will be adjusted accordingly. This will allow you to build your Invacare® Top End® Eliminator™ OSR Wheelchair exactly how you want it!

How to use:

- Use any Adobe Acrobat application to electronically fill out the order form
- Use the default selections to order a product with the most common configuration
- Use the **"Reset"** button located in the forms footer to reset selections to their defaults
- Use the **"Clear All" button** located in the forms footer to clear selections, including defaults. This feature allows for a blank order form to be printed and completed by hand
- Use the **"E-mail to Top End order box"** button located on the last page of the form to e-mail the completed form directly to the topendorders@invacare.com e-mail inbox
- Forms can be completed electronically, printed, and e-mailed or faxed to topendorders@invacare.com or 727-522-1007 if access to Acrobat Reader for Windows is unavailable

(Acrobat Reader for Windows with java-script enabled is required to save and attach to e-mail)

Features:

- Tabulated options for an organized and standardized work-flow
- New more specific option codes which will appear on confirmations to improve customer / manufacturing communication
- Common special requests have been given option codes
- Built in exclusions to help eliminate errors
- Ability to save the form after it is filled out
- Ability to e-mail the form inside Acrobat Reader rather than saving a copy and attaching to an e-mail
- Ability to clear the entire form so it can be printed and completed by hand in the field when a computer is not available
- Ability to reset form to default selections



Invacare® Top End® Eliminator™ OSR Racing Wheelchair



Order Request Type:			
☐ Quote	☐ CAD/Quote	☐ Order	PO #
Provider / Customer Information:			
Provider Information:		End User Information:	
Date:	Name:	Height:	
IVC Acct #:	Ship to Address:	Weight:	
Contact:	Address Line 2:	Leg Length:	
Phone #:	City:	Hip Width:	
Fax/Email:	State:	Zip:	Phone #:
☐ Ship to Provider	☐ Ship to End User		Email

Complete chair starts at: \$3,395.00

Frame: (see definitions)

Item #:	Description:	Price:
☐ 10015	Eliminator OSR frame - Open V cage	\$3,130
☐ 10020	Eliminator OSR frame - U cage	\$3,130
☐ 10021	Eliminator OSR frame - I cage	\$3,130

Camber Angle:

Item #:	Description:	Price:
☐ 45005	11 degree camber	\$0
☐ 45006	12 degree camber	\$0
☐ 45007	13 degree camber	\$0
☐ 45008	15 degree camber	\$0

Steering: (see diagram 4)

Item #:	Description:	Price:
Steering Handle:		
☐ 65120	Horizontal Steering Handle	\$0
☐ 65121	Vertical Steering Handle	\$0
Track Control:		
☐ 65115	Triangle Track Control	\$0
☐ 65116	Oval Track Control	\$50
Track Control Mount:		
☐ 65110	Top Mount Steering	\$0
☐ 65111	Bottom Mount Steering	\$0

Frame Options:

Item #:	Description:	Price:
☐ 75010	Aluminum Solid Seat (see definitions)	\$75
☐ 75015	2" Cut-out Axle (see definitions)	\$250
☐ 75016	3" Cut-out Axle (see definitions)	\$250
☐ 75017	4" Cut-out Axle (see definitions)	\$250

Sizing for all cage styles:

Item #:	Description:	Price:
Upper Frame Width:		
☐ 15195	10.0" Upper frame width	\$0
☐ 15196	10.5" Upper Frame Width	\$0
☐ 15197	11.0" Upper Frame Width	\$0
☐ 15198	11.5" Upper Frame Width	\$0
☐ 15199	12.0" Upper Frame Width	\$0
☐ 15200	12.5" Upper Frame Width	\$0
☐ 15201	13.0" Upper Frame Width	\$0
☐ 15202	13.5" Upper Frame Width	\$0
☐ 15203	14.0" Upper Frame Width	\$0
☐ 15204	14.5" Upper Frame Width	\$0
☐ 15205	15.0" Upper Frame Width	\$0
☐ 15206	15.5" Upper Frame Width	\$0
☐ 15207	16.0" Upper Frame Width	\$0
☐ 15208	16.5" Upper Frame Width	\$0
☐ 15209	17.0" Upper Frame Width	\$0
Lower Frame Width:		
☐ 15217	11.0" Lower frame width	\$0
☐ 15218	11.5" Lower frame width	\$0
☐ 15219	12.0" Lower frame width	\$0
☐ 15220	12.5" Lower frame width	\$0
☐ 15221	13.0" Lower frame width	\$0
☐ 15222	13.5" Lower frame width	\$0
☐ 15223	14.0" Lower frame width	\$0
☐ 15224	14.5" Lower frame width	\$0
☐ 15225	15.0" Lower frame width	\$0
☐ 15226	15.5" Lower frame width	\$0
☐ 15227	16.0" Lower frame width	\$0
☐ 15228	16.5" Lower frame width	\$0
☐ 15229	17.0" Lower frame width	\$0
☐ 15230	17.5" Lower frame width	\$0
☐ 15231	18.0" Lower frame width	\$0



Invacare® Top End® Eliminator™ OSR Racing Wheelchair



Sizing for all cage styles:

Item #:	Description:	Price:
Rear Seat Height:		
☐ 15280	15.0" Rear Seat Height	\$0
☐ 15281	15.5" Rear Seat Height	\$0
☐ 15282	16.0" Rear Seat Height	\$0
☐ 15283	16.5" Rear Seat Height	\$0
☐ 15284	17.0" Rear Seat Height	\$0
☐ 15285	17.5" Rear Seat Height	\$0
☐ 15286	18.0" Rear Seat Height	\$0
☐ 15287	18.5" Rear Seat Height	\$0
☐ 15288	19.0" Rear Seat Height	\$0
☐ 15289	19.5" Rear Seat Height	\$0
☐ 15290	20.0" Rear Seat Height	\$0
Overall Length:		
☐ 15245	64" Overall Length	\$0
☐ 15247	66" Overall Length	\$0
☐ 15249	68" Overall Length	\$0
☐ 15251	70" Overall Length	\$0
☐ 15253	72" Overall Length	\$0
☐ 15255	74" Overall Length	\$0
☐ 15257	76" Overall Length	\$0
Axle Position:		
☐ 15236	5.0" Axle Position	\$0
☐ 15237	5.5" Axle Position	\$0
☐ 15238	6.0" Axle Position	\$0
☐ 15239	6.5" Axle Position	\$0
☐ 15240	7.0" Axle Position	\$0
☐ 15241	7.5" Axle Position	\$0
☐ 15242	8.0" Axle Position	\$0

Addition Sizing for Open V and I cage styles:

Item #:	Description:	Price:
Front Seat Height:		
☐ 15295	16.0" Front Seat Height	\$0
☐ 15296	16.5" Front Seat Height	\$0
☐ 15297	17.0" Front Seat Height	\$0
☐ 15298	17.5" Front Seat Height	\$0
☐ 15299	18.0" Front Seat Height	\$0
☐ 15300	18.5" Front Seat Height	\$0
☐ 15301	19.0" Front Seat Height	\$0
☐ 15302	19.5" Front Seat Height	\$0
☐ 15303	20.0" Front Seat Height	\$0

Additional Sizing for Open V cage:

Item #:	Description:	Price:
Cage depth:		
☐ 15265	15" Cage Depth	\$0
V depth:		
☐ 15260	24" V Depth	\$0
☐ 15261	25" V Depth	\$0
☐ 15262	26" V Depth	\$0

Additional Sizing for I cage:

Item #:	Description:	Price:
Depth:		
☐ 15270	15" Depth	\$0

Sizing for U cage style:

Item #:	Description:	Price:
Knee Height:		
☐ 15305	16" Knee Height	\$0
☐ 15306	17" Knee Height	\$0
☐ 15307	18" Knee Height	\$0
☐ 15308	19" Knee Height	\$0
☐ 15309	20" Knee Height	\$0
Knee depth:		
☐ 15275	18" Knee Depth	\$0
Knee Width:		
☐ 15315	7.0" Knee Width	\$0
☐ 15316	7.5" Knee Width	\$0
☐ 15317	8.0" Knee Width	\$0
☐ 15318	8.5" Knee Width	\$0
☐ 15319	9.0" Knee Width	\$0
☐ 15320	9.5" Knee Width	\$0
☐ 15321	10.0" Knee Width	\$0
☐ 15322	10.5" Knee Width	\$0
☐ 15323	11.0" Knee Width	\$0
☐ 15324	11.5" Knee Width	\$0
☐ 15325	12.0" Knee Width	\$0

Wheel Options:

Item #:	Description:	Tire:	Price:
Front wheels:			
☐ 50146	20" Corima Carbon Fiber	Tubular	\$600
☐ 50145	20" Metal Spoked	Tubular	\$0
☐ 50150	Omit Front Wheel	N/A	\$0
Rear wheels:			
☐ 50139	26" High Performance	Clincher	\$265
☐ 50147	700C Corima Carbon Fiber Disk	Tubular	\$2,515
☐ 50148	700C Corima Carbon Fiber 4 spoke	Tubular	\$2,515
☐ 50149	700C Zipp Carbon Fiber Disk	Tubular	\$2,515
☐ 50143	700C High Performance	Clincher	\$265
☐ 50151	Omit Rear Wheels	N/A	\$0

Sizing:

Item #:	Description:	Price:
Rubber Coated Pushrim Diameter:		
☐ 50330	12" Diameter Pushrims	\$0
☐ 50332	13" Diameter Pushrims	\$0
☐ 50334	14" Diameter Pushrims	\$0
☐ 50336	15" Diameter Pushrims	\$0
☐ 50338	16" Diameter Pushrims	\$0

Upholstery:

Item #:	Description:	Price:
Material/Color:		
☒ 30001	Nylon, Black Upholstery	\$0
☐ 30002	Nylon, Woodland Camo Upholstery	\$100
☐ 30003	Nylon, Digital Camo Upholstery	\$100
☐ 30018	Seat Cushion	\$75
Seat Style:		
☐ 30016	6" Seat Upholstery Length	\$0
☐ 30017	9" Amputee Seat Upholstery Length	\$0



Invacare® Top End® Eliminator™ OSR Racing Wheelchair

Options and Accessories:		
Item #:	Description:	Price:
Safety:		
☐ 25015	Safety Light, Front (White)	\$25
☐ 25016	Safety Light, Rear (Red)	\$25
☐ 25073	Safety Flag	\$30
☐ 25018	Helmet, Small/Medium	\$75
☐ 25019	Helmet, Medium/Large	\$75
Positioning:		
☐ 60014	Clamp-on Footrest, Adjustable	\$200
☐ 25068	Additional 2" Strap	\$20
☐ 25001	Click Strap set, 8" Small	\$75
☐ 25002	Click Strap set, 10" Medium	\$80
☐ 25003	Click Strap set, 12" Large	\$86
Performance Enhancement:		
☐ 25022	Wheel Alignment Gauge	\$75
☐ 25024	Cordless Computer	\$125
Transportation:		
☐ 25007	Bike Rack, Hitch 1-1/4"	\$450
☐ 25008	Bike Rack, Hitch 2"	\$450
☐ 25031	Wheel bag	\$180
Training:		
☐ 25033	Heartrate Monitor	\$100
☐ 25069	Indoor Racing Roller	\$800
Apparel:		
☐ 25009	T-Shirt, Large	\$10
☐ 25010	T-Shirt, Extra Large	\$10
☐ 25026	Racing Jersey, Large	\$70
☐ 25027	Racing Jersey, Extra Large	\$70

Color Options:			
Item #:	Color:	Hue:	Price:
Frame Color:			
☐ 55002	Red	Tangerine Red	\$0
☐ 55003	Red	Raspberry Red	\$0
☐ 55005	Orange	Simple Orange	\$0
☐ 55007	Yellow	Safety Yellow	\$0
☐ 55008	Green	Grasshopper Green	\$0
☐ 55010	Blue	Midnight Blue	\$0
☐ 55011	Blue	Royal Blue	\$0
☐ 55012	Blue	Electric Blueberry	\$0
☐ 55014	Purple	Grape Madness	\$0
☐ 55015	Pink	Bubble-Gum Pink	\$0
☐ 55016	White	White Snow	\$0
☐ 55017	Black	Stealth Black	\$0
☐ 55018	Black	Wet Black	\$0
☐ 55019	Black	Black Prism	\$0
☐ 55020	Silver	Silver Metallic	\$0
☐ 55023	Grey	Stealth Grey	\$0

Optional Alternative Fork Color:			
Item #:	Color:	Hue:	Price:
Fork Color: (Pick one only when two-tone is chosen)			
☐ 55025	Red	Tangerine Red	\$100
☐ 55026	Red	Raspberry Red	\$100
☐ 55028	Orange	Simple Orange	\$100
☐ 55030	Yellow	Safety Yellow	\$100
☐ 55031	Green	Grasshopper Green	\$100
☐ 55033	Blue	Midnight Blue	\$100
☐ 55034	Blue	Royal Blue	\$100
☐ 55035	Blue	Electric Blueberry	\$100
☐ 55037	Purple	Grape Madness	\$100
☐ 55038	Pink	Bubble-Gum Pink	\$100
☐ 55039	White	White Snow	\$100
☐ 55040	Black	Stealth Black	\$100
☐ 55041	Black	Wet Black	\$100
☐ 55042	Black	Black Prism	\$100
☐ 55043	Silver	Silver Metallic	\$100
☐ 55046	Grey	Stealth Grey	\$100

Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	

E-mail to Top End order box

Invacare Corporation Top End Division

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

©2015 Invacare Corporation

PN: 1181415 REV A KPIPARG 07-01-2016

[Clear All](#) [Reset Form](#)



Invacare® Top End® Eliminator™ OSR Racing Wheelchair

Definitions:	
Open V cage:	Recommended for beginner level athletes with T-10 SCI and above as that level of injury requires more of a bucket. Designed for people who prefer the seated position and/or may also want to kneel. This cage does require good transfer skills.
U cage:	Recommended for experienced customers who exclusively use the kneeling position. Customers tend to have good balance and are T-10 injury level or below. U cage includes an aluminum knee plate for extra rigidity.
I cage:	I Cage is the easiest transfer style to use but cannot be converted to a kneeling position. This style is typically used by higher level injury athletes who sit in the traditional feet
Aluminum solid seat:	Provides a solid surface. Not recommended if customers intend to adjust the seat upholstery
Cut out axle:	Recommended when a customer who desires a kneeling position requires additional space for their lower legs. This additional space is accomplished by removing the center portion of the camber tube after welding. Creates more room for legs with low seat heights.

