

Top End® Excelerator™ Handcycle

PROVIDER / CUSTOMER INFORMATION:		
Provider Information:		End User Information:
Date:	Contact Name:	Athlete Name:
Acct #:	Ship to Address:	
Contact:	City:	State:
Phone #:	Zip/Postal Code:	Phone:
Fax/Email:	Email:	
<input type="radio"/> Ship to Provider		<input checked="" type="radio"/> Ship to End User

ORDER REQUEST TYPE:		
<input type="radio"/> Quote	<input type="radio"/> Order	PO # <input type="text"/>

BASE MODEL:		
<input checked="" type="radio"/> 10005	Top End Excelerator Handcycle	\$4,899

Assembly: <i>(see definitions)</i>		
<input checked="" type="radio"/> 10048	Assembled	\$0
Shipping:		
Shipping Fee <i>(Continental US Only - See definitions)</i>		\$325

SEAT WIDTH:		
<input type="radio"/> 15017	16" Seat Width	\$0
<input type="radio"/> 15021	18" Seat Width	\$0
<input type="radio"/> 15025	20" Seat Width	\$0

CAMBER ANGLE:		
<input checked="" type="radio"/> 45002	3 Degree Camber	\$0

WHEELBASE: <i>(Front to Rear Axle)</i>		
<input checked="" type="radio"/> 75085	40" Wheelbase	\$0

WHEEL SETS:		
<input checked="" type="radio"/> 50121	24" Front Wheels & 24" Rear Wheels with Quick Release Axles	\$0

TIRES:		
<input checked="" type="radio"/> 50213	24" X 1.75" Treaded Tires	\$0

FOOTREST:		
<input checked="" type="radio"/> 60002	Dual, Height/Angle Adjustable	\$0



DRIVE CONTROL:		
HAND PEDALS: <i>(see definitions)</i>		
<input type="radio"/> 65137	Medium 3D Printed with Power Plates	\$0
<input type="radio"/> 65134	Tri-Pin Hand Pedals	\$361
SHIFTING:		
<input checked="" type="radio"/> 65060	Twist Activated Shifter with Internal Gearing	\$0
<input type="radio"/> 65061	Quad Twist Shifting Adaptor <i>(see definitions)</i>	\$181
BRAKING:		
<input checked="" type="radio"/> 65080	Reverse Activated Drum Brake	\$0
<input checked="" type="radio"/> 65090	Locking Parking Brake	\$0
CHAINRING:		
<input checked="" type="radio"/> 65100	Single 38 Tooth Chainring	\$0
CRANK TYPE: <i>(see definitions)</i>		
<input type="radio"/> 65054	Crank Arms: 175mm Long x 12" Wide	\$0
<input type="radio"/> 65141	V - Crankset	\$360

Top End Sports LLC

All rights reserved. Trademarks are identified by the symbols "™" and "®". All trademarks are owned by or licensed to Top End unless otherwise noted. Specifications and prices are subject to change without notification.
©2023 Top End Sports LLC

Top End®

Excelerator™ Handcycle

UPHOLSTERY:		
MATERIAL & COLOR:		
<input checked="" type="radio"/>	30001	Nylon, Black Upholstery \$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/>	30014	Adjustable Tension Seat Upholstery \$0
<input checked="" type="radio"/>	25060	Hook and Loop Positioning Strap \$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/>	30029	Fold - Over Back Upholstery \$0
SEAT CUSHION:		
<input checked="" type="radio"/>	25063	Seat Cushion \$0

FRAME COLOR:			
<input type="radio"/>	55018	Black	Wet Black \$0
<input type="radio"/>	55019	Black	Black Prism \$0
<input type="radio"/>	55017	Black	Stealth Black \$0
<input type="radio"/>	55010	Blue	Midnight Blue \$0
<input type="radio"/>	55011	Blue	Royal Blue \$0
<input type="radio"/>	55012	Blue	Electric Blueberry \$0
<input type="radio"/>	55023	Grey	Stealth Grey \$0
<input type="radio"/>	55053	Green	Agriculture Green \$0
<input type="radio"/>	55008	Green	Grasshopper Green \$0
<input type="radio"/>	55054	Green	Neon Green \$0
<input type="radio"/>	55006	Orange	Sunkissed Orange \$0
<input type="radio"/>	55005	Orange	Simple Orange \$0
<input type="radio"/>	55051	Orange	Neon Orange \$0
<input type="radio"/>	55015	Pink	Bubblegum Pink \$0
<input type="radio"/>	55057	Pink	Neon Pink \$0
<input type="radio"/>	55055	Purple	Lilac Purple \$0
<input type="radio"/>	55001	Red	Simple Red \$0
<input type="radio"/>	55002	Red	Tangerine Red \$0
<input type="radio"/>	55003	Red	Raspberry Red \$0
<input type="radio"/>	55050	Red	Neon Red \$0
<input type="radio"/>	55020	Silver	Silver Metallic \$0
<input type="radio"/>	55016	White	White Snow \$0
<input type="radio"/>	55007	Yellow	Safety Yellow \$0
<input type="radio"/>	55052	Yellow	Neon Yellow \$0
<input type="radio"/>	55070	Raw	Buffed Aluminum \$0
<input type="radio"/>	55090	Bronze	Las Vegas Bronze \$0
<input type="radio"/>	55091	Green	Mint Green \$0
<input type="radio"/>	55092	Pink	Light Pink \$0
<input type="radio"/>	55093	Grey	Fighter Jet Grey \$0
<input type="radio"/>	55094	Blue	Tiffany Blue \$0
<input type="radio"/>	55095	Blue	Carolina Blue \$0

ACCESSORIES:		
STANDARD:		
<input checked="" type="radio"/>	25061	Safety Flag \$0
<input checked="" type="radio"/>	25062	Chain Guard \$0
UPGRADABLE:		
SAFETY:		
<input type="radio"/>	25097	Mirror, Rear View \$37
HYDRATION:		
<input checked="" type="radio"/>	25021	Water Bottle and Cage \$60



Note: All Special requests are subject to approval, additional charges, and additional lead time.

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	

Top End Sports LLC

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End unless otherwise noted.
Specifications and prices are subject to change without notification.
©2023 Top End Sports LLC.

Top End® Excelerator™ Handcycle

HANDCYCLE SPECIFICATIONS:	
Weight:	50 pounds
Weight Capacity:	350 pounds
Turning Radius:	12 feet for 360 degree circle, 8 feet on a 3-point turn
Transfer Height:	20" from ground with 2" cushion
Frame Style:	Single main tube, upright position, 1-3/4" X .065" ovalized carbon steel tubing
Fork:	Carbon Steel
Propulsion / Shifting:	Internal front hub, Twist activated shifter mounted on crank assembly
Cranks:	165mm X 12" Adjustable height
Back:	Adjustable height from 14"-18"
Footrest:	Dual adjustable height and angle, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Reverse Activated Drum and Locking Parking Brake mounted on fork
Wheelbase:	40" Axle to Axle
Camber:	3 or 9 degree available
Wheels:	24" metal spoked wheels
Axles:	Quick Release
Tires / Tubes:	24" X 1.75" Treaded Tires (45 PSI)
Positioning Straps:	Hook and loop positioning strap included
Cushions:	Seat cushion included
Safety / Visibility:	Flag and reflectors included
Chain Guard:	Full length shock resistant composite material except when Mountain Drive is chosen. With Mountain Drive selection a half moon shaped black composite chainguard is used.

DEFINITIONS:	
Assembly and Shipping:	Arrives mostly assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Quad Twist Shifter Adaptor:	Three pronged adaptor which allows user to tap for shifting.
Tri-pin Hand Pedals:	3-Point contact hand pedal similar to hand controls used for vehicles. Allows for the hands to be kept in place for people with limited grip strength.
V - Crankset:	Option for people who need more clearance due to long legs or large midsections. Cranks are 180mm long and 20" wide.
Accessories:	All Add-on accessories may be boxed separately and require self installation. Additionally, all add-on options may be shipped separately.

Warranty Info:

Top End warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty.