

Top End® Force™ RX Handcycle

PROVIDER / CUSTOMER INFORMATION:

Provider Information:		End User Information:	
Date:	Contact Name:	Athlete Name:	
Acct #:	Ship to Address:		
Contact:	City:	State:	
Phone #:	Zip/Postal Code:	Phone:	
Fax/Email:	Email:		
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	

ORDER REQUEST TYPE:

☐ Quote
 ☒ Order
 PO #

BASE MODEL:

<input checked="" type="radio"/> 10009	Force RX Handcycle	\$9,499
--	--------------------	---------

Assembly: (see definitions)

<input checked="" type="radio"/> 10048	Assembled	\$0
--	-----------	-----

Shipping:

Shipping Fee (Continental US Only - See definitions)	\$325
---	-------

WHEEL SETS:

WHEEL SIZES:

<input checked="" type="radio"/> 50171	Front, 26" and Rear, 20"	\$0
--	--------------------------	-----

26" FRONT and 20" REAR:

<input type="radio"/> 50176	High Performance, 26" Front High Performance, 20" Rear (Clincher)	\$0
<input type="radio"/> 50177	C7 Carbon, 26" Front C7 Carbon, 20" Rear (Clincher)	\$1,530




WHEELBASE: (see definitions)

<input checked="" type="radio"/> 75014	Standard Length	\$0
--	-----------------	-----

SEAT WIDTH: (see diagram)

<input type="radio"/> 15011	13" Seat Width	\$0
<input type="radio"/> 15013	14" Seat Width	\$0
<input type="radio"/> 15015	15" Seat Width	\$0
<input type="radio"/> 15017	16" Seat Width	\$0
<input type="radio"/> 15019	17" Seat Width	\$0
<input type="radio"/> 15021	18" Seat Width	\$0



CAMBER ANGLE:

<input checked="" type="radio"/> 45002	3 Degree Camber	\$0
--	-----------------	-----


Top End Sports LLC.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End Sports LLC unless otherwise noted. Specifications and prices are subject to change without notification.

©2023 Top End Sports LLC

Top End® Force™ RX Handcycle

FORK OPTIONS:		
FORK STYLE:		
<input type="radio"/> 65181	Adjustable Height Fixed Length <i>(see definitions)</i>	\$0
<input type="radio"/> 65180	Adjustable Height Adjustable Length	\$361
LOWER FORK STYLE:		
<input type="radio"/> 75025	Aerodynamic	\$0
<input type="radio"/> 75026	Wide <i>Required with Disk Brake</i>	\$0



BRAKING:		
<input type="radio"/> 65084	Rim Brakes, Tektro <i>(see definitions)</i>	\$0
<input type="radio"/> 65082	Front disk brake <i>Must select Wide Lower Fork</i>	\$361

CHAINRINGS:		
<input type="radio"/> 65096	Double, 48-34T <i>Not available with eTap shifting package</i>	\$0
<input type="radio"/> 65118	Single, 36T <i>Required for eTap shifting Not available with 2x Rival shifting</i>	\$0

SHIFTER OPTIONS: <i>(see diagram)</i>		
<input type="radio"/> 65095	Shimano 105 2x Trigger Shifting Package <i>Shifter & brake mounted on right pedal only</i> Shimano 105 Shifter Shimano 105 Derailleurs Avid Brake Lever 11-34T 11 Speed Cassette Cable Driven chainring derailleur	\$0
<input type="radio"/> 65088	SRAM eTap 1x Eagle GX Electronic Shifting Package <i>Shifter & brake mounted on right pedal</i> SRAM eTap Shifter SRAM eTap Eagle GX 12 Speed Derailleur SRAM Time Trail Brake Lever 11-50T 12 speed cassette	\$1,415
<input type="radio"/> 65094	SRAM GX 1x Trigger Shifting Package <i>Shifter & brake mounted on right pedal only</i> SRAM GX Shifter SRAM GX Cable 12 Speed Derailleur Avid Brake Lever 11-50T 12 speed Cassette	\$0


Speed Chart for Shifting Packages				
Shifting Package TE65094				
	Cadence (Pedal RPM)			
Gear Combination	70	80	90	100
34x11	15.6	17.9	20.1	22.3
34x34	5.1	5.8	6.5	7.2
48x11	22.1	25.2	28.4	31.5
48x34	7.1	8.2	9.2	10.2

Shifting Package TE65094 and TE65088				
	Cadence (Pedal RPM)			
Gear Combination	70.0	80.0	90.0	100.0
36x11	16.5	18.9	21.3	23.6
36x50	3.6	4.2	4.7	5.2

Top End®

Force™ RX Handcycle


CRANKS: (see diagram below)		
CRANK STYLE:		
<input checked="" type="radio"/> 65055	V Cranks	\$0
CRANK LENGTH:		
<input type="radio"/> 65037	165mm Long	\$0
<input type="radio"/> 65146	Custom Crank Length	\$54
CRANK WIDTH:		
<input type="radio"/> 65009	14" Wide <i>N/A with 25 degree handpedals</i>	\$0
<input type="radio"/> 65147	Custom Crank Width	\$54



HAND PEDALS:		
SIZE:		
<input checked="" type="radio"/> 65137	Medium (1-3/4") 3D Printed W/Power Plates <i>(15 Degree only)</i>	\$0



FOOTREST OPTIONS: (see diagram below)		
<input type="radio"/> 60008	Standard Footrests, 3-3/4" Wide	\$0
<input type="radio"/> 60009	Standard Footrests, 4-3/4" Wide	\$0



The diagram shows a side view of a footrest with a double-strap design. A vertical arrow indicates the 'Width' of the footrest at its widest point.

UPHOLSTERY:		
MATERIAL:		
<input checked="" type="radio"/> 30001	Nylon, Black Upholstery	\$0
BACK UPHOLSTERY:		
<input checked="" type="radio"/> 25064	Back Cushion	\$0
SEAT UPHOLSTERY:		
<input checked="" type="radio"/> 25063	Seat Cushion	\$0
<input checked="" type="radio"/> 25060	Hook and Loop Positioning Strap	\$0

FRAME COLOR:			
<input type="radio"/> 55018	Black	Wet Black	\$0
<input type="radio"/> 55019	Black	Black Prism	\$0
<input type="radio"/> 55017	Black	Stealth Black	\$0
<input type="radio"/> 55010	Blue	Midnight Blue	\$0
<input type="radio"/> 55011	Blue	Royal Blue	\$0
<input type="radio"/> 55012	Blue	Electric Blueberry	\$0
<input type="radio"/> 55023	Grey	Stealth Grey	\$0
<input type="radio"/> 55053	Green	Agriculture Green	\$0
<input type="radio"/> 55008	Green	Grasshopper Green	\$0
<input type="radio"/> 55054	Green	Neon Green	\$0
<input type="radio"/> 55006	Orange	Sunkissed Orange	\$0
<input type="radio"/> 55005	Orange	Simple Orange	\$0
<input type="radio"/> 55051	Orange	Neon Orange	\$0
<input type="radio"/> 55015	Pink	Bubblegum Pink	\$0
<input type="radio"/> 55057	Pink	Neon Pink	\$0
<input type="radio"/> 55055	Purple	Lilac Purple	\$0
<input type="radio"/> 55001	Red	Simple Red	\$0
<input type="radio"/> 55002	Red	Tangerine Red	\$0
<input type="radio"/> 55003	Red	Raspberry Red	\$0
<input type="radio"/> 55050	Red	Neon Red	\$0
<input type="radio"/> 55020	Silver	Silver Metallic	\$0
<input type="radio"/> 55016	White	White Snow	\$0
<input type="radio"/> 55007	Yellow	Safety Yellow	\$0
<input type="radio"/> 55052	Yellow	Neon Yellow	\$0
<input type="radio"/> 55070	Raw	Buffed Aluminum	\$0
<input type="radio"/> 55090	Bronze	Las Vegas Bronze	\$0
<input type="radio"/> 55091	Green	Mint Green	\$0
<input type="radio"/> 55092	Pink	Light Pink	\$0
<input type="radio"/> 55093	Grey	Fighter Jet Grey	\$0
<input type="radio"/> 55094	Blue	Tiffany Blue	\$0
<input type="radio"/> 55095	Blue	Carolina Blue	\$0

Top End Sports LLC.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End Sports LLC unless otherwise noted. Specifications and prices are subject to change without notification.

©2023 Top End Sports LLC

1181419~N


Top End®

Force™ RX Handcycle

FORK COLOR:			
<input type="radio"/>		Same as Frame Color	
<input type="radio"/> 55067	N/A	Not applicable with 55066 (buffed frame), must select 55044 (buffed fork)	\$0
<input type="radio"/> 55041	Black	Wet Black	\$120
<input type="radio"/> 55042	Black	Black Prism	\$120
<input type="radio"/> 55040	Black	Stealth Black	\$120
<input type="radio"/> 55033	Blue	Midnight Blue	\$120
<input type="radio"/> 55034	Blue	Royal Blue	\$120
<input type="radio"/> 55035	Blue	Electric Blueberry	\$120
<input type="radio"/> 55046	Gray	Stealth Grey	\$120
<input type="radio"/> 55061	Green	Agriculture Green	\$120
<input type="radio"/> 55031	Green	Grasshopper Green	\$120
<input type="radio"/> 55062	Green	Neon Green	\$120
<input type="radio"/> 55029	Orange	Sunkissed Orange	\$120
<input type="radio"/> 55028	Orange	Simple Orange	\$120
<input type="radio"/> 55059	Orange	Neon Orange	\$120
<input type="radio"/> 55038	Pink	Bubblegum Pink	\$120
<input type="radio"/> 55065	Pink	Neon Pink	\$120
<input type="radio"/> 55063	Purple	Lilac Purple	\$120
<input type="radio"/> 55024	Red	Simple Red	\$120
<input type="radio"/> 55025	Red	Tangerine Red	\$120
<input type="radio"/> 55026	Red	Raspberry Red	\$120
<input type="radio"/> 55058	Red	Neon Red	\$120
<input type="radio"/> 55043	Silver	Silver Metallic	\$120
<input type="radio"/> 55039	White	White Snow	\$120
<input type="radio"/> 55030	Yellow	Safety Yellow	\$120
<input type="radio"/> 55060	Yellow	Neon Yellow	\$120
<input type="radio"/> 55096	Bronze	Las Vegas Bronze	\$120
<input type="radio"/> 55097	Green	Mint Green	\$120
<input type="radio"/> 55098	Pink	Light Pink	\$120
<input type="radio"/> 55099	Grey	Fighter Jet Grey	\$120
<input type="radio"/> 55100	Blue	Tiffany Blue	\$120
<input type="radio"/> 55101	Blue	Carolina Blue	\$120
<input type="radio"/> 55044	Raw	Buffed Aluminum	\$241

ACCESSORIES:		
STANDARD:		
<input checked="" type="radio"/> 25061	Safety Flag	\$0
<input checked="" type="radio"/> 25067	Drafting Bumper	\$0
<input checked="" type="radio"/> 25071	Adjustable Headrest	\$0
<input checked="" type="radio"/> 25038	Mirror, Chain Guard Mounted	\$0
UPGRADABLE:		
<input type="radio"/> 25039	Leg Guard	\$90
<input type="radio"/> 60015	Leg Support, Curved (Pair)	\$157
HYDRATION:		
<input type="radio"/> 25021	Water Bottle and Cage	\$0

DECAL COLORS: (decals are reflective)		
FRAME DECAL COLOR:		
<input type="radio"/> 55071	Black Frame Decals	\$0
<input type="radio"/> 55072	White Frame Decals	\$0
FORK DECAL COLOR:		
<input type="radio"/> 55073	Black Fork Decals	\$0
<input type="radio"/> 55074	White Fork Decals	\$0
FORK COLOR		
<input type="radio"/> 55106	Fork Color: Bomber Orange	\$120
<input type="radio"/> 55076	Fork Color: Custom	\$224
<input type="radio"/> 55068	Fork Color: Flat Black	\$0
<input type="radio"/> 55104	Fork Color: Wrinkle Black	\$0
FRAME COLOR		
<input type="radio"/> 55105	Frame Color: Bomber Orange	\$0
<input type="radio"/> 55066	Frame Color: Buffed Aluminum	\$421
<input type="radio"/> 55102	Frame Color: Flat Black	\$0
<input type="radio"/> 55056	Frame Color: Purple Death	\$0
<input type="radio"/> 55017	Frame Color: Stealth Black	\$0
<input type="radio"/> 55103	Frame Color: Wrinkle Black	\$0
<input type="radio"/> 25102	Water Bottle and Cage	\$0



Top End Sports LLC.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End Sports LLC unless otherwise noted. Specifications and prices are subject to change without notification.

©2023 Top End Sports LLC

Top End®

Force™ RX Handcycle

HANDCYCLE SPECIFICATIONS:	
Weight:	Approximately 28 pounds
Turning Radius:	18 feet for 360 degree circle, 12 feet on a 3-point turn
Transfer Height:	5" from ground with 2" cushion
Frame Style:	Ovalized 7005 Aluminum tubing which is heat treated for increased strength and stiffness.
Fork:	Removable, Travel ready, with adjustable cranks
Propulsion / Shifting:	11 - 22 speeds with choice of shifting components
Fork Style:	Adjustable length and height options
Back:	Carbon Fiber, adjustable angle from 50 to 70 degrees
Footrest:	Dual adjustable length, includes straps
Hand Pedals:	Ovalized ergonomic with frictionless bearings
Brake:	Caliper brakes standard, optional front disk brake
Wheelbase:	60" to 63" depending on wheel selection. See definitions
Camber:	3 or 0 degree available
Wheels:	20", 26" or 700c wheels with optional upgrades
Axles:	Threaded
Tires:	Standard tires provided. Exact brand varies based on availability
Positioning Straps:	Hook and loop positioning strap included
Cushions:	Seat and back cushion included
Safety / Visibility:	Flag and reflectors included
Chain Guard:	Carbon Fiber chain guard included

Warranty Info:

Top End Sports LLC. warrants the frames and forks when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty

EXCLUSIONS AND DEFINITIONS:	
Top End Force RX Frame and Fork:	The Top End Force RX handcycle is made of 7005 series aluminum. An adjustable carbon fiber backrest and an adjustable bolt-on headrest are included. Speeds vary between 11 and 22 depending on shifting selections.
Assembled vs Unassembled and shipping:	If TE10048 (Assembled) is selected the bike will be mostly assembled upon delivery. Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
Wheel Sets:	Wheel sets include a front wheel, two rear wheels with axles, tires and tubes.
Tires:	Standard tires included.
Fork Style:	Fork adjustment range is about: 19" to 21.5" with a 26" front wheel & 20" to 22.5" with a 700c front wheels
eTap Electronic Shifting:	ETap is a wireless electronic shifting system. Power Plate upgrade required. Limited chainring selections available and only available with C7 wheels.
Power Plates:	Required with eTap shifting package
Omit Wheels:	When an omit wheels option is selected the following will NOT be included. Front or rear wheels, tires, tubes, axles. When a shifting package is selected the cassette from that package will be included.
Wheelbase:	Wheelbase is measured from the middle of the rear axle to the middle of the front hub. When Standard length is chosen wheelbase will be 60" with a 26" front wheel and 61" with a 700c front wheel. When extended length is chosen wheelbase will be 61" with 26" front wheel and 63" with a 700c front wheel.
Braking Options:	The Tektro braking package includes a caliper and cables for the hand brake lever included in the shifting package and flip lever parking brake assembly. The disc brake braking package includes a disc brake rotor, caliper, cables and a Tektro flip lever parking brake assembly.
Omit Shifting:	When omit shifting option is selected the shifter, chainring derailleur, cassette derailleur, cassette, friction shifter for brake caliper, chain and shifter cables will NOT be included. The handpedals and fork will be made at accept trigger shifter components mounted on the right handpedal. A hand brake lever, caliper and cables as well as a flip-lever parking brake will be included. The brand/style is determined by the braking option selected.
25 Degree Handpedal Angle:	Only available with eTap shifting package due to interference. Not available with crank widths of 12"-15".

Note: All Special requests are subject to approval, additional charges, and additional lead time

Special Request 1)	
Special Request 2)	
Special Request 3)	
Special Request 4)	
Special Request 5)	

Top End Sports LLC.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End Sports LLC unless otherwise noted. Specifications and prices are subject to change without notification.

©2023 Top End Sports LLC


Top End®

Force™ RX Handcycle

FORK COLOR:			
<input type="radio"/>		Same as Frame Color	
<input type="radio"/> 55067	N/A	Not applicable with 55066 (buffed frame), must select 55044 (buffed fork)	\$0
<input type="radio"/> 55041	Black	Wet Black	\$120
<input type="radio"/> 55042	Black	Black Prism	\$120
<input type="radio"/> 55040	Black	Stealth Black	\$120
<input type="radio"/> 55033	Blue	Midnight Blue	\$120
<input type="radio"/> 55034	Blue	Royal Blue	\$120
<input type="radio"/> 55035	Blue	Electric Blueberry	\$120
<input type="radio"/> 55046	Gray	Stealth Grey	\$120
<input type="radio"/> 55061	Green	Agriculture Green	\$120
<input type="radio"/> 55031	Green	Grasshopper Green	\$120
<input type="radio"/> 55062	Green	Neon Green	\$120
<input type="radio"/> 55029	Orange	Sunkissed Orange	\$120
<input type="radio"/> 55028	Orange	Simple Orange	\$120
<input type="radio"/> 55059	Orange	Neon Orange	\$120
<input type="radio"/> 55038	Pink	Bubblegum Pink	\$120
<input type="radio"/> 55065	Pink	Neon Pink	\$120
<input type="radio"/> 55063	Purple	Lilac Purple	\$120
<input type="radio"/> 55024	Red	Simple Red	\$120
<input type="radio"/> 55025	Red	Tangerine Red	\$120
<input type="radio"/> 55026	Red	Raspberry Red	\$120
<input type="radio"/> 55058	Red	Neon Red	\$120
<input type="radio"/> 55043	Silver	Silver Metallic	\$120
<input type="radio"/> 55039	White	White Snow	\$120
<input type="radio"/> 55030	Yellow	Safety Yellow	\$120
<input type="radio"/> 55060	Yellow	Neon Yellow	\$120
<input type="radio"/> 55096	Bronze	Las Vegas Bronze	\$120
<input type="radio"/> 55097	Green	Mint Green	\$120
<input type="radio"/> 55098	Pink	Light Pink	\$120
<input type="radio"/> 55099	Grey	Fighter Jet Grey	\$120
<input type="radio"/> 55100	Blue	Tiffany Blue	\$120
<input type="radio"/> 55101	Blue	Carolina Blue	\$120
<input type="radio"/> 55044	Raw	Buffed Aluminum	\$241

ACCESSORIES:		
STANDARD:		
<input checked="" type="radio"/> 25061	Safety Flag	\$0
<input checked="" type="radio"/> 25067	Drafting Bumper	\$0
<input checked="" type="radio"/> 25071	Adjustable Headrest	\$0
<input checked="" type="radio"/> 25038	Mirror, Chain Guard Mounted	\$0
UPGRADABLE:		
<input type="radio"/> 25039	Leg Guard	\$90
<input type="radio"/> 60015	Leg Support, Curved (Pair)	\$157
HYDRATION:		
<input type="radio"/> 25021	Water Bottle and Cage	\$0

DECAL COLORS: (decals are reflective)		
FRAME DECAL COLOR:		
<input type="radio"/> 55071	Black Frame Decals	\$0
<input type="radio"/> 55072	White Frame Decals	\$0
FORK DECAL COLOR:		
<input type="radio"/> 55073	Black Fork Decals	\$0
<input type="radio"/> 55074	White Fork Decals	\$0
FORK COLOR		
<input type="radio"/> 55106	Fork Color: Bomber Orange	\$120
<input type="radio"/> 55076	Fork Color: Custom	\$224
<input type="radio"/> 55068	Fork Color: Flat Black	\$0
<input type="radio"/> 55104	Fork Color: Wrinkle Black	\$0
FRAME COLOR		
<input type="radio"/> 55105	Frame Color: Bomber Orange	\$0
<input type="radio"/> 55066	Frame Color: Buffed Aluminum	\$421
<input type="radio"/> 55102	Frame Color: Flat Black	\$0
<input type="radio"/> 55056	Frame Color: Purple Death	\$0
<input type="radio"/> 55017	Frame Color: Stealth Black	\$0
<input type="radio"/> 55103	Frame Color: Wrinkle Black	\$0
<input type="radio"/> 25102	Water Bottle and Cage	\$0



Top End Sports LLC.

All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Top End Sports LLC unless otherwise noted. Specifications and prices are subject to change without notification.

©2023 Top End Sports LLC