

# Top End® Little Excelerator™ Handcycle



## PROVIDER / CUSTOMER INFORMATION:

Provider Information:		End User Information:	
Date:	Contact Name:	Athlete Name:	
Acct #:	Ship to Address:		
Contact:	City:	State:	
Phone #:	Zip/Postal Code:	Phone:	
Fax/Email:	Email:		
<input type="radio"/> Ship to Provider		<input type="radio"/> Ship to End User	

## ORDER REQUEST TYPE:

☐ Quote
 ☐ Order
 PO #

### BASE MODEL:

☒ 10008 Top End Little Excelerator Handcycle \$4,399

### Assembly: *(see definitions)*

☒ 10048 Assembled \$0

### Shipping:

Shipping Fee  
*(Continental US Only - See definitions)* \$325

### SEAT WIDTH:

☒ 15013 14" Seat Width \$0

### CAMBER ANGLE:

☒ 45004 9 Degree Camber \$0

### WHEELBASE: *(Front to Rear Axle)*

☒ 75074 35" Wheelbase \$0

### WHEEL SETS:

☒ 50120 20" Front Wheel & 20" Rear Wheels with Quick Release Axles \$0

### TIRES:

☒ 50212 20" Treaded Tires \$0

### FOOTREST:

☒ 60002 Dual, Height/Angle Adjustable \$0

### UPHOLSTERY:

#### MATERIAL & COLOR:

☒ 30001 Nylon, Black Upholstery \$0

#### SEAT UPHOLSTERY:

☒ 30014 Adjustable Tension Seat Upholstery \$0

☒ 25060 Hook and Loop Positioning Strap \$0

#### BACK UPHOLSTERY:

☒ 30029 Fold - Over Back Upholstery \$0

#### SEAT CUSHION:

☒ 25063 Seat Cushion \$0



### DRIVE CONTROL:

#### CRANK TYPE:

☒ 65074 Crank Arms: 155mm Long by 12" Wide \$0

#### HAND PEDAL TYPE: *(see definitions)*

☐ 65137 Medium 3D Printed with Power Plates \$0

☐ 65134 Tri-Pin Hand Pedals \$361

#### SHIFTING:

☒ 65060 Twist Activated Shifter with Internal Gearing \$0

☐ 65061 Quad Twist Shifting Adaptor  
*(see definitions)* \$181

#### BRAKING:

☒ 65080 Reverse Activated Drum Brake \$0

☒ 65090 Locking Parking Brake \$0

#### CHAINRING:

☒ 65100 Single 38 Tooth Chainring \$0

### ACCESSORIES:

#### STANDARD:

☒ 25061 Safety Flag \$0

☒ 25062 Chain Guard \$0

#### UPGRADABLE:

#### SAFETY:

☐ 25097 Mirror, Rear View \$37

#### HYDRATION:

☒ 25021 Water Bottle and Cage \$0

# Top End® Little Excelerator™ Handcycle



## FRAME COLOR:

<input type="radio"/>	55018	Black	Wet Black	\$0
<input type="radio"/>	55019	Black	Black Prism	\$0
<input type="radio"/>	55017	Black	Stealth Black	\$0
<input type="radio"/>	55010	Blue	Midnight Blue	\$0
<input type="radio"/>	55011	Blue	Royal Blue	\$0
<input type="radio"/>	55012	Blue	Electric Blueberry	\$0
<input type="radio"/>	55023	Grey	Stealth Grey	\$0
<input type="radio"/>	55053	Green	Agriculture Green	\$0
<input type="radio"/>	55008	Green	Grasshopper Green	\$0
<input type="radio"/>	55054	Green	Neon Green	\$0
<input type="radio"/>	55006	Orange	Sunkissed Orange	\$0
<input type="radio"/>	55005	Orange	Simple Orange	\$0
<input type="radio"/>	55051	Orange	Neon Orange	\$0
<input type="radio"/>	55015	Pink	Bubblegum Pink	\$0
<input type="radio"/>	55057	Pink	Neon Pink	\$0
<input type="radio"/>	55055	Purple	Lilac Purple	\$0
<input type="radio"/>	55001	Red	Simple Red	\$0
<input type="radio"/>	55002	Red	Tangerine Red	\$0
<input type="radio"/>	55003	Red	Raspberry Red	\$0
<input type="radio"/>	55050	Red	Neon Red	\$0
<input type="radio"/>	55020	Silver	Silver Metallic	\$0
<input type="radio"/>	55016	White	White Snow	\$0
<input type="radio"/>	55007	Yellow	Safety Yellow	\$0
<input type="radio"/>	55052	Yellow	Neon Yellow	\$0
<input type="radio"/>	55070	Raw	Buffed Aluminum	\$0
<input type="radio"/>	55090	Bronze	Las Vegas Bronze	\$0
<input type="radio"/>	55091	Green	Mint Green	\$0
<input type="radio"/>	55092	Pink	Light Pink	\$0
<input type="radio"/>	55093	Grey	Fighter Jet Grey	\$0
<input type="radio"/>	55094	Blue	Tiffany Blue	\$0
<input type="radio"/>	55095	Blue	Carolina Blue	\$0



## Custom Color Options



Colors may vary due to monitor settings.  
Color variation in sample represent depth from lighting.

# Top End® Little Excelerator™ Handcycle

HANDCYCLE SPECIFICATIONS:	
<b>Weight:</b>	45 pounds
<b>Turning Radius:</b>	10 feet for 360 degree circle, 7 feet on a 3-point turn
<b>Transfer Height:</b>	17.5" from ground with 2" cushion
<b>Frame Style:</b>	Single main tube, upright position, 1-3/4" X .065" ovalized carbon steel tubing
<b>Seat:</b>	Seat adjust forward and rearward with sliding mechanism
<b>Fork:</b>	Carbon Steel
<b>Propulsion / Shifting:</b>	Internal front hub, twist activated shifter mounted on crank support
<b>Cranks:</b>	155mm adjustable height
<b>Back:</b>	Adjustable height
<b>Footrest:</b>	Dual adjustable height and angle, includes straps
<b>Hand Pedals:</b>	Ovalized ergonomic with frictionless bearings
<b>Brake:</b>	Reverse activated drum and Locking Parking Brake mounted on fork
<b>Wheelbase:</b>	56" axle to axle
<b>Camber:</b>	9 degree
<b>Wheels:</b>	20" metal spoked wheels
<b>Axles:</b>	Quick Release
<b>Tires / Tubes:</b>	20" X 1.75" tires (40 PSI)
<b>Positioning Straps:</b>	Hook and loop strap included
<b>Cushions:</b>	Seat cushion included
<b>Safety / Visibility:</b>	Flag and wheel reflectors included
<b>Chain Guard:</b>	Full length, shock resistant composite material

DEFINITIONS:	
<b>Assembly and Shipping:</b>	Arrives factory assembled, rear wheels and some accessories will need to be installed. Shipping rate is \$325 for all shipments in the continental US. Freight quotes for shipments outside the continental US are available upon request.
<b>Quad Twist Shifter Adaptor:</b>	Three pronged adaptor which allows user to tap for shifting.
<b>Tri-pin Hand Pedals:</b>	3-Point contact hand pedal similar to hand controls used for vehicles. Allows for the hands to be kept in place for people with limited grip strength.

## Warranty Info:

Top End warrants the frames when purchased new and unused to be free from defects in materials and workmanship for a period of three (3) years from the date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty. Top End warrants the upholstered materials (seat and back) and remaining components of this product when purchased new and unused to be free from defects in materials and workmanship for a period of thirteen (13) months from date of purchase from Top End or a dealer, with a copy of the seller's invoice required for coverage under this warranty.

