

T R A D I T I O N A L S E R I E S

Rise



Cambridge
PR-401M




Cambridge
PR-401L




Williamsburg
PR-747

Softouch
with STAIN DEFENSE



Hazelnut




Sterling




Copper



Calypso



Shiraz



Island*

*Std. on
PR-747 only

Traditional Series

The Traditional Series features two models in classic styles that will subtly match any decor. The Cambridge and Williamsburg feature plush, overstuffed backrests for hours of comfort.



Rise.

Roll.

Ride.



GOLDEN
technologies

YOUR LIFE IN MOTION



TRADITIONAL SERIES



Cambridge Cambridge Williamsburg

Model Number	PR-401M	PR-401L	PR-747
Size	Medium	Large	Medium
Weight Capacity	375 Lb.	375 Lb.	375 Lb.
Recline Positions	3	3	3
Standard Fabrics	5	5	6
Fabric Options	All	All	All
Chaise Seat	YES	YES	YES
Back Type	Biscuit	Biscuit	Button
Head Rest Cover	YES	YES	YES
Arm Covers	YES	YES	YES
Battery Back Up	YES	YES	YES
Overall Width	33"	35"	32.5"
Overall Height	41.5"	43.5"	44"
Width Between Arms	20"	22"	21"
Floor to Top of Seat	19"	21.5"	19.5"
Seat Depth	21"	22"	21.5"
Seat to Top of Back	28"	30"	28"
Seat to Top of Arm	7"	8.5"	6"
Seat Angle Raised	20°	20°	15°
Floor to Top of Foot Rest Reclined	22"	23"	23"
Floor to Top of Back Reclined	33"	34"	36"
Back Angle Closed	110°	110°	100°
Back Angle Reclined	145°	145°	130°
Foot Rest Angle Closed	90°	90°	90°
Foot Rest Angle Opened	175°	175°	180°
Distance Required from Wall Reclined	17"	19"	18"
Extension reclined from Back of Seat to end of Foot Rest	41.5"	44"	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.
Removable Back Standard	YES	YES	YES
Smart Tek System Standard	YES	YES	YES
Customizable*	NO	NO	NO
Left Hand Control Option Available	YES	YES	YES
Extra Pocket Option Available	YES	YES	YES
Heat & Massage Option Available	YES	YES	YES
Foot Rest Extension Option Available	YES	YES	YES
Straight Lift Option Available	YES	YES	YES

* See Customize Chart for Dimension limits

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness.
All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68
Listing # E112417
FNAE 80-214A



Golden Technologies • 401 Bridge Street, Old Forge PA 18518 • Tel: (800) 624-6374 • Fax: (800) 628-5165 • www.goldentech.com

SHT-TRAD-082813