

BRAND	MODEL				
-------	-------	--	--	--	--

Ki	Catalyst 5 (*)				
	Catalyst 5Ti (*)				
	Catalyst 5Vx (*)				
	Rogue (+)				
	Rogue ALX (+)				
	Rogue XP (*)				
(+) Chairs have a camber tube diameter of 1-1/8"/28mm					

Quickie	Helium (+)				
	Helium Freestyle (+)				
	Nitrum (+)				
	5R (+)				
	Neon (+)				
	Helix 2 (+)				
	Neon 2 (+)				
	Q55 X				
	Quickie 2				
	QX/QXi				
	QRi (+)				
	Xenon (all models +)				
(+) Chairs have a camber tube diameter of 1-1/4"/31.75mm					

RGK	Tiga (all models)				
	Octane (all models)				
	Hi Lite (all models)				
<i>When ordering "single post fixed center of gravity," it is necessary to request RGK's build-up on the camber tube to accommodate rear camber tube-mounted power assist systems.</i>					

Motion Composite	Apex A (+)				
	Apex C (+)				
	Helio C2				
	Helio C2 Carbon				
	Veloce Carbon				
	Move				
	Helio A6				
	Helio A7				
(+) Chairs have a camber tube diameter of 1-1/4"/31.75mm					

TiLite	ZR (+)				
	ZRA (+)				
	TR (+)				
	TRA (+)				
	ZGX				
	TX (until 2017)				
	Z (+)				
	AEROT (+)				
	X				
	TWIST (+)				
(+) Chairs have a camber tube diameter of 1-1/4"/31.75mm					

Zippie	Zone (+)				
(+) Chairs have a camber tube diameter of 1-1/4"/31.75mm					